


Rišo tempimo sistema su priedais

Patient Positioning
Limb Positioners



Wrist Tower Kit and Accessories

Reference #	Description
72202056	Wrist Tower System

System includes

72202114	Wrist Tower	1. Rišo tempimo stovas
7210577	Traction Accessory	
72200660	Wrist Arthroscopy Accessory with Wrist Plate	2.2. Rišo plokštelė

Rišo stabilizavimo rinkinys.
Finger Traps (Disposable)

72201550	Small/Medium Wrist Stabilization Kits (5 per case)	2.3
72201551	Medium/Large Wrist Stabilization Kits (5 per case)	
72201552	Large/XLarge Wrist Stabilization Kits (5 per case)	2.4

Additional Accessories

72200643	US Bed Rail (2)
72200675	Euro Bed Rail (2)
72200666	Wrist Plate replacement (radiolucent – shows a shadow)
72200665	Finger Trap Sizing Ring
72202309	Wrist Tower Tray
72202310	Wrist Accessory Tray

SPIDER2 and SPIDER Limb Positioners – Wrist



Reference #	Description
72203299	SPIDER2 Limb Positioner
7210569	SPIDER Limb Positioner
7210570	Piggy Back Connector, 2 per box, (required for 7210569)
72200660	Wrist Arthroscopy Accessory With Wrist Plate
72200664	Wrist Counter Traction Post
7210577	Traction Accessory

Finger Traps (Disposable)

72201550	Small/Medium Wrist Stabilization Kits (5 per case)
72201551	Medium/Large Wrist Stabilization Kits (5 per case)
72201552	Large/XLarge Wrist Stabilization Kits (5 per case)

TENET Wrist Tower



Wrist distraction,
without the heavy lifting




The TENET Wrist Tower is easy to set up and adjust without using weights to apply distraction. It also provides excellent stability, allowing you to operate with confidence.

Flexible design

- Low-profile design ensures excellent C-arm access
- Includes removable hand plate
- Offers three sizes of finger traps (trys dydžiai pirštų fiksatorius)
- Allows for ulnar and radial deviation and full rotation of patient's arm

Easy to use

- Fine traction system eliminates weights
- Allows intraoperative positioning within the sterile field



Convenient

- Simple to set up, break down, and store
- Easily sterilized and fully autoclavable
- Compatible with most OR tables

Riešo tarpinio stovas fiksavimas
priė bet kokių operacinių stalo.
Secure positioning


- Provides stable distraction, up to 30 lbs

13 lb.

Western Elbow Positioner

The Western Elbow Positioner fits into a Clark Clamp and provides positioning options for lateral and prone elbow procedures.

- Low profile device provides exposure to the operative elbow. The slim design of the pad allows for easy draping
- Support Bar can be positioned close to the patient's chest which provides more room around the elbow to work
- The Bicep Strap secures the patient's arm to the support pad





Western Elbow Positioner

Wrist Tower

The Wrist Tower is designed to be used as a standalone tower that attaches to most OR tables and is accessible from within the sterile field.²¹⁻²³ With adjustable traction and flexible positioning options, the Wrist Tower provides access and stability.^{20, 22, 23}

- **Flexible Design**
Low-profile design to help provide proper C-arm access
- **Traction**
Measurable traction without the interference of cumbersome weights¹⁸
- **Positioning**
Ulnar and radial deviation and full rotation of patient's arm
- **Convenient**
Easily sterilized: fully autoclavable²²
Compatible with our Wrist Stabilization Kits offered in small, medium and large sizes

sterilized in pillow sterilization




Wrist Tower

72

TENET Wrist Tower

Gaining access to the operative hand and wrist is simple with the TENET Wrist Tower.

Tempimo sistema fiksuojama prie riešo tempimo stovo.



Ordering Information

72202056 TENET Wrist Tower System

System Includes:

<u>72202114</u>	<u>TENET Wrist Tower</u>
<u>7210577</u>	<u>Traction Accessory</u>
<u>72200660</u>	<u>Wrist Arthroscopy Accessory with Wrist Plate</u>


Finger Traps (Disposable)

<u>72201550</u>	<u>Small (case of 5)</u>
<u>72201551</u>	<u>Medium (case of 5)</u>
<u>72201552</u>	<u>Large (case of 5)</u>

Additional Accessories Sold Separately

<u>72200643</u>	<u>US Bed Rail (2)</u>
<u>72200675</u>	<u>EUR Bed Rail (2)</u>

1-1.
 Fitavimmo
 prije operacijo
 stola komponenta.



Tempinuo sistema

7210577

RiSO ortros-
 kopijos prietas
 su RiSO
 plotele.
 72200660

Manufactured by TENET Medical Engineering and distributed by Smith & Nephew.

Endoscopy
 Smith & Nephew, Inc.
 Andover, MA 01810
 USA

www.smith-nephew.com
 T +1 978 749 1000
 US Customer Service: +1 800 343 5717
 International Customer Service: +1 978 749 1140

®Trademarks of Smith & Nephew. Certain marks Reg. US Pat & TM Off. All other trademarks acknowledged.
 ©2009 Smith & Nephew. All rights reserved.
 Printed in USA. 09/09 10600434 Rev. A

74