

Scatole portavetrini

Slide boxes **Dėžutės objektyvinių stiklelių laikymui**

| | Cod. | Descrizione Description | Q.ta Q.ty |
|--------------|------|--|-------------|
| 11 PD | 0544 | scatole portavetrini plastica da 25 posti slide boxes for 25 slides | 10 pcs |
| 10 PD | 0545 | scatole portavetrini plastica da 50 posti slide boxes for 50 slides | 10 pcs |
| 9 PD | 0516 | scatole portavetrini plastica da 100 posti slide boxes for 100 slides | 1 pcs |
| | 0520 | scatole portavetrini security da 25 posti slide boxes for 25 slides - new model | 10 pcs |
| | 0562 | scatole portavetrini security da 50 posti slide boxes for 50 slides | 10 pcs |
| | 0521 | scatole portavetrini security da 100 posti slide boxes for 100 slides - new model | 1 pcs |
| | 0350 | scatole portavetrini security da 12 posti slide boxes for 12 slides - new model | 1 pcs |



Cartella portavetrini in cartone

Cardboard slide mailers

| | Cod. | Descrizione Description | Q.ta Q.ty |
|--|------|---------------------------|-------------|
| | 3341 | 1 posto/1 place | 250 pcs |
| | 0549 | 2 posti/2 places | 50 pcs |
| | 0551 | 3 posti/3 places | 100 pcs |
| | 0669 | 4 posti/4 places | 50 pcs |
| | 3149 | 5 posti/5 places | 50 pcs |
| | 3378 | 6 posti/6 places | 10 pcs |
| | 0552 | 10 posti/10 places | 50 pcs |
| | 0550 | 20 posti/20 places | 50 pcs |



Vassoio in cartone

Cardboard trays

| | Cod. | Descrizione Description | Q.ta Q.ty |
|--|------|--|-------------|
| | 3151 | 10 posti/10 places | 10 pcs |
| | 0755 | 20 posti/20 places | 10 pcs |
| | 0575 | 20 posti senza divisorio/20 places without partition | 10 pcs |
| | 0815 | 20 posti con divisorio centrale/ 20 places with central partition | 10 pcs |



Tissue Processing Cassettes

32.1 PD Mega kasetēs su papildoma vieta standartinei kasetei

Supa Mega Mothership Cassette



With an internal cavity a third of the depth of the conventional Supa Mega, CellPath's Supa Mega Mothership provides an innovative solution for the processing, embedding and tracking of large tissue specimens.

In accordance with recommendations of HER2 testing in the UK¹ and ASCO/CAP in the USA^{2,3}, by reducing specimen thickness and shortening processing time, whilst ensuring high quality tissue processing.

The Supa Mega Mothership cassette uses an innovative permanent docking bay for the majority of standard sized tissue processing cassettes* printed with ID bar codes and/or test information. This enables the specimen to be entered onto a LIMS guaranteeing its accurate and reliable tracking.


Dimensions: 52.5mm x 75mm x 18.3mm



Matmenys

KEY BENEFITS

- 1 Enables specimens to be entered onto Laboratory Information Management System (LIMS)
- 2 Better reagent flow and improved turnaround times
- 3 Casting done in 30 minutes at -8°C
- 4 Mothership Cassette can accommodate the majority of available standard sized cassettes
- 5 Large colour range available with bespoke colours available on request



The new innovative latch design allows the Supa Mega Mothership cassette to house most available standard sized tissue cassettes.

The Supa Mega Mothership Cassettes is available in a range of colours



| Item Code | Description (Packaged Loose) | Quantity | Colour |
|--------------|-------------------------------|----------|--------|
| EAO-0102-02A | Supa Mega Mothership Cassette | 100 | White |
| EAO-0104-02A | Supa Mega Mothership Cassette | 100 | Green |
| EAO-0106-02A | Supa Mega Mothership Cassette | 100 | Blue |
| EAO-0108-02A | Supa Mega Mothership Cassette | 100 | Pink |
| EAO-0109-02A | Supa Mega Mothership Cassette | 100 | Yellow |
| EAO-0110-02A | Supa Mega Mothership Cassette | 100 | Red |

1 HER2 Testing in the UK: Further update to Recommendations. R.A. Walker et al. J. Clin. Pathol. 2008; 61: 818-824.
 2 American Society of Clinical Oncology/College of American Pathologists Guideline Recommendations for Immunohistochemical Testing of Estrogen and Progesterone Receptors in Breast Cancer. Hammond et al. Arch Pathol Lab Med-2010
 3 ASCO/CAP Guideline Recommendations for Human Epidermal Growth Factor Receptor 2 Testing in Breast Cancer. A.C. Wolff et al. Arch Pathol Lab Med - Vol 131, January 2007

* The Supa Mega Mothership has been verified for use with CellPath's System III, System II+, System II Hex and PI cassettes. Alternative cassettes can also now be housed in the docking system, although these require validation by the user/laboratory.

FiloSlide 100 Supa Mega Slides

The sturdy FiloSlide plastic box offers safe and cost-effective microscope slide storage for up to 100 Supa Mega slides per box. They can be stacked for easy storage and feature a snap closure that seals out dirt and dust. They are also splash-proof with an index card included for easy specimen identification.

Dimensions: 231mm x 191mm x 58mm



KEY BENEFITS

- 1 Store up to 100 Supa Mega slides
- 2 Splash-proof
- 3 Index card included
- 4 Stackable for easy storage

| Item Code | Description | Quantity | Colour |
|--------------|--|----------|--|
| WEQ-0906-00A | FiloSlide 100 Box for Supa Mega Slides | 1 | Blue  |

Cassette Boxes

32.2 PD Archyvavimo dėžutės Mega blokų kasetėms

BlocFile 1 Supa Mega Cassette Box

BlocFile 1 Supa Mega Cassette boxes are an economical drawer system for Supa Mega cassette storage. It is manufactured from high-density fibreboard which ensures its strength, rigidity and durability. The box can be supplied flat-packed or ready assembled.

Dimensions: 400mm x 240mm x 95mm - holds 80 Supa Mega cassettes

Talpina apie 80vnt. kasečių



KEY BENEFITS

- 1 Supplied assembled or flat-packed
- 2 Manufactured from recycled high-density fibreboard

| Item Code | Description | Quantity |
|--------------|---|----------|
| WCB-1100-08F | BlocFile 1 - 400mm Deep Supa Mega Cassette — Flatpack | 25 |



32.3 PD Archyvavimo dėžutės Mega stikleliams

BlocStor 3 Supa Mega Box

The BlocStor 3 Supa Mega Cassette Box is manufactured from recycled fibreboard and provides a cost-efficient archiving and storage solution. It is designed for strength and rigidity and the long-term archiving of Supa Mega tissue blocks and slides. It also maximises the utilisation of storage space.

The BlocStor 3 is supplied with a pre-printed archive label for clear identification and rapid retrieval. BlocStor boxes are available assembled or flat-packed.

Dimensions: 395mm x 240mm x 95mm - holds 70 Supa Mega cassettes

| Item Code | Description | Quantity |
|--------------|--|----------|
| WCB-0500-08F | BlocStor 3 Supa Mega Cassette — Flatpack | 25 |



KEY BENEFITS

- 1 Manufactured from recycled material
- 2 Supplied with a pre-printed archive label
- 3 Cost effective archive and storage

Archive Boxes

CellNass SmartStor Block and Slide Archive Boxes

Manufactured from recycled fibreboard, the SmartStor Block and Slide Archive Boxes are an economical solution to the long-term archiving of routine slides.

They are compatible with the CellNass SmartStor 2 cabinet and SmartStor 8 cabinet solution.

SmartStor Supa Mega Block Boxes accommodate 120 Supa Mega blocks; 256 Supa Mega Slim Blocks or 160 Supa Mega Mothership blocks.

SmartStor Supa Mega Slide Boxes accommodate 3,360 Supa Mega™ Slides.

| Item Code | Description | Quantity |
|--------------|---|----------|
| WSD-0400-08A | CellNass SmartStor Supa Mega Block Archive Box — Flatpack | 25 |
| WSD-0500-08A | CellNass SmartStor Supa Mega Slide Archive Box — Flatpack | 25 |



KEY BENEFITS

- 1 Supplied assembled or flat-packed
- 2 Made from recycled material
- 3 Cost-effective archiving solution

Slide Rack Adapter - Stainless Steel

The Slide Rack Adapter enables the user to stain both Mega and Supa Mega slides on automated staining equipment bringing convenience, efficiency and cost savings all at the same time. This makes the Slide Rack Adapter ideal for all laboratories looking to improve their process.




KEY BENEFITS

- 1 Stainless steel for a lifetime of use **nerūdijančio plieno**
- 2 Enables a LEAN staining process
- 3 Use with Mega, Supa Mega and potentially any bespoke sized slide

skirtas talpinti didelius stiklelius

| Item Code | Description | Quantity |
|--------------|--------------------------------------|----------|
| RMC-2000-63A | Slide Rack Adapter - Stainless Steel | 2 |

Slide Rack not included.




Slide Mailers

SlideRite Mega and Supa Mega Slide Cardboard Mailer

The robust cardboard mailers are constructed from one-piece moulding that provides a high level of protection for Mega and Supa Mega slides.

KEY BENEFITS

- 1 Economical
- 2 Reusable
- 3 Tough and durable



Economical and reusable

| Item Code | Description | Quantity |
|--------------|--|----------|
| WEX-0400-01A | SlideRite Mega Slide Cardboard Mailer - fits slides up to 76mm x 39mm | 10 |
| WEX-0500-01A | SlideRite Supa Mega Slide Cardboard Mailer - fits slides up to 76mm x 52mm | 10 |



GŁUCHOŁASKIE ZAKŁADY PAPIERNICZE

Aligninas 36 PD

- Be chloro
- Bekvapis
- Nedulkantis
- Baltumas 80%
- Tankis 28 g/m²
- pH 7
- Balintas;
- Higroskopiškas 8g/g pagal EN ISO 12625-8:2011
- Skendimo laikas 3s
- Lapai 40x60cm
- Pakuotėse po 1 kg, 5 kg





GŁUCHOŁASKIE ZAKŁADY PAPIERNICZE

Cellulose wadding in sheets

36 PD

- Chlorine free
- Scentless
- Lint free
- Whiteness is 80%
- Grammage is 28 g / m²
- pH 7
- Bleached;
- Hygroscopic 8g/g regarding EN ISO 12625-8:2011
- Sinking time 3s
- Lapai 40x60cm
- Packages: 1 kg, 5 kg



Contenitori per campioni con tappo a vite e sottotappo *Specimen containers with screw cap and inner stopper*

CE IVD

Materiale/Material: HDPE (Polietilene alta densità/High density polyethylene)

| Cod. | Descrizione Description | Q.ta Q.ty |
|------|---------------------------|-------------|
| 2434 | 50/60 ml | 500 pcs |
| 0143 | 70 ml | 500 pcs |
| 0144 | 120 ml | 500 pcs |
| 0145 | 250 ml | 250 pcs |
| 0146 | 500 ml | 240 pcs |
| 0147 | 1000 ml | 150 pcs |
| 0148 | 2000 ml | 74 pcs |



46PD Indas 250ml

Contenitori per biopsie

Biopsy containers with pressure cap

CE IVD

Contenitori con tappo a pressione di piccole dimensioni per il trasporto di biopsie o pezzi operatori di dimensioni ridotte.

Materiale/Material: PS (Polistirolo/Polystyrene) **Pagamintas iš PE;**

Histologinių mėginių transportavimo indas
Be etiketės su sandariai užspaudžiamu/užsukamu dangteliu.


| Cod. | Descrizione Description | Q.ta Q.ty |
|------------------|-----------------------------|-------------|
| 0130 | 10 ml | 100 pcs |
| 0131 | 20 ml | 100 pcs |
| 0132 | 30 ml | 100 pcs |
| 0133 | 50 ml | 600 pcs |
| 0134 | 75 ml | 100 pcs |
| 39PD 0135 | 100 ml 100 ml talpos | 100 pcs |
| 0136 | 150 ml | 100 pcs |
| 0137 | 200 ml | 100 pcs |
| 46PD 0138 | 250 ml 250 ml talpos | 100 pcs |





MEDICAL

Features



Microtome Blade Ultimate Tip at 5000 Power

- Evenly matched and aligned blade edges
- Unique coating and uniform hone on blade edges provide consistent and reliable performance
- Good for general paraffin and cryostat use
- Low Profile .010"
- High Profile .012"
- Round corner blades
- New improved dispenser
- Trial packs available

Target Market

Histology Lab

42 PD Mikrotominiai peiliukai
MicroCoat[®] Technology



Plonų parafininių audinių pjūvių kokybiškam pjovimui

★ **Made in U.S.A.**

| ITEM CODE | PRODUCT DESCRIPTION | PACKING | SHIP UNIT | QTY PER SHIP UNIT | SHIP WT. lbs. | DIMENSIONS LxWxH in. |
|-----------|---|------------------------|-----------|-------------------|---------------|----------------------|
| 72-0110 | Blade, Personna Plus low profile .010" | 50 BL/DP 10 DP/CASE | CS | 500 | 2.50 | 5.25 x 4.75 x 7.5 |
| 72-0114 | Blade, Personna Plus high profile .012" | 50 BL/DP 10 DP/CASE | CS | 500 | 2.63 | 5.25 x 4.75 x 7.5 |
| 72-0133 | Blade, Personna Plus low profile .010" | 10 BL/DP 10 DP/CASE | CS | 100 | 2.38 | 5.25 x 4.75 x 7.5 |
| 72-0134 | Blade, Personna Plus high profile .012" | 10 BL/DP 10 DP/CASE | CS | 100 | 2.41 | 5.25 x 4.75 x 7.5 |

To see the full line of Personna medical products: <http://medical.atblades.com>

Email: david.piper@durham-duplex.com

Features

Benefits

universalūs, žemo profilio, 80mm ilgio

Low Profile
 .010 x .315 x 3.15 in (.25 x 8 x 80mm)

High Profile
 .012 x .554 x 3.15 in (.30 x 14.1 x 80mm)

- Longer length to maximize cutting area
- For general paraffin and cryostat use

ašmuo suapvalintais galais

Round corner blade

- Improves operator handling safety

pagaminti iš nerūdijančio plieno

High quality stainless steel

- Smoother cutting edge
- Greater resistance to corrosion

PTFE polymer coating

- Smoother cuts with less compression
- Greater section clarity
- Increased lubricity and longevity

Chrome hard coating

- Greater durability

3 Facet Edge

- Finer blade tip
- Optimal cutting angles

Made in USA

- High quality, precision blades

Uniquely designed dispenser:
 Blade Follower
 Transparent construction
 Rubber grip pusher
 Impact resistant materials
 Built in blade disposal

- One hand dispensing
- Stabilizes blades for consistent dispensing
- Visibility of blade supply
- Improved comfort and safety
- Operator safety
- Convenient disposal of used blades

Warning: Blades are extremely sharp. Handle with care.
 Used sharps are contaminated. Disease may occur from
 contact or injury. Dispose of in puncture resistant container.
 Refer to OSHA or local standards setting organizational
 guidelines for handling and disposal of sharps.



SHEFFIELD • ENGLAND

Razor sharp since 1910

Durham Duplex
312-314 Petre Street
Sheffield, S4 8LT, UK

Tel: +44 (0)114 243 2313
Export: +44 (0)114 251 8468
export@durham-duplex.com

23 May 2019

Confidential: not to be disclosed without written permission of Durham Duplex.

We hereby advise technical details of the blade listed below

B575 002 A01

- o AccuTec Blades' 72-0110
- o Low profile Personna Plus disposable microtome blade
- o 3.15"/80mm long
- o 0.25"/8mm high
- o 0.010"/0.0254mm thick
- o Triple facet with inclusive bevels approx. 11°, 21° & 35° ašmens kampas-35°
- o Coated stainless steel
- o 50 blades per dispenser, 10 dispensers per case



Sales Manager Europe
Durham Duplex

44 PD Stikleliai imunohistocheminėms reakcijoms



TOMO® Adhesive Microscope Slide

ISO 13485 Certified. | IVDR CE marked.

Matsunami microscope slides are manufactured in Japan using optically clear, exceptionally flat, low fluorescence glass and finished with 90° and clipped corners.

90° kampas; Baltos spalvos šlifluotas galas



Matmenys: 75x25x1mm



Hidrofiliniai objektiniai stikleliai, užtikrinantys stiprų sukibimą su audiniu ir vienodą reagentų pasiskirstymą

Matsunami TOMO® hydrophilic Adhesive slides offer superior tissue adhesion, staining performance and slide surface fluidics. Proprietary slide coating

ensures consistent performance even after exposure to heat stress and humid environmental conditions. Optimized for use with thermal printers.

Recommended for use with automated IHC platforms.

Atsparūs karščio poveikiui ir drėgnai aplinkai;
Optimizuoti naudojimui su šiluminiais spausdintuvais;

- Using optically clear, exceptionally flat, low fluorescence glass
pagaminti iš optiškai švaraus, išskirtinai plokščio stiklo
- Superior tissue adhesion, staining performance and slide surface fluidics
- Optimized for use in various IHC applications
rekomenduojami naudoti Ventana Benchmark automatinėmis IHC sistemomis

Dokumentą elektroniniu parašu
pasirašė TOMA ŠREIDERĖ
Data: 2024-06-28 12:25:37
Paskirtis: 721598
Vieta: Energetikų g. 8
Kontaktinė informacija:
Administratorė