

EMOSTAN®

Highly absorbing
medical pad
 designed to plug
 and stop bleeding
 at the end of
 the infusional
 treatment



Description

Highly absorbent pad with compressive hemostatic action, specific for post-dialysis arteriovenous access care to be used after the removal of the fistula needles.
 Also suitable for patients who require intravenous treatment and blood transfusion.

Intended use

Post dialysis dressing to stop the bleeding after the removal of the fistula needles.

Indication

Designed to stop bleeding and to reduce the hemorrhagic risk, thanks to the compression on the needle holes and to work as a barrier against the risk of infection.
 The non-adherent, biocompatible layer on the pad does not stick to the skin and guarantees a non-traumatic removal.

Available Models

CODE	PAD DIMENSION	PAD THICKNESS	PCS/PACK	PACKS/BOX
S0210	∅ 2.6 cm	0.7 cm ± 10%	2	500
S0210-1			1	700
S0211	∅ 3.5 cm		2	500
S0211-1			1	700
S0212	3.5 x 4.5 cm		2	300
S0212-1			1	500