

Olympus EndoTherapy - Your Picks

HEMOSTASIS

SINGLE-USE APPLICATOR




- Eliminates the reprocessing procedure after use
- Cost-efficient single-use clip-fixing device for multiple clipping
- Fast and easy clip reloading (push-pull action)
- Clips reloadable during one clinical case
- Rotation function for easy targeting
- Compatible with all the different clip types of EZ Clip (HX-610 series)
- Three working lengths available

Article Number	Article Name	Quantity	Min. Working Channel ø	Working Length	For Use With	Area of Compatibility	MR Icon (x)
N5432930	HX-810UR	10	● 2.8 mm	2300 mm	HX-610 series		X

Olympus EndoTherapy - Your Picks

CLIPS FOR SINGLE-USE APPLICATOR



- Fast and effective mechanical hemostasis
- Eight types of EZ Clips to suit all applications
- Rotatable
- Sterile cartridge-housed clips
- 90° and 135° jaw angles available
- Short tail
- Colored cartridges to easily identify each clip
- Additional uses include closure of mucosal defects, marking, and anchoring
- MR conditional

Article Number	update	Article Name	Quantity	Clip Arm Length	Jaw Angle (A)	Features	Area of Compatibility	MR Icon (x)
A,B N6012430	x	HX-610-135	40	~ 7,5 mm	135°	For pathological mucosa, deep ulcerous lesions, cancerous tissue		X
A,B N6012830	x	HX-610-135L	40	~ 9 mm	135°	Wide clip opening; for post-EMR and post-polypectomy clipping in polyps with big stalk; for tangential approach		X
A,B N6012330	x	HX-610-090	40	~ 7,5 mm	90°	For healthy, elastic mucosa; Mallory-Weiss syndrome; open vessels; small perforations; marking		X
A,B N6012730	x	HX-610-090L	40	~ 9 mm	90°	For healthy, elastic mucosa; large-scale lesions; large adenomatous pedicles before and after polypectomy		X