

Olympus EndoTherapy - Your Picks


ELECTROSURGERY

ESG-150 SYSTEM



A

B



- Basic and advanced user levels for adjusted usability
- Touch screen
- Storage function for settings
- 29 languages
- Aligned with ESG-300: User interface and HF settings are equal

	Article Number	update	Article Name	Quantity	Description	For Use With	Area of Compatibility
2.2.7	B A02915A		MAJ-860	1	Active cord for EndoTherapy instruments, 8 mm	ESG-150	
2.2.8	A WA99510A	x	Cable for neutral electrode	1	4.5 m cable	ESG-150	