

Scalpel Blades

Supplied in packs of 100 pieces, sterile



10
46.014.10



11
46.014.11



12
46.014.12



12d
46.014.12D



15
46.014.15



15c
46.014.15C



Safety Blade Remover


46.005.65

This practical tool allows safe replacing of scalpel blades.




1

Place the blade remover



2

Press the blade remover



3

Lever the blade upwards



4

Pull the blade forward and throw it in a containment



46.007.00 Scalpel Handle **ZEPF-Line**, straight, for blades # 10 -15, 14.5 cm



46.007.01 Scalpel Handle, **ZEPF-Line**, angled, for blades # 10 -15, 14.5 cm



46.007.40 Parallel Blade Holder, **ZEPF-Line**, distance 1.5 mm

46.007.45 Parallel Blade Holder, **ZEPF-Line**, distance 2.0 mm