

DRQ 14/21 REV.02 2024/09/20

Description:	Sterilization control tape
Final dimensions:	19 cm x50 m
Composition:	Crepe Paper,Indicator strips.
Latex:	Free
PVC:	Free
Re use:	No
Sterilize by:	Intended for: steam
CMR:	Free
Endocrine Disruptors:	Free
Thickness:	0,11±0,01 mm
Shelf life:	60 months
CFU/g	N.A.
Type ISO 11 140-1	Type I
Note	Free heavy metals
Peel Adhesion	450±80 g/25 mm



Nr and Revision Date:	21 -- 2024/10/23
Applicable standards:	EN ISO 11140-1:2014 ;

Important Notes: Suitable for steam sterilization 121°C or 134°C

PACKING					
Level	Units	Type	Dimension	Stamp	Labels
1	1	Tape 19mmx50m			
2	48	Carton box			

EAN Carton Box/Caixa:05605622007724

Storage instructions	Store in a dry place at room temperature and protect from sunlight. The product has a good stability and therefore, normal conditions of devices storage shall apply.
Waste instructions	Waste: No special requirements. The waste must be disposed according to local legal regulations for the disposable of medical waste.
Intended purpose	Process indicator type 1 for moist heat sterilization process. After an effective full sterilization cycle, the indicator ink will change from the initial color (pink) to black colour. It should be exposed to steam sterilization conditions for 134°C-3.5 minutes or 121°C -15 minutes.
<i>Any modification of the product as presented or any use out of its intended purpose are not of the manufacturer responsibility</i>	

Non-harmonized symbols used in article labeling

Environment - Packaging waste	Weight (gr) of packing material used in 100 pcs	12500
--------------------------------------	---	-------

ISSUED BY: Quality Management Department, 2024/12/06

CONFIDENTIAL: The information contained herein is confidential and belongs to Bastos Viegas, s.a. It is provided for specific purposes and in connection with established agreements. It is not permitted disclosure, reproduction, copying or distribution, in whole or in part, for other purposes or to any third parties.