

	 Search	Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Diameter inch	Length cm	Taper Length cm

[Home](#) / [PRODUCTS](#) / [AMPLATZ SUPPORT WIRE GUIDE](#)

[Share this page](#)

Amplatz Support Wire Guide

[SPECIFICATIONS](#) [IMAGES](#)

Intended for diagnostic and interventional procedures.

This product line is serviced by the following clinical division(s):
[Aortic Intervention »](#)
[Interventional Radiology »](#)
[Peripheral Intervention »](#)

Instructions for Use (IFUs) linked below are the most current version associated with the product. To find an IFU associated with a specific Lot Number, please visit <https://ifu.cookmedical.com/>.

[LEARN HOW TO ORDER »](#)
Set up a customer account and order this or other Cook products

[ASK A SUPPORT QUESTION »](#)
Ask Customer Support about our products or services

[FIND PRODUCT TRAINING »](#)
Connect with our Vista® Education & Training program

[VIEW YOUR ACCOUNT INFORMATION »](#)
Register your customer number on our Customer Portal

Feedback

Specifications

[COLLAPSED](#) [EXPANDED](#)

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Diameter inch	Length cm
Stiff					
Straight - PTFE-Coated Stainless Steel					

Cookie Settings

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Diameter inch	Length cm	Taper Length cm
G27032	THSF-35-145-ASG		N/A	.035	145	
G28781	THSF-35-180-ASG		N/A	.035	180	
G27035	THSF-35-260-ASG		N/A	.035	260	
Extra-Stiff Straight - PTFE-Coated Stainless Steel						
G04047	THSF-25-260-AES		N/A	.025	260	
G04003	THSF-35-80-AES		N/A	.035	80	
ADDITIONAL SPECS						
Description	-					
Floppy Tip Length cm	3					
Tip Curve Radius mm	—					
G03095	THSF-35-145-AES		N/A	.035	145	
G23638	THSF-35-145-AES-ST		N/A	.035	145	
G03330	THSF-35-180-AES		N/A	.035	180	
G23639	THSF-35-180-AES-ST		N/A	.035	180	
G03562	THSF-35-260-AES		N/A	.035	260	
G23640	THSF-35-260-AES-ST		N/A	.035	260	
G04001	THSF-35-300-AES		N/A	.035	300	
G23641	THSF-35-300-AES-ST		N/A	.035	300	

Extra-Stiff


[Cookie Settings](#)

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Diameter inch	Length cm	Taper Length cm
Curved - PTFE-Coated Stainless Steel						
G05067	THSCF-25-260-3-AES		N/A	.025	260	
G03750	THSCF-32-260-3-AES		N/A	.032	260	
G04237	THSCF-35-80-3-AES		N/A	.035	80	
ADDITIONAL SPECS						
Description	-					
Floppy Tip Length cm	3					
Tip Curve Radius mm	3					
G03203	THSCF-35-145-3-AES		N/A	.035	145	
G03565	THSCF-35-180-3-AES		N/A	.035	180	
G03460	THSCF-35-260-3-AES		N/A	.035	260	
Ultra-Stiff						
Straight - PTFE-Coated Stainless Steel						
G28786	THSF-35-90-AUS1		N/A	.035	90	
G28787	THSF-35-90-AUS2		N/A	.035	90	
G19850	THSF-35-145-AUS1		N/A	.035	145	
G19852	THSF-35-145-AUS2		N/A	.035	145	
G19865	THSF-35-180-AUS1		N/A	.035	180	
G20131	THSF-35-180-AUS2		N/A	.035	180	
G19851	THSF-35-260-AUS1		N/A	.035	260	

[Cookie Settings](#)

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Diameter inch	Length cm	Taper Length cm
Ultra-Stiff						
Curved - PTFE-Coated Stainless Steel						
G20171	THSCF-35-60-3-AUS1		N/A	.035	60	
G20302	THSCF-35-90-3-AUS1		N/A	.035	90	
ADDITIONAL SPECS						
Description		-				
Floppy Tip Length cm		2.5				
Tip Curve Radius mm		3				
G28761	THSCF-35-90-3-AUS2		N/A	.035	90	
G19867	THSCF-35-145-3-AUS1		N/A	.035	145	
G19853	THSCF-35-180-3-AUS2		N/A	.035	180	
G19868	THSCF-35-180-3-AUS1		N/A	.035	180	
G20170	THSCF-35-260-3-AUS1		N/A	.035	260	

 Image Library





Order Number

Reference Part Number

Instructions
for Use (IFU)

MR
Status

Diameter
inch

Length
cm

Taper
Length
cm

Products
Patient Resources
Support
About
Divisions
Contact
Education & Training
Customer Portal

Privacy Statement
Terms of Use
Human Rights Responsibility Statement
Cookie Compliance
Sitemap

Sign Up. Get Email about News, Products, and Events from Cook Medical.

SIGN UP

Not all products shown on this website may be approved or available in all jurisdictions.
Consult with your local Cook representative, distribution company or customer support centre for details.