

Search for products, documents, pages, or articles



HOME / PRODUCTS / FIXED-CORE WIRE GUIDE

Share this page

Fixed-Core Wire Guide

SPECIFICATIONS IMAGES

18 p.d Biliarinis drenavimo rinkinys

Used to assist in inserting catheters, introducer sheaths, and other medical devices in vascular applications.

This product line is serviced by the following clinical division(s):
Interventional Radiology »
Peripheral Intervention »

Instructions for Use (IFUs) linked below are the most current version associated with the product. To find an IFU associated with a specific Lot Number, please visit <https://ifu.cookmedical.com/>.

LEARN HOW TO ORDER »
Set up a customer account and order this or other Cook products

ASK A SUPPORT QUESTION »
Ask Customer Support about our products or services

FIND PRODUCT TRAINING »
Connect with our Vista® Education & Training program

VIEW YOUR ACCOUNT INFORMATION »
Register your customer number on our Customer Portal

Feedback

Specifications

COLLAPSED EXPANDED

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Diameter inch	Length cm	Tap cm
--------------	-----------------------	----------------------------	-----------	---------------	-----------	--------

Safe-T-J® Tip
PTFE-Coated Stainless Steel

Cookie Settings

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Diameter inch	Length cm	Tip Length cm
G02160	TSCF-18-180-3		N/A	.018	180	7
G00455	TSCF-21-260-3		N/A	.021	260	7
G00474	TSCF-25-50-3		N/A	.025	50	7
G00466	TSCF-25-145-3		N/A	.025	145	7
G00471	TSCF-25-260-3		N/A	.025	260	7
G00488	TSCF-32-145-3		N/A	.032	145	6
G00524	TSCF-35-50-3		N/A	.035	50	6
G00529	TSCF-35-80-3		N/A	.035	80	6

ADDITIONAL SPECS

Description	-
Floppy Tip Length cm	2.5
Tip Curve Radius mm	3

G00511	TSCF-35-145-3		N/A	.035	145	6
G00509	TSCF-35-145-15		N/A	.035	145	6
G00902	TSCF-35-180-3		N/A	.035	180	6
G02481	TSCF-35-180-15		N/A	.035	180	6
G01505	TSCF-35-220-3		N/A	.035	220	6
G00517	TSCF-35-260-3		N/A	.035	260	6
G00541	TSCF-38-145-3		N/A	.038	145	6

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Diameter inch	Length cm	Tip Length cm
Heavy Duty with a Safe-T-J Tip						
PTFE-Coated Stainless Steel						
G04069	THSCF-35-80-3		N/A	.035	80	7

ADDITIONAL SPECS


Description	Value
Floppy Tip Length cm	2.5
Tip Curve Radius mm	3

G00412	THSCF-35-145-3		N/A	.035	145	7
G01197	THSCF-35-260-3		N/A	.035	260	7
Straight Tip						
PTFE-Coated Stainless Steel						
G00593	TSF-18-50		N/A	.018	50	7
G00587	TSF-18-145		N/A	.018	145	7
G00590	TSF-18-260		N/A	.018	260	7
G00609	TSF-21-50		N/A	.021	50	7
G00602	TSF-21-145		N/A	.021	145	7
G00617	TSF-25-145		N/A	.025	145	7
G00638	TSF-32-145		N/A	.032	145	6
G00661	TSF-35-50		N/A	.035	50	6
G00664	TSF-35-80		N/A	.035	80	6

[Cookie Settings](#)

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Diameter inch	Length cm	Tip cm
G00650	TSF-35-145		N/A	.035	145	6
G00652	TSF-35-180		N/A	.035	180	6
G00655	TSF-35-260		N/A	.035	260	6
Heavy Duty with a Straight Tip PTFE-Coated Stainless Steel						
G01057	THSF-35-145		N/A	.035	145	7
G01237	THSF-35-180		N/A	.035	180	7
G00423	THSF-35-260		N/A	.035	260	7

Image Library



© Cook 2024

 +353 61 334 440

- Products
- Patient Resources
- Support
- About
- Divisions
- Contact
- Education & Training
- Customer Portal

- Privacy Statement
- Terms of Use
- Human Rights Responsibility Statement
- Cookie Compliance
- Sitemap

[Cookie Settings](#)

Sign Up. Get Email about News, Products, and Events from Cook Medical.

SIGN UP

Not all products shown on this website may be approved or available in all jurisdictions.
Consult with your local Cook representative, distribution company or customer support centre for details.