



## Medoject® Hypodermic needles

Colour Code	Size Metric	Size G <sup>17</sup>	Article No.
Light yellow	0.30 x 12 mm	30G x 1/2"	CH30012
Grey	0.40 x 12 mm	27G x 1/2"	CH27012
Grey	0.40 x 20 mm	27G x 4/5"	CH27045
Grey	0.40 x 20 mm	27G x 3/4"	CH27034
Grey	0.40 x 40 mm	27G x 1 1/2"	CH27112
Brown	0.45 x 12 mm	26G x 1/2"	CH26012
Brown	0.45 x 16 mm	26G x 5/8"	CH26058
Brown	0.45 x 25 mm	26G x 1"	CH26100
Orange	0.50 x 16 mm	25G x 5/8"	CH25058
Orange	0.50 x 20 mm	25G x 4/5"	CH25045
Orange	0.50 x 25 mm	25G x 1"	CH25100
Orange	0.50 x 40 mm	25G x 1 1/2"	CH25112
Purple	0.55 x 25 mm	24G x 1"	CH24100
Purple	0.55 x 38 mm	24G x 1 1/2"	CH24112
Blue	0.60 x 16 mm	23G x 5/8"	CH23058
Blue	0.60 x 25 mm	23G x 1"	CH23100
Blue	0.60 x 30 mm	23G x 1 1/4"	CH23114
Blue	0.60 x 40 mm	23G x 1 1/2"	CH23112
Blue	0.60 x 60 mm	23G x 2 1/4"	CH23214
Black	0.70 x 25 mm	22G x 1"	CH22100
Black	0.70 x 30 mm	22G x 1 1/4"	CH22114
Black	0.70 x 35 mm	22G x 1 2/5"	CH22125
Black	0.70 x 40 mm	22G x 1 1/2"	CH22112
Black	0.70 x 50 mm	22G x 2"	CH22200
Green	0.80 x 16 mm	21G x 5/8"	CH21058
Green	0.80 x 25 mm	21G x 1"	CH21100
Green	0.80 x 30 mm	21G x 1 1/4"	CH21114
Green	0.80 x 40 mm	21G x 1 1/2"	CH21112
Green	0.80 x 50 mm	21G x 2"	CH21200
Yellow	0.90 x 25 mm	20G x 1"	CH20100
Yellow	0.90 x 30 mm	20G x 1 1/4"	CH20114
Yellow	0.90 x 40 mm	20G x 1 1/2"	CH20112
Yellow	0.90 x 50 mm	20G x 2"	CH20200
Ivory	1.10 x 25 mm	19G x 1"	CH19100
Ivory	1.10 x 30 mm	19G x 1 1/4"	CH19114
Ivory	1.10 x 40 mm	19G x 1 1/2"	CH19112
Ivory	1.10 x 50 mm	19G x 2"	CH19200
Pink	1.20 x 25 mm	18G x 1"	CH18100
Pink	1.20 x 40 mm	18G x 1 1/2"	CH18112
Pink	1.20 x 50 mm	18G x 2"	CH18200
Red purple	1.50 x 50 mm	17G x 2"	CH17200
White	1.60 x 25 mm	16G x 1"	CH16100
White	1.60 x 40 mm	16G x 1 1/2"	CH16112
Blue grey	1.80 x 40 mm	15G x 1 1/2"	CH15112

## Medoject®

Nepirogeniškos

Netoksiškos

Medoject needles are made of polypropylene hub and stainless steel tube. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. The needles are intended for intramuscular, subcutaneous, intravenous and intra-arterial application. They are individually packed in boxes of 100 pcs, cartons of 4000 pcs.

**USER CONVENIENCE:** Colour-coded transparent hub for flow visual control, metric and gauge measurements on packaging and mainly thin wall tube material.

**USER SAFETY:** Safety needle is available.

**COMFORT:** Grinded, precision lancet point and silicone treated for smooth penetration.

**PATIENT SAFETY:** Single-use only. Sterilised by ETO.

Sterilizuotos Etileno oksido dujomis

## Medoject® Hypodermic needles

Medoject Hypodermic needles are intended for intramuscular, subcutaneous, intravenous and intra-arterial application. There is possibility of short bevel lancet point grinding.

According to EN, ISO standards.

Quality management certified according to EN ISO 13485:2016. CE marked.

Shelf life 5 years



### Medoject® Safety Hypodermic needles

Medoject Safety Hypodermic needles have an integrated safety feature for one-handed manipulation.

Colour Code	Size Metric	Size G <sup>1/2</sup>	Article No.
Blue ●	0.60 x 30 mm	23G x 1 1/4"	CHS23114
Black ●	0.70 x 40 mm	22G x 1 1/2"	CHS22112
Green ●	0.80 x 40 mm	21G x 1 1/2"	CHS21112
Yellow ●	0.90 x 40 mm	20G x 1 1/2"	CHS20112

**Shelf life 5 years**

### Medoject® Blunt fill needles

Medoject Blunt fill needles are also available with filter of 5µm for aspiration filtration procedure, which risks particle elimination during procedure.

Colour Code	Size Metric	Size G <sup>1/2</sup>	Article No.
Red ●	1.20 x 40 mm Blunt fill	18G x 1 1/2"	CH18112BF
Pink ●	1.20 x 40 mm Blunt fill w/ filter	18G x 1 1/2"	CH18112F
Pink ●	1.20 x 50 mm Blunt fill w/ filter	18G x 2"	CH18200F

**Shelf life 5 years**



According to EN, ISO standards.  
Quality management certified according to EN ISO 13485:2016. CE marked.

### Medoject® Special Hypodermic needles Specialios hipoderminės adatos P. d. 3.1

Medoject Special Hypodermic needles are intended for a wide special application range. There is possibility of short bevel lancet point grinding.

Colour Code	Dydis (mm) Dydis (G)		Article No.
	Size Metric	Size G <sup>1/2</sup>	
Blue ●	0.60 x 80 mm	23G x 3 1/8"	CH23318
Green ●	0.80 x 100 mm	21G x 4"	CH21400
Green ●	0.80 x 120 mm	21G x 4 3/4"	CH21434
Yellow ●	0.90 x 70 mm	20G x 2 4/5"	CH20245
Pink ●	1.20 x 70 mm	18G x 2 4/5"	CH18245
Pink ●	1.20 x 100 mm	18G x 4"	CH18400
White ○	1.60 x 70 mm	16G x 2 4/5"	CH16245
Blue grey ●	1.80 x 70 mm	15G x 2 4/5"	CH15245
Pale green ●	2.00 x 70 mm	14G x 2 4/5"	CH14245
Pale green ●	2.00 x 80 mm	14G x 3 1/8"	CH14318
Pale green ●	2.00 x 100 mm	14G x 4"	CH14400
Pale green ●	2.00 x 120 mm	14G x 4 3/4"	CH14434

Other sizes available upon request.

**Shelf life 5 years**

### Inox® Ophthalmic needles

Inox Ophthalmic needles has a blunt end for ophthalmic application, which includes scratch-free edges for comfort use..

Colour Code	Size Metric	Size G <sup>1/2</sup>	Article No. STRAIGHT	Article No. BENT
Grey ●	0.40 x 20 mm	27G x 4/5"	CHOP27045S	CHOP27045B
Blue ●	0.60 x 20 mm	23G x 4/5"	CHOP23045S	CHOP23045B

**Shelf life 5 years**

# The systems at a glance

Delta Cut biopsinė sistema (šaudyklė ir adatos)

## DeltaCut

Product	Art. No.	PU
DeltaCut biopsy system	304B000001	1



### Dydis


### Kodas

### Suitable coaxial cannula

Product	Size	Art. No.	PU
<b>DeltaCut biopsy cannulas, sterile</b>	18 G x 100 mm (1,20 mm)	315S120100	10
	18 G x 150 mm (1,20 mm)	315S120150	10
	18 G x 200 mm (1,20 mm)	315S120200	10
	18 G x 250 mm (1,20 mm)	315S120250	10
	18 G x 280 mm (1,20 mm)	315S120280	10
	16 G x 100 mm (1,60 mm)	315S160100	10
	16 G x 150 mm (1,60 mm)	315S160150	10
	16 G x 200 mm (1,60 mm)	315S160200	10
	14 G x 100 mm (2,10 mm)	315S210100	10
	14 G x 150 mm (2,10 mm)	315S210150	10

Size	Art. No.	PU
17 G x 44 mm	313S150014	25
17 G x 94 mm	313S150064	10
17 G x 144 mm	313S150114	10
14 G x 44 mm	313S200014	25
14 G x 94 mm	313S200064	10
13 G x 44 mm	313S250014	25

## CoreCut



Product	Size	Art. No.	PU
<b>CoreCut, sterile</b>	18 G x 150 mm (1,20 mm)	515S120150	1
	18 G x 200 mm (1,20 mm)	515S120200	1
	18 G x 250 mm (1,20 mm)	515S120250	1
	16 G x 150 mm (1,60 mm)	515S160150	1
	16 G x 200 mm (1,60 mm)	515S160200	1
	14 G x 100 mm (2,10 mm)	515S210100	1
	14 G x 150 mm (2,10 mm)	515S210150	1



### 3 Way Stop Cock



- 360° tap rotation.
- Red and blue pegs/stubs for arterial and venous line identification.
- Pressure tested to 5 bars.
- Minimum priming volume & low dead space.
- Polycarbonate material.
- Sterile, transparent, single use.
- Male/male/female connectors.
- It can be open three channels at the same time.



### 3-jų krypčių kranelis



- 360° kranelio sukimasis.
- Raudoni ir mėlyni kaiščiai / smeigtukai arterijų ir venų identifikavimui.
- Slėgis išbandytas iki 5 barų.
- Minimalus užpildymo tūris ir maža tuščioji erdvė.
- Polikarbonato medžiaga.
- Sterilus, skaidrus, vienkartinis.
- Vyriškos/vyriškos/moteriškos jungtys.
- Vienu metu galima atidaryti tris kanalus.

# SAFECATH

Disposable Central Venous Catheter Kit 

## Product Information


### Features

- Soft blue tip
- Patented connection design
- Flexible polyurethane material
- Latex-free

Be latekso

### Ranges

- Single/double/triple lumen
- 5-50cm multi-catheter lengths



# Package Componets



styga-pravedējas 0,032“ vienas galas tiesus-minkštas, kitas J formas

## Standard Package:

- Catheter (1)
- Guidewire (1)
- Dilator (1)
- Introducer Needle (0/1)
- Y-Introducer Needle (0/1)
- Catheter Clamp (1)
- Rigid Fastener (1)
- Blue Introducer Syringe (0/1)
- Heparin Cap (1/2/3)
- Stanch Clip (1/2/3)
- Syringe (1/2)
- Injection Needle (1)
- Scalpe (1)

## Ordering Information

Model	guidewire kink-proof nitinol guidewire, one end straight, the other “J” shaped;		dilator material: PP Plētiklis	blue introducer syringe	introducer needle		Y-introducer needle		fixatio n clip	number syringes (5ml)	Symbol on thr package
	OD/inch	length/mm			OD/G	length/mm	OD/G	length/mm			
MMCVCBJ 1-14-15	0.035	600	8F	5 mL	18	70	/	/	Ø2.1	1	
MMCVCBJ 1-14-20	0.035	600	8F	5 mL	18	70	/	/	Ø2.1	1	
MMCVCBJ 1-14-30	0.035	700	8F	5 mL	18	70	/	/	Ø2.1	1	
MMCVCBJ 1-16-15	0.035	600	7F	5 mL	18	70	/	/	Ø1.5	1	
<b>MMCVCBJ 1-16-20</b>	<b>0.035</b>	<b>600</b>	<b>8F</b>	<b>5 mL</b>	<b>18</b>	<b>70</b>	/	/	Ø1.5	1	
MMCVCBJ 1-16-30	0.035	700	7F	5 mL	18	70	/	/	Ø1.5	1	
MMCVCBJ2-70-15	0.035	600	8.5F	5 mL	18	70	/	/	Ø2.3	1	
MMCVCBJ2-70-20	0.035	600	8.5F	5 mL	18	70	/	/	Ø2.3	1	
MMCVCBJ2-70-30	0.035	700	8.5F	5 mL	18	70	/	/	Ø2.3	1	
MMCVCBJ2-70-50	0.035	700	8.5F	5 mL	18	70	/	/	Ø2.3	1	
MMCVCBJ3-70-15	0.035	600	8.5F	5 mL	18	70	/	/	Ø2.3	1	
MMCVCBJ3-70-20	0.035	600	8.5F	5 mL	18	70	/	/	Ø2.3	1	
MMCVCBJ3-70-30	0.035	700	8.5F	5 mL	18	70	/	/	Ø2.3	1	
MMCVCBJY1-14-15	0.035	600	8F	/	/	/	18	70	Ø2.1	2	
MMCVCBJY1-14-20	0.035	600	8F	/	/	/	18	70	Ø2.1	2	
MMCVCBJY1-14-30	0.035	700	8F	/	/	/	18	70	Ø2.1	2	
MMCVCBJY1-16-15	0.035	600	7F	/	/	/	18	70	Ø1.5	2	
MMCVCBJY1-16-20	0.035	600	7F	/	/	/	18	70	Ø1.5	2	
MMCVCBJY1-16-30	0.035	700	7F	/	/	/	18	70	Ø1.5	2	
MMCVCBJY1-18-15	0.021	600	5F	/	/	/	20	38	Ø1.3	2	
MMCVCBJY1-18-20	0.021	600	5F	/	/	/	20	38	Ø1.3	2	
MMCVCBJY1-18-30	0.021	600	5F	/	/	/	20	38	Ø1.3	2	
MMCVCBJY1-20-13	0.021	600	5F	/	/	/	20	38	Ø1.0	2	
MMCVCBJY1-20-20	0.021	600	5F	/	/	/	20	38	Ø1.0	2	
MMCVCBJY2-40-05	0.021	600	5F	/	/	/	20	38	Ø1.3	2	
MMCVCBJY2-40-08	0.021	600	5F	/	/	/	20	38	Ø1.3	2	
MMCVCBJY2-40-13	0.021	600	5F	/	/	/	20	38	Ø1.3	2	
MMCVCBJY2-50-08	0.021	600	6F	/	/	/	20	38	Ø1.7	2	
MMCVCBJY2-50-13	0.021	600	6F	/	/	/	20	38	Ø1.7	2	
MMCVCBJY2-50-20	0.021	600	6F	/	/	/	20	38	Ø1.7	2	
MMCVCBJY2-70-15	0.035	600	8.5F	/	/	/	18	70	Ø2.3	2	
MMCVCBJY2-70-20	0.035	600	8.5F	/	/	/	18	70	Ø2.3	2	
MMCVCBJY2-70-30	0.035	700	8.5F	/	/	/	18	70	Ø2.3	2	
MMCVCBJY2-70-50	0.035	700	8.5F	/	/	/	18	70	Ø2.3	2	
MMCVCBJY3-55-08	0.021	600	6F	/	/	/	20	38	Ø1.7	2	
MMCVCBJY3-55-13	0.021	600	6F	/	/	/	20	38	Ø1.7	2	
MMCVCBJY3-70-15	0.035	600	8.5F	/	/	/	18	70	Ø2.3	2	
MMCVCBJY3-70-20	0.035	600	8.5F	/	/	/	18	70	Ø2.3	2	
MMCVCBJY3-70-30	0.035	700	8.5F	/	/	/	18	70	Ø2.3	2	


8 p.d.

# PRODUCT TECHNICAL DATA SHEET

## DESCRIPTION

**TORASET – Thoracentesis drainage set with 14G-16G-19G needles**

**Pleuros punkcijos rinkinys su trimis adatomis**



2l surinkimo maišas

CODE	BAG VOLUME ml	NEEDLES trys adatos, 80mm	CARTON pcs.	GTIN-UDI SINGLE PRODUCT	GTIN-UDI CARTON
M032010S	2000	14G/16G/19G – 80 mm	30	18053676294524	38053676294528
M032000S	2000	14G/16G/19G – 55 mm	30	18053676294494	38053676294498

## ● PRODUCT COMPOSITION:

The device consists of:

1. **2000 ml collection bag** in PVC, with non-return valve, "push-pull" drain tap, **sujungimo vamzdelis** connection tube I.D. 5.2 mm x O.D. 6.5 mm and **length cm 90: ilgis- 90cm**
2. **Male Luer Lock connector** in ABS;
3. **3-way stopcock** (360° rotation) made in **Polyethylene** and **Polycarbonate**, with fluid direction indicators to identify the different procedure phases (aspiration of pleural effusion from the patient; discharge into collection bag);  
**3-ju krypčių kranelis**
4. **Protective cap** in blue **ABS**;
5. **Female Luer lock connector** in **ABS**;
6. **Tube** in **DEHP FREE PVC**, I.D. 4.8 mm x O.D. 6.8 mm and length cm 40;  
**prailginimo linija**
7. **Vented cap** in **Polyethylene**;
8. **60 ml syringe** to be connected to the 3-way stopcock, made in **Polypropylene** and **Synthetic Isoprene rubber**;  
**luer lock švirkštas 60ml**
9. **Stainless Steel 14G Needle** with a length of 55 or 80 mm (see table), packaged in microperforated bag;
10. **Stainless Steel 16G Needle** with a length of 55 or 80 mm (see table), packaged in microperforated bag;
11. **Stainless Steel 19G Needle** with a length of 55 or 80 mm (see table), packaged in microperforated bag.

## ● INTENDED USE:

Closed circuit system thoracentesis, particularly suitable for the aspiration of pleural effusion liquid.

Following an evaluation by the healthcare professional, and under his direct control, the device can be used for paracentesis procedures.

## ● PACKAGING: **steriliame įpakavime**

Primary packaging: single device in micro-perforated **PE** inner pouch bag + blister in medical grade paper/PP-PE film.

Secondary packaging: carton with 30 pieces.

## ● PRODUCTION PROCESS:

The device is manufactured in accordance with HMC Premedical S.p.A. Quality System and in comply with the requirements of the standard EN ISO 13485.

## ● CONTROL ON THE PRODUCT:

In all stages of processing, according to internal procedures and to sampling plans defined by norm ISO 2859-1.

## ● CLASSIFICATION:

Class IIa sterile.

## ● STERILIZATION: **sterilizuota EO**

Product sterilized by ethylene oxide according to a validated course of treatment in accordance with current regulations.

Shelf life: **5 years** from the date of sterilization.

**For single use** only and non re-sterilizable.

**Vienkartiniam naudojimui**

## ● STORAGE:

Standard storage procedures and conditions.

## ● DISPOSAL:

The user must follow the legal regulations regarding disposal of hospital waste.

## ● WARNINGS:

The device must be used exclusively by healthcare professionals.

## ● REGISTRATION TO ITALIAN MEDICAL DEVICES REPERTOIRE:

CND: **A060204**

RDM: **1575159**

GMDN: **36787**

## ● UNIT OF SALE:

Carton with **30 pieces**.

## ● GTIN-UDI:

See table.



HMC PREMEDICAL S.p.A.  
Via Bosco, 1/3 - 41037 Mirandola (MO) Italy  
TEL. +39 0535 22704 - FAX. +39 0535 609546  
[info@hmcgroup.it](mailto:info@hmcgroup.it)

MOD. 07 POS. 002 REV.00  
DATA 23/06/2021

# PRODUCT TECHNICAL DATA SHEET



REV.	CHANGES	ISSUED	VERIFIED AND APPROVED	DATE
00	First issue – english language	E. Benassi	D. Bosetti	24/10/2017
01	Data sheet revision; drawing and product composition update	E. Benassi	S. Tralli	31/01/2020
02	Intended use updated	E. Benassi	S. Tralli	29/09/2020
03	Leggl manufacturer's address update	E. Benassi	S. Tralli	28/04/2021
08	Data sheet update; jump to rev. 08 to align with italian version	E. Benassi	S. Tralli	09/09/2021
09	Table update with GTIN-UDI	E. Benassi	S. Tralli	02/02/2022
10	Drawing update	E. Benassi	S. Tralli	11/02/2022

## BPDS-2000 Pleural Drainage System

- Underwater seal drainage
- 2000 ml collecting capacity
- Firm stand to avoid accidental overturn
- Handle for patient convenience
- Easy to write on drainage marking scale
- Funnel to facilitate priming
- Easy to milk, soft and kink resistant silk finished tubing
- Corrugated support tubing that prevents drainage obstruction due to kinking at collector inlet
- Length of tubing: 150 cm
- Tube connector specially designed to accept different thoracic catheter sizes
- Optional luer lock port cap on the connector helps to empty coagulum collected in the tubing
- Priming volume: 180cc




Ref	
175 2001 1	w/ Luer Lock Port
175 2002 1	w/ Antifoam
175 2005 1	w/ Latex Tubing
017 2004 1	w/o Luer Lock Port

## BPDS-700 Pleural Drainage System

Pleuros drenažo sistema

- Underwater seal drainage
- 700 ml collecting capacity
- 700 ml graduota plastikinė surinkimo talpa
- Firm stand to avoid accidental overturn
- Handle for patient convenience
- Easy to write on drainage marking scale
- Cap
- Easy to milk, soft and kink resistant silk finished tubing
- Corrugated support tubing that prevents drainage obstruction due to kinking at collector inlet
- Length of tubing: 150 cm
- 150cm vamzdelis
- Tube connector specially designed to accept different thoracic catheter sizes
- Optional luer lock port cap on the connector helps to empty coagulum collected in the tubing
- Priming volume: 50cc



vamzdelio galas  
pritaikytas  
skirtingiems  
krūtinės kateterių  
dydžiams

vienetiniame  
steriliame  
įpakavime

Ref	
175 7001 1	w/ Luer Lock Port
175 7002 1	w/ Antifoam
017 7004 1	w/o Luer Lock Port

## Ultra Cath Novo

### Thoracic Drainage Catheters

90080-90089

#### Ultra Cath Novo - Straight

#### Thoracic Drainage Catheter - Straight (DEHP free)

- Suitable for effective drainage after cardio-thoracic & thoracic surgery
- Atraumatic & rounded open distal end with smooth eyes
- Proximal end is fitted with pull through tapered tongue
- With radio-opaque line and markings at every 2 cm from the last eye
- Cross side eyes to prevent tissue aspiration
- Matching size connector for easy connection to the drainage system
- Soft, frosted and kink resistant PVC tubing
- Length: 45 cm
- Sizes: 8, 10, 12, 14, 16, 18, 20, 22, 24, 28, 32, 36 & 40 FG



90540-90549

#### Ultra Cath Novo - Angled

#### Thoracic Drainage Catheter - Angled (DEHP free)

- Catheter is provided with curved (90°) angle for easy aspiration of drainage

90090-90178 P.d. 17.1

#### Ultra Cath Novo - Trocar

#### Thoracic Drainage Catheter with Trocar (DEHP free)

- Catheter is provided with trocar for easy penetration, placement and handling
- Markings provided at 5, 10, 15 & 20 cm from the last eye to help ascertain depth of placement
- Two lateral eyes for efficient drainage
- Proximal end fitted with tapered connector for easy connection to the drainage bottle
- Size : 8, 10, 12, 14, 16, 18, 20, 22, 24, 28, 32, 36, 40 FG



90090 - 90178

**TORAKALINIS KATETERIS SU TROKARU** (be DEHP)

- Kateteris pateikiamas kartu su trokaru lengvam įsiskverbimui, įtaisymui ir valdymui
- Žymos kas 5, 10, 15 ir 20 cm nuo paskutinės akies, padedančios nustatyti įtaisymo gylį
- Dvi šoninės akys efektyviam drenažui
- Proksimalinis galiukas pritaikytas siaurėjančiam konektoriui lengvam susijungimui su drenažo buteliu
- Dydžiai: 8, 10, 12, 14, 16, 18, 20, 22, 24, 28, 32, 36, 40 FG

# Vaginal Speculum P. d. 20.1 Vaginální spekule

**Description:**

It can be used for the checking and treatment of gynecology disease in hospitals or clinics. Before using it, user must read the usage manual, and then take the speculum, push the "duckbill" of speculum that is closed into the vaginal slowly. Open the "duckbill" of speculum according to demand, fix the nut of speculum, then it can be used for checking and treatment of gynecology disease. After all, loose the nut of speculum, and take "duckbill" of speculum out of the vaginal slowly.

**Specification:**

Made of polystyrene;  
It is a transparent substance, non-toxic, medical grade  
Size: Large/Middle/Small  
Supplied sterile in individual polybag or Film+ Polybag or blister pack  
sterilized by ethylene oxide gas

**From of packing:**

sterilized in individual polybag or Film+ Polybag or blister pack; 100pcs/ carton;

**Size of packing:**

Export carton:520x380x360mm

**G.Weight.:7kgs**

**N.Weight: 6kgs**

Other Medical Products



**GT159-100** Disposable Vaginal Speculum with side screw




**GT159-103** Disposable Vaginal Speculum with middle screw



**GT159-106** American type Disposable Vaginal Speculum



**GT159-101** French type 1 Disposable Vaginal Speculum  
Vienkartinė vaginální spekule



**GT159-104** Australia type Disposable Vaginal Speculum



**GT159-102** Spanish type Disposable Vaginal Speculum



**GT159-105** South Asian type Disposable Vaginal Speculum



**GT207-P22 Cervical Rambrush**

Usage: Sampling Examination For Male And Female


Material: Stainless Steel Thread + PVC Handle + Nylon Tip, PVC

Model: Slid Type, Injection Type

Packing: 1Pc/Pack, 100Packs/box,


50box/ 41 X 35 X 14cm, 13/12Kgs

# ECG Paper




- Single-lead ECG Paper
- 3-lead ECG Paper
- 6-lead ECG Paper
- 12-lead ECG Paper

**GT148-101**  
ECG Paper




50mmX20m-30m


Single-lead ECG Paper



6-lead ECG Paper



3-lead ECG Paper




12-lead ECG Paper

- Single-lead ECG Paper: 50mm X 20m, 200pcs/ctn
- 3-lead ECG Paper: 90mm X 90mm-200P, 100pcs/ctn
- 6-lead ECG Paper: 112mm X 90mm-200P, 100pcs/ctn
- 12-lead ECG Paper: 210mm X 300mm-200P, 20pcs/ctn

Other Medical Products

1. Sketch map of the structure of product:



2.1 Vaginal Speculum French Type



**Basic Size:**

Dydis Size	Basic Dimensions							
	L	L1	L2	B	B1	C	C1	H
L	109.5±1	109±1	116.5±2	9.5±0.5	9.5±0.5	30±0.5	30±5	83±3
M	108.5±1	108±1	116.5±2	9.2±0.5	9.2±0.5	26.5±0.5	26.5±0.5	82±3
S	108.5±1	108±1	116.5±2	9.2±0.5	9.2±0.5	24.2±0.5	24.2±0.5	81.5±3
XS	108.5±1	108±1	116.5±2	9±0.5	9±0.5	20.5±0.5	20.5±0.5	81±3
XXS	108.5±1	108±1	116.5±2	7.7±0.5	7.7±0.5	15±0.5	15±0.5	81±3

## Rectal Tube


- (1) Made of non-toxic, non-irritant PVC (Medical grade) or DEHP free.
- (2) Color: Transparent or translucent (Surface fog form).
- (3) Closed tip with two lateral eyes.
- (4) Connector (color coded)-Size (Fr/Ch): F18, F20, F22, F24, F26, F28, F30, F32, F34, F36.
- (5) Package: Single in PE or blister package.
- (6) Sterilized by EO GAS.



## GT020-100

### Rektalinis zondas

- (1) Pagamintas iš netoksiško, nedirginančio PVC (medicinės paskirties) ir be DEHP.
- (2) Spalva: skaidrus arba permatomas.
- (3) Uždaras galiukas su dviem šoninėm angelėms.
- (4) Konektorius (koduotas spalva)- **Dydžiai (Fr/Ch): F18, F20, F22, F24, F26, F28, F30, F32, F34, F36.**
- (5) Pakuotė: Individualus įpakavimas PE pakuotėje arba Blister pakuotėje
- (6) **Sterilizuota EO duiomis.**







## GENERAL SURGERY

WWW.ZARYS.COM

atsiurbimo sujungimo vamzdelis

### surgical suction tube

sterile

#### features:

- made of soft and flexible PVC of optimal thickness to prevent the tube from collapsing and facilitate operating at the same time
- available types of connectors:
  - funnel - funnel
  - funnel - kapkon
  - funnel-funnel cut-to-fit
- available sizes:
  - CH 24 (outer diameter: 8.0 mm, inner diameter: 5.6 mm)
  - CH 30 (outer diameter: 10.0 mm, inner diameter: 7.1 mm)
- grooves on the surface of the tube prevent it from bending, twisting and collapsing
- smooth inner surface prevents secretions and body fluids from sticking to the walls of the tube and ensures correct and optimal flow
- all types of connectors are compatible with standard suction cannulas
- single-use **sterilus**
- EO sterilized **vienkartinis**
- double packaging: internal- foil, external- paper-foil

#### indication:


- to evacuate secretions and body fluids



ref.	dydis	connector	vamzdelio ilgis	bulk packaging
	size		tube length	
P. d. 24.1 DO-24201	CH 24	funnel-funnel	210 cm	50 pcs
DO-24301	CH 24	funnel-funnel	300 cm	50 pcs
P. d. 24.2 DO-30201	CH 30	funnel-funnel	210 cm	50 pcs
DO-30301	CH 30	funnel-funnel	300 cm	50 pcs
DO-24202	CH 24	funnel-kapkon	210 cm	50 pcs
DO-24204	CH 24	funnel-funnel cut-to-fit	210 cm	50 pcs

## Silicone drain sterile

Article number: 21xxx | Packing unit (VPE): 25 | Shipping unit (VSE): 25




- Flexible transparent silicone drain;
- Open atraumatic nozzle with 6 side holes;
- With contrast strips;
- Individually packaged and sterilized;
- Length: 100-500 mm

CH	PU/SU	REF
8	25/100	21899
10	25/100	21900
12	25/100	21901
14	25/100	21130
15	25/100	21902
18	25/100	21903
20	25/100	21131
21	25/100	21904
24	25/100	21905
26	25/100	21132
27	25/100	21906
30	25/100	21907
33	25/100	21908
36	25/100	21909
39	25/100	21133



## Silikoninis drenas

Article number: 21xxx | Packing unit (VPE): 25 | Shipping unit (VSE): 25



- Lankstus skaidraus silikono drenas;
- Atviras atraumatinis antgalis su 6 šoninėmis skylutėmis;
- Su kontrastinėmis juostelėmis;
- Individualiai supakuotas ir sterilizuotas;
- Visas ilgis: 500 mm

CH	Ilgis	PU/SU	REF
8	500 mm	25/100	21899
10	500 mm	25/100	21900
12	500 mm	25/100	21901
14	500 mm	25/100	21130
15	500 mm	25/100	21902
18	500 mm	25/100	21903
20	500 mm	25/100	21131
21	500 mm	25/100	21904
24	500 mm	25/100	21905
26	500 mm	25/100	21132
27	500 mm	25/100	21906
30	500 mm	25/100	21907
33	500 mm	25/100	21908
36	500 mm	25/100	21909
39	500 mm	25/100	21133

\ Treatment \ Post-operative Redon systems \ Low vacuum system

**500ml**



**low vacuum system for routine use in operating theatre**  
 product features:

- starting with a vacuum of 120 mbar
- content: 500 ml
- connecting tube, length: 125 cm
- universal stepped connector for safe connection of ch. 6 to ch. 18 redon drains
- set with 2 drains and guiding needle
- length of the drain: 800 mm
- perforation of the drain: 150 mm
- individually packed and sterilized
- special designs available on request
- UMDNS-Code: 11-307
- Klasse/Regel: IIa/7

**Art.-Nr.: 21586, PRI-LOW-VAC Redon set, 500 ml with drain (500/150) Ch. 16** P.d. 29.1

Packing unit (VPE): 35, Shipping unit (VSE): 35



**low vacuum system for routine use in operating theatre**  
 product features:

- starting with a vacuum of 120 mbar
- content: 500 ml
- connecting tube, length: 125 cm
- **universal stepped connector for safe connection of ch. 6 to ch. 18 redon drains**
- length of the drain: 500 mm
- perforation of the drain: 150 mm
- sets with drain
- individually packed and sterilized
- special designs available on request

UMDNS-Code: 11-307


Klasse/Regel: IIa/7

**Art.-Nr.: 21587, PRI-LOW-VAC Redon set, 500 ml with drain (500/150) Ch. 18**

Packing unit (VPE): 35, Shipping unit (VSE): 35

## Redon drains (PVC) without trocar Redon drenai

Redon drains are suitable for use with high-vacuum systems and low-vacuum systems for postoperative wound drainage.



### Product features:

- Made of soft, flexible polyvinyl chloride (PVC)
- With radiopaque stripe
- Implant tested in accordance with USP XXIII
- Alternating perforation to avoid ingrowth of tissue
- Rounded drainage openings for removal of secretion
- Length markings allow simple repositioning

Ch.	Length	PU/SU	REF
6	500 mm	100/300	21860
8	500 mm	100/300	21861
10	500 mm	100/300	21862
12	500 mm	100/300	21863
14	500 mm	100/300	21864
16	500 mm	100/300	21865
18	500 mm	100/300	21866

Ch.	Length	PU/SU	REF
6	800 mm	100/300	21880
8	800 mm	100/300	21881
10	800 mm	100/300	21882
12	800 mm	100/300	21883
14	800 mm	100/300	21884
16	800 mm	100/300	21885
18	800 mm	100/300	21886

## Žemo vakuumo Redon, operacinio lauko sistema su universalium konektoriumi, CH 6-18



### Žemo vakuumo sistema kasdieniam naudojimui operacinėse

Produkto savybės:

- ▲ pradedant 120 mbar vakuumu;
- ▲ talpa: 20, 50, 200, 250, 500 ml;
- ▲ jungiamasis vamzdelis, ilgis: 125 cm;
- ▲ universalus pakopinis konektorius saugiam sujungimui, CH 6-18;
- ▲ rinkinys su dreneu ir pagrindine adata;
- ▲ atskirai supakuoti, sterilūs;
- ▲ galimi specialūs modeliai;

# LOW-VACUUM SYSTEMS

## PRI-LOW-VAC® 500 ml



### PRI-LOW-VAC® 500 ml complete set

- Redon bellows with connecting tube, Redon drain (PVC) with trocar

Ch.	Drain length	PU/SU	REF
10	800 mm	35/ 35	<b>21605</b>
12	800 mm	35/ 35	<b>21601</b>
14	800 mm	35/ 35	<b>21599</b>



### Talpa 500ml

### PRI-LOW-VAC® 500 ml OR system

- Redon bellows with connecting tube

PU/SU	REF
35/ 35	<b>21631</b>



### PRI-LOW-VAC® 500 ml ward system (replacement bellows)

- Redon bellows with closed lid attachment

PU/SU	REF
35/ 35	<b>21621</b>

**KARMAN ŞİRINGA Karman Syringe**

Tek ve çift valfli  
Manuel vakumlu  
5 ml Medikal grade silikon  
Yedek malzemeler  
Reusable  
4, 5, 6, 7 kanüller tek valfli  
8, 10, 12, kanüller çift valfli  
şırınga ile kullanılır.

Single and double valve  
Manual vacuum  
5 ml Medical grade silicon  
Spare parts  
Reusable  
4, 5, 6, 7 karmans - can be use with single  
valve syringe  
8, 10, 12 karmans - can be use with double  
valve syringe



Tip/Type	Ref No
Tek Valfli/Single Valve	240 102
Çift Valfli/Double Valve	240 101

**P.d. 31.1 KARMAN KANUL Karman Canula**

Tip/Type	Ref No
No : 4	240 111
No : 5	240 112
No : 6	240 113
No : 7	240 114
No : 8	240 115
No : 9	240 116
No : 10	240 117
No : 11	240 119
No : 12	240 118

**ENDOMETRİAL ÖRNEK ALMA KANULÜ  
Endometrial Cell Sampler Cannula**

Yuvarlatılmış uç  
Hassas derecelendirilmiş silindir  
Sızdırmaz steril conta  
3,1 mm dış çap

Rounded tip  
Precisely graduated cylinder  
Sterile  
3,1 mm outside diameter



Tip/Type	Ref No
Standart Steril/Sterile	240 201



## KARMAN Kaniulė



vienkartinis,  
sterilus

Viengubo arba dvigubo vožtuvo variantai

Rankinis vakuumas

Medicininės paskirties silikonas 5ml

Atsarginės dalys

Dydis	Katalogo nr.
Nr:4	240 111
Nr:5	240 112
Nr:6	240 113
Nr:7	240 114
Nr:8	240 115
Nr:9	240 116
Nr:10	240 117
Nr:11	240 119
Nr:12	240 118

	<b>TECHNICAL FILE</b>	<b>Document No</b>	KYS.TD.023/4
	<b>KARMAN SYRINGE &amp; KARMAN CANNULA</b>	<b>Revision No</b>	02
	Ref.MDD 93/42/EEC (2007/47/EC)	<b>Revision Date</b>	31.01.2011
	Product Ref. No: 240XXX	<b>Issue Date</b>	15.03.2005
		<b>Page Number</b>	7/8

#### 1.4.2. **KARMAN CANNULA** Karman kaniülë


<b>THE STANDARD SPECIFICATION of THE PRODUCT –TS 6074:1988</b>									
<i>PRODUCT SPECIFICATIONS</i>			<i>VALUE</i>						TOLERANCE
VISUAL SPECIFICATION and STRUCTURAL SECURITY			External and external surfaces must be smooth. It must be equally circular cross-section. Distal tip must be closed and atraumatic and have two lateral eyes for no: 4, 5, 6, 7, 8-single eye for no: 9,10,12						-
THE TIGHTNESS of THE CONNECTION BETWEEN ADAPTOR and CANNULA for 4 and 5.			Adaptor has length in 4mm and external diameter must be 6mm. Cannula must be inserted to the adaptor as not get out.						
THE TIGHTNESS of THE CONNECTION BETWEEN ADAPTOR and CANNULA for 6, 7, 8, 9 and 10.			Adaptor must be harmonies with the tip in 12mm of the syringe with double valve. Connection between adaptor and cannula must be tight.						
SCALA LINES of THE CANNULA			There must be at least 5 pieces centimeter signs that must to begin, where beyond 6 cm from the distal end of the cannula.						
THE TIGHTNESS of THE CONNECTION BETWEEN KARMAN CANNULA and CARMAN SYRINGE			Leakage must not be at the connection points between and karman cannula and carman syringe. It must be controlled the suction power of the syringe.						-
THE STRENGTH of THE NEGATIVE PRESSURE of THE CARMAN CANNULA			It must not be form deformation in shape of the carman cannula.						-
STERILITY TEST			<b>Privalo büti sterilus</b> <b>IT MUST BE STERILE</b>						-
PYROGENITY TEST			IT MUST BE APYROGEN						-
<b>DIMENSIONAL SPECIFICATION of THE CARMAN CANNULA</b>									
<b>CANNULA NO:</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>12</b>	TOLERANCE
NOMINAL EXT.DIAMETER	4	5	6	7	8	9	10	12	±%2.5mm
MINIMUM INT.DIAMETER	3,0	3,5	4,5	5,5	6,5	7,5	8,5	10,5	±%2.5mm
THE COLORS of THE ADAPTOR	SARI	YEŞİL	MAVİ	AÇIK KAHVE RENGİ	TURUNCU	KAHVE RENGİ	KOYU YEŞİL	LACIVERT	
THE LENGTH of THE LATERAL EYE	6,0	7,5	9,0	10,0	12,0	22,0	23,0	26,0	±%2.5mm
RADIUS of THE LATERAL EYE	1,80	2,25	2,70	3,10	3,60	4,5	5,0	6,0	±%2.5mm
DIST.BETWEEN TWO LAT.EYES	2,0	2,5	3,0	3,5	4,0	-	-	-	±%2.5mm
DIST.BETWEEN DIST.END and LATERAL EYE	4	4	4	4	4,5	5,5	5,5	5,5	±%2.5mm
DIST.BETWEEN DISTAL END and SCALA LINES	60	60	60	60	60	60	60	60	±%2.5mm
THE LENGTH of THE CANNULA	200	200	200	200	200	230	230	230	±%2.5mm

PREPARED BY: Coskun DURMAYUKSEL

APPROVED BY: Muammer Berksöz


**SURGICAL DRAPES**  
**2-LAYER DRAPES**
**STERILE DRAPES WITH CELLULOSE LAYER - 2-LAYERED**


ref.	dimensions	packaging unit	 bulk packaging
AT-PF-S 1	50cm x 50cm	1 pc.	700 pcs
AT-PF-S 3	50cm x 70cm	1 pc.	450 pcs
AT-PF-S 5	80cm x 80cm	1 pc.	200 pcs
AT-PF-S 6	75cm x 90cm	1 pc.	250 pcs
AT-PF-S 8	90cm x 150cm	1 pc.	125 pcs



135

**STERILE SURGICAL DRAPES - 2-LAYERED**


ref.	dimensions	packaging unit	 bulk packaging
AT-NF-S 21 E	37cm x 45cm	100 pcs	4 x 100 pcs
AT-NF-S 1 E	45cm x 45cm	100 pcs	4 x 100 pcs
AT-NF-S 3 E	45cm x 60cm	100 pcs	4 x 100 pcs
AT-NF-S 4 E	45cm x 75cm	100 pcs	4 x 100 pcs
AT-NF-S 5 E	75cm x 75cm	50 pcs	4 x 50 pcs
AT-NF-S 9 E	75cm x 90cm	50 pcs	4 x 50 pcs
AT-NF-S 13 E	100cm x 150cm	50 pcs	4 x 50 pcs
AT-NF-S 15 E	150cm x 180cm	20 pcs	4 x 20 pcs



economic  
version

**STERILE SURGICAL DRAPES - 2-LAYERED**

ref.	dimensions	packaging unit	 bulk packaging
AT-NF-S 1	45cm x 45cm	1 pc.	240 pcs
AT-NF-S 2	50cm x 50cm	1 pc.	200 pcs
AT-NF-S 3	50cm x 60cm	1 pc.	280 pcs
AT-NF-S 4	50cm x 70cm	1 pc.	160 pcs
AT-NF-S 5	75cm x 75cm	1 pc.	150 pcs
AT-NF-S 9	75cm x 90cm	1 pc.	150 pcs
AT-NF-S 10	90cm x 120cm	1 pc.	80 pcs
AT-NF-S 13	100cm x 150cm	1 pc.	50 pcs
AT-NF-S 14	120cm x 150cm	1 pc.	50 pcs
AT-NF-S 15	150cm x 180cm	1 pc.	20 pcs
AT-NF-S 16	150cm x 200cm	1 pc.	20 pcs
AT-NF-S 17	180cm x 180cm	1 pc.	20 pcs
AT-NF-S 20	150cm x 240cm	1 pc.	20 pcs
AT-NF-S 22	170cm x 300cm	1 pc.	20 pcs


**STERILE SURGICAL DRAPES WITH ADHESIVE EDGE - 2-LAYERED**

ref.	dimensions	packaging unit	 bulk packaging
AT-NFA-S 1 E	45cm x 45cm	100 pcs	4 x 100 pcs
AT-NFA-S 4 E	45cm x 75cm	100 pcs	4 x 100 pcs
AT-NFA-S 9 E	75cm x 90cm	50 pcs	4 x 50 pcs
AT-NFA-S 10 E	90cm x 120cm	50 pcs	4 x 50 pcs
AT-NFA-S 13 E	90cm x 150cm	50 pcs	4 x 50 pcs
AT-NFA-S 15 E	150cm x 175cm	20 pcs	4 x 20 pcs
AT-NFA-S 16 E	150cm x 200cm	20 pcs	4 x 20 pcs
AT-NFA-S 20 E	150cm x 240cm	20 pcs	4 x 20 pcs



adhesive tape  
on the longer  
side of drape



economic  
version

	<b>LIMB COVER</b> <b>sterile</b>	att.no.3 pt.2 subpt.2.3
		revision 2
	<b>TD-11</b>	date: 2017-10-19
		amendmend of: 2022-09-08

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02, fax. +48 32 370 38 08

#### INDICATION:

- Patient and equipment drapes for surgical procedures
- Intended for use in the operating room
- Intended for use during surgical procedures
- As an impermeable barrier for liquids and microorganisms
- Absorb excess fluids during surgical procedures
- Elimination of dusting due to the use of cotton-free materials

#### FEATURES:

Kojos dangalas pagamintas iš dv sluoksniu laminuoto audinio, chirurginė lipni juosta iš neaustinės medžiagos

- **Limb cover made of two-layer laminate, surgical adhesive tape made of spunlace non-woven**
- **Fully absorbent and impermeable to fluids** Visiškai absorbuojanti ir nepralaidi skysčiams
- Providing an optimal barrier against microbial migration
- Fulfill the requirements of the current standard EN 13795-1
- Sterile, single-use
- Packed with one piece of surgical adhesive tape
- Non-flammable
- Dust-free
- Latex-free
- **Paper-foil packaging** Popieriaus folijos pakuotė. Yra ir kartoninė skirta prekių gabenimui į operacinę, tretinė - skirta transportavimui.
- Label with four adhesive tag labels for data archiving. The following information is available on each adhesive tag label: ref number, expiry date, batch number, manufacturer's data and barcode



10x50cm

#### SIZE AND REFERENCE

CATALOGUE NUMBER	DIMENSIONS	INDIVIDUAL PACKAGING	BULK PACKAGING
AT-SD-S 6	30cm x 60cm	1 pc.+ adhesive tape 10x50cm	200 sets
AT-SD-S 7	30cm x 80cm	1 pc.+ adhesive tape 10x50cm	120 sets
<b>AT-SD-S 8</b>	<b>35cm x 120cm</b>	<b>1 pc.+ adhesive tape 10x50cm</b>	<b>100 sets</b>

	<b>LIMB COVER</b> <b>sterile</b>	att.no.3 pt.2 subpt.2.3
		revision 2
	<b>TD-11</b>	date: 2017-10-19
		amendmend of: 2022-09-08

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02, fax. +48 32 370 38 08

#### INDICATION:

- Patient and equipment drapes for surgical procedures
- Intended for use in the operating room
- Intended for use during surgical procedures
- As an impermeable barrier for liquids and microorganisms
- Absorb excess fluids during surgical procedures
- Elimination of dusting due to the use of cotton-free materials

#### FEATURES:

Kojos dangalas pagamintas iš dv sluoksniu laminuoto audinio, chirurginė lipni juosta iš neaustinės medžiagos

- **Limb cover made of two-layer laminate, surgical adhesive tape made of spunlace non-woven**
- **Fully absorbent and impermeable to fluids** Visiškai absorbuojanti ir nepralaidi skysčiams
- Providing an optimal barrier against microbial migration
- Fulfill the requirements of the current standard EN 13795-1
- Sterile, single-use
- Packed with one piece of surgical adhesive tape
- Non-flammable
- Dust-free
- Latex-free
- **Paper-foil packaging** Popieriaus folijos pakuotė. Yra ir kartoninė skirta prekių gabenimui į operacinę, tretinė - skirta transportavimui.
- Label with four adhesive tag labels for data archiving. The following information is available on each adhesive tag label: ref number, expiry date, batch number, manufacturer's data and barcode



10x50cm

#### SIZE AND REFERENCE

CATALOGUE NUMBER	DIMENSIONS	INDIVIDUAL PACKAGING	BULK PACKAGING
AT-SD-S 6	30cm x 60cm	1 pc.+ adhesive tape 10x50cm	200 sets
<b>AT-SD-S 7</b>	<b>30cm x 80cm</b>	<b>1 pc.+ adhesive tape 10x50cm</b>	<b>120 sets</b>
AT-SD-S 8	35cm x 120cm	1 pc.+ adhesive tape 10x50cm	100 sets

**KAMEROS UŽDANGALAS**

Tipas	Nr.
Standartinis tipas	170 101

- Dydis - 15x250 cm;
- Pagamintas iš polietileno 0,035mm su lipnia juosta fiksavimui;
- Teleskopiškai supakuotas;
- Vamzdelio apvalkalo optikos jungtis baigiasi trikampio forma, perforacija patogiam nuplėšimui;
- Sterilus

## Camera cover



Type	No.
Standart	170 101

- Size – 15 x 250 cm;
- Made of polyethylene 0,035mm with adhesive tape for fixing;
- Telescope folding;
- The tube shell optic connectors end up in a triangular, perforated, comfortable cut;
- Sterile

# Surgicare® Dual

Dual Protection **Powder Free** Sterile Surgical Gloves for High Risk & Trauma

Surgicare Dual is made from Natural Rubber Latex. The glove is anatomically shaped with thumb ball effect, reducing hand fatigue during surgical procedures. **Surgicare Dual provides primary and secondary protection for both doctor and patient.**

## EXCELLENT DONNING AND FIT

The specially designed mould produces gloves of high comfort.

## OTHER FEATURES

- ▶ Visual warning of a puncture on over glove in presence of fluid, indicator appears.
- ▶ **Recommended for High Risk, Orthopaedic & Trauma surgery.**
- ▶ Protects against pathogen transmission and is a barrier to cross infection.
- ▶ **Inner Polymer Coating provides smooth and easy donning of gloves. specialiu polimeru padengtas vidinis paviršius**
- ▶ The powderfree Under Glove protects the user against latex protein allergies. Protein content no more than 50µg of protein per gram and dm<sup>2</sup>.
- ▶ Manufactured in dust free environment with low bioburden and therefore reduced pyrogen.
- ▶ Sterilized in a validated process cycle as per ISO 11135:2014 norms SAL 10<sup>-6</sup>.

## STANDARDS

- ▶ ASTM : D-3577-2009 & EN 455-1:2000, EN 455-2:2015, EN 455-3:2015 & ISO 10282:2014
- ▶ CE Mark - MDD Class IIa
- ▶ Passes Biocompatibility test as per ISO 10993-10

## QUALITY CONTROL

- ▶ 100% Inspected
- ▶ **Factory Standard of 0.65AQL for Pinholes**

Recommended for:

- ▶ **High Risk, Orthopaedic & Trauma Surgery**



**Kanam**  
KANAM LATEX  
INDUSTRIES PVT. LTD.

The trusted name the world over

Ooppotttil Buildings, K.K. Road  
Kottayam-686 001, Kerala, India  
Tel: 0481-2560108, 2563513  
Fax: 0481-2563614  
E-mail: customercare@kanamlatex.com  
ho@kanamlatex.com

## Product Information

### TYPE

Sterile Latex **Dual** Surgical Gloves.

### PRIMARY MATERIAL

Natural Rubber Latex (NR).

### COLOUR

Under Glove: Green.

Over Glove: Creamy White.

### DESIGN

Hand Specific, Curved finger, Beaded cuff, Micro rough textured.

### POWDER FREE RESIDUE

<2mg/glove.

### STERILIZATION

Ethylene Oxide / Gamma as per Customer requirement.

### SHELF LIFE

3 Years in original packaging if stored as per storage instructions.

### PACKAGING

Material	Quantity	Dimension (cm)
Paper Wallet	2 Pairs (1 Unit)	20(L) x 11(W)
Pouch	2 Pairs (1 Unit)	27(L) x 15(W)
Inner Box	25 Units	27(L) x 15(W) x 20(H)
Master Carton	200 Units	61(L) x 28(W) x 41(H)
Container Quantities	20FCL: 75,600 units/40FCL: 1,59,600 units	

### PHYSICAL DIMENSIONS – Under Glove

Size	Length (mm)	Palm Width (mm)	Thickness (mm)		
			Cuff	Palm	Finger
5 ½	270	72 ± 4	0.14	0.17	0.19
6		77 ± 5			
6 ½		83 ± 5			
7	285	89 ± 5			
7 ½		95 ± 5			
8		102 ± 6			
8 ½		108 ± 6			
9		114 ± 6			

### PHYSICAL DIMENSIONS - Over Glove

Size	Length (mm)	Palm Width (mm)	Thickness (mm)		
			Cuff	Palm	Finger
5 ½	270	72 ± 4	0.14	0.17	0.19
6		77 ± 5			
6 ½		83 ± 5			
7	285	89 ± 5			
7 ½		95 ± 5			
8		102 ± 6			
8 ½		108 ± 6			
9		114 ± 6			

### PHYSICAL PROPERTIES - Under Glove & Over Glove

Characteristics	Before Ageing	After Ageing 70±2°C for 166±2hrs
Tensile Strength (Mpa) min.	24	18
Ultimate Elongation (%) min.	750	560
Stress at 500% Elongation (Mpa) Max.	5.5	-
Force at break (N) Min.	9.0	9.0

### PRE-SHIPMENT QUALITY INSPECTION

Sampling Procedure: ISO 2859 Part 1

Sampling Plan: Single Sampling, Normal Inspection

Characteristics	Inspection Level	AQL
Freedom from Holes	G-1	1.5
Visual attributes – Major Defects	G-1	1.5
Visual attributes – Minor Defects	G-1	2.5
Physical Dimension	S-2 / N=13	4.0
Physical Properties	S-2	4.0
Force at break	N=13	N/A
Extractable Protein Content	N=3	N/A
Powder free residue	N=5	N/A
Sterility	N=5	N/A

KLIP/G17/Rev#04



# Surgicare®

## Dual

Dual Protection Sterile Surgical Gloves for High Risk & Trauma

**KANAM LATEX INDUSTRIES PVT. LTD.****SURGICARE DUAL**

Surgicare Dual yra pagamintos iš natūralios gumos latekso. Pirštinė yra anatominės konfigūracijos bei su specialiu efektu sumažinančiu rankų nuovargį chirurginių procedūrų metu. **Surgicare Dual suteikia pirminę bei antrinę apsaugą ir pacientui ir gydytojui.** Specialiai sukurta forma padeda gaminti didelio komforto pirštines.

**Kitos savybės:**

- Vizualinio pradūrimo įspėjimo sistema, esant skysčiams.
- **Rekomenduojamos aukštos rizikos, ortopedinei bei traumų chirurgijai.**
- Apsaugo nuo patogenų perdavimo bei yra barjeras infekcijai
- Apatinė pirštinė **(be pudros)** apsaugo nuo latekso proteinų alergijos. Proteino kiekis ne daugiau nei  $50\mu\text{g}/\text{dm}^2$ .
- Pagaminta aplinkoje be dulkių su žemu biologinės taršos lygiu, dėl to sumažinamas pirogenų kiekis.
- Sterilizuotas patvirtintu proceso ciklu, pagal ISO 11135:2014

**Standartai:**

- ASTM : D-3577-2009 ir EN 455-1:2000, EN455-2:2015, EN 455-3:2015, ISO 10282:2014
- Turi CE sertifikatą – Medicinos prietaisų klasė IIa / Apsaugos priemonių kategorija III.
- Praeina biologinio suderinamumo testą pagal ISO 10993-10

**Kokybės kontrolė**

- 100% patikrintos
- **Gamyklos standartas AQL 0.65**

**Produkto Informacija**

- **Tipas:** Sterilios **Dvigubos** Chirurginės Pirštinės. Sterilizuotos ETO.
- **Pagrindinė medžiaga:** Natūralios gumos lateksas
- **Spalva:** Vidinė pirštinė – žalia, Išorinė pirštinė – kreminė/balta
- **Dizainas:** **dešinei ir kairei rankoms,** lenktos ties pirštais, susuktu rankogaliu, su mikrotekstūra.
- **Pudros likutis:** <2mg/pirštinei
- **Sterilizacija:** Etileno Oksidas/Gamma spinduliuotė (pagal kliento pageidavimą)

**Vidinė Pirštinė**

Dydis	Ilgis (mm) Min	Delno plotis (mm) Min	Storis (mm) Min		
			Rankogalis	Delnas	Pirštai
5½	270	72±4	0.14	0.17	0.19
6		77±5			
6½		83±5			
7	280	89±5			
7½		95±5			
8		102±6			
8½		108±6			
9		114±6			

### Išorinė piršinė

Dydis	Ilgis (mm) Min	Delno plotis (mm) Min	Storis (mm) Min		
			Rankogalis	Delnas	Pirštai
5½	270	72±4	0.14	0.17	0.19
6		77±5			
6½		83±5			
7	280	89±5			
7½		95±5			
8		102±6			
8½		108±6			
9		114±6			

### Fizinės Savybės

Savybės	Prieš Sendinimą	Po Sendinimo 70±2°C, 166±2 valandas
Stiprumas tempiant (Mpa) min.	24	18
Prailgėjimas (%) min.	750	560
Stiprumas prie 500% prailgėjimo (Mpa) Max.	5.5	-
Plyšimo riba (N) Min.	9.0	9.0

### Kokybės Patikrinimas Prieš Transportavimą

Tikrinimo procedūra: ISO 2859 Part 1

Tikrinimo planas: vienas pvz., įprastas patikrinimas


Savybės	Patikrinimo lygis	AQL
Skylių nebuvimas	G-1	1.5
Vizualiniai atributai – Pagrindiniai defektai	G-1	1.5
Vizualiniai atributai – Nedideli defektai	G-1	1.5
Fiziniai išmatavimai	S-2 / N=13	4.0
Fizinės savybės	S-2	4.0
Plyšimo riba	N=13	Nėra
Išgaunamas proteino likutis	N=3	Nėra
Pudros likutis	N=5	Nėra
Sterilumas	N=5	Nėra

**I  
p  
a  
k  
a  
v  
i  
m  
a**

<b>Medžiaga</b>	<b>Kiekis</b>	<b>Išmatavimai (cm)</b>
Popierinis vokas	2 poros (1 pakuotė)	20(L) x 11.5(W)
Atplėšiama pakuotė	2 poros (1 pakuotė)	27(L) x 15(W)
Vidinė dėžutė	25 Poros	27(L) x 15(W) x 22(H)
Dėžė	200 Poros	62(L) x 28.5(W) x 46(H)
Konteinerio kiekiai	20FCL: 72,000 pak/40FCL: 1,52,000 pakuotės	

# ESVI sensor

EKG ELEKTRODAI



## ES GS48

**NAUDOJIMAS: STRESO  
TESTUI, STEBĖJIMUI, EKG**

1. **Apdangalas** – viena pusė silikonizuota, PET dangalas
2. **Pedengimas** – tirštas gelis
3. **Laikiklis** - PE putos 1 mm storio
4. **Žyma** - PE, apvalus Ø24 mm, su užspaudu
5. **Jutiklis** - Ag/AgCl
6. **Smeigė** – nerūdijantis plienas
7. **Nuėmimas**– Popieriaus lapas, kuris palengvina elektrodo nuėmimą nuo dangalo

### PAPILDOMA INFORMACIJA :

Toksikologijos Direktyva EN ISO10993

Be latekso, vienakrtiniai, skirti suaugusiems

Tinkamumo laikas 24 mėnesiai

### Pakavimo metodas

Pakuotės kiekis (vnt) 50 Dėžės  
kiekis (vnt) 2000

Paletės kiekis (vnt) 112,000      Pakuotės  
medžiaga: PA / ALU / PE

**esvico**

DESIGN, PRODUCTION AND SALE OF MEDICAL DEVICES CLASS I  
- DISPOSABLE ECG ELECTRODES





OMS Kodas	Produkto aprašymas	Pavadinimas	Vidinė pakuotė	Išorinė pakuotė
1123	2,7g pakelis	OptiLube	144 vnt dėžutė	1440 vnt dėžė
1120	5g pakelis	OptiLube	150 vnt dėžutė	1500 vnt dėžė
<b>1125</b>	<b>6ml švirkštuose</b>	<b>OptiLube</b>	<b>25 vnt dėžutė</b>	<b>250 vnt dėžė</b>

OptiLube yra svarbus universalus gelis, rekomenduojamas šioms procedūroms:

- kateterių įvedimui
- spekulijų įvedimui
- rektaliniams termometrams
- klizmoms
- ginekologinėms apžiūroms
- tamponų įvedimui
- skaitmeninei ir instrumentinei apžiūrai
- prezervatyvų sutepimui
- bendroms hospitalinėms procedūroms



OptiLube™ sterilus lubrikantas pakeliai skirti individualiam naudojimui. Šios pakuotės panaikina užkrėtimo riziką ir yra idealiai tinkamos patogiam naudojimui. Galimas įpakavimas 2,7 g arba 5 g.



# DIGITAL THERMOMETER

## KFT-01/KFT-02



- Model No.: KFT-01
- Display: LCD,3.5digits, 19mmx7.5mm
- Temperature measurement: thermistor
- Apply for: Oral / Axillary / Rectal
- 40 seconds response
- sound warning when peak temperature is reached
- Fever Alarm

- Accuracy:  $\pm 0.1^{\circ}\text{C}$ ; ( $\pm 0.2^{\circ}\text{F}$ )
- Resolution:  $0.1^{\circ}\text{C}$  ( $0.1^{\circ}\text{F}$ )
- Measurement scope:  $32^{\circ}\text{C}$ - $42.9^{\circ}\text{C}$  ( $89.6^{\circ}\text{F}$ - $109.2^{\circ}\text{F}$ )  
Temperature  $< 32^{\circ}\text{C}$  ( $89.6^{\circ}\text{F}$ ) display  $\text{Lo}^{\circ}\text{C}$  ( $\text{Lo}^{\circ}\text{F}$ )  
Temperature  $> 42.9^{\circ}\text{C}$  ( $109.2^{\circ}\text{F}$ ) display  $\text{Hi}^{\circ}\text{C}$  ( $\text{Hi}^{\circ}\text{F}$ )

- Memory records: 1 sets
- Low voltage Indication
- Automatically turn-off function (8 min)

- Validity: 3 years, (Batteries are not included and the product can be used after calibration by a qualified body)
- Power source and voltage: DC 1.5V, LR41 or AG3
- Battery life: Approx.1200 times or 2 years
- Dimension:  $125\text{mm} \times 18\text{mm} \times 9.2\text{mm}$
- Unit weight: Approx. 10g including battery
- Work environment:  
Temperature  $10^{\circ}\text{C}$ - $40^{\circ}\text{C}$  ; Humidity: 15%RH~85%RH  
Atmospheric pressure: 70KPa~106KPa
- Transport/Storage condition:  
Temperature  $-20^{\circ}\text{C}$ - $60^{\circ}\text{C}$ ; Humidity: 10%RH~90%RH  
Atmospheric pressure: 70KPa~106KPa

- Model No.: KFT-02
- Display: LCD,3.5digits, 15.6mmx6mm
- Temperature measurement: thermistor
- Apply for: Oral / Axillary / Rectal
- 40 seconds response

- Accuracy:  $\pm 0.1^{\circ}\text{C}$ ; ( $\pm 0.2^{\circ}\text{F}$ )
- Resolution:  $0.1^{\circ}\text{C}$  ( $0.1^{\circ}\text{F}$ )
- Measurement scope:  $32^{\circ}\text{C}$ - $42.9^{\circ}\text{C}$  ( $89.6^{\circ}\text{F}$ - $109.2^{\circ}\text{F}$ )  
Temperature  $< 32^{\circ}\text{C}$  ( $89.6^{\circ}\text{F}$ ) display  $\text{Lo}^{\circ}\text{C}$  ( $\text{Lo}^{\circ}\text{F}$ )  
Temperature  $> 42.9^{\circ}\text{C}$  ( $109.2^{\circ}\text{F}$ ) display  $\text{Hi}^{\circ}\text{C}$  ( $\text{Hi}^{\circ}\text{F}$ )

- Waterproof level: IPX7
- Memory records: 1 sets
- Low voltage Indication
- Automatically turn-off function (8 min)
- Validity: 3 years

- Power source and voltage: DC 1.5V, LR41 or AG3
- Battery life: Approx.1500 times or 3 years
- Dimension:  $119\text{mm} \times 15.6\text{mm} \times 9.3\text{mm}$
- Unit weight: Approx. 7g including battery
- Work environment:  
Temperature  $10^{\circ}\text{C}$ - $40^{\circ}\text{C}$  ; Humidity: 15%RH~85%RH  
Atmospheric pressure: 70KPa~106KPa
- Transport/Storage condition:  
Temperature  $-20^{\circ}\text{C}$ - $60^{\circ}\text{C}$ ; Humidity: 10%RH~90%RH  
Atmospheric pressure: 70KPa~106KPa



COCET 可事达  
KANGFU MEDICAL

## SKAITMENINIS TERMOMETRAS

### KFT-01

- Modelio numeris.: KFT-01
- Ekranas: LCD,3.5 skaičių, 19mmx7.5mm
- Temperatūros matavimas: termistorius
- Galimas matavimas: burnje / pažastyje / tiesiojoje žarnoje
  
- rezultatas po 40 sekundžių
- garsinis įspėjimas, kai pasiekama aukščiausia temperatūra
- Karščiavimo signalas
  
- Tikslumas:  $\pm 01^{\circ}\text{C}$ ; ( $\pm 02^{\circ}\text{F}$ )
- Skyra:  $0.1^{\circ}\text{C}$  ( $0.1^{\circ}\text{F}$ )
- Matavimo apimtis:  $32^{\circ}\text{C}$ - $42.9^{\circ}\text{C}$  ( $89.6^{\circ}\text{F}$ - $109.2^{\circ}\text{F}$ )  
 Temperatura  $<32^{\circ}\text{C}$  ( $89.6^{\circ}\text{F}$ ) ekranas Lo $^{\circ}\text{C}$  (Lo $^{\circ}\text{F}$ )  
 Temperatura  $>42.9^{\circ}\text{C}$  ( $109.2^{\circ}\text{F}$ ) ekranas Hi $^{\circ}\text{C}$  (Hi $^{\circ}\text{F}$ )
  
- Atminties įrašymas: 1 atminties kartas
- Žemos įtampos indikacija
- Automatinio išsijungimo funkcija po 8 minučių
  
- Garantija: 3 metai, (Baterijos nėra įtrauktos ir gaminį galima naudoti kvalifikuotai institucijai sukalibravus)
- Energijos šaltinis ir galingumas: DC 1.5V, LR41 or AG3
- Baterijos gyvavimas: Vidutiniškai 1200 kartų arba 2 metus
- Išmatavimai: 125mm × 18mm × 9.2mm
- Vieneto svoris: Vidutiniškai 10g įskaitant bateriją
- Darbinė aplinka :  
 Temperatura  $10^{\circ}\text{C}$ ~ $40^{\circ}\text{C}$  ; Drėgmė: 15%RH~85%RH  
 Atmosferos slėgis : 70KPa~106KPa
- Transportavimo/laikymo sąlygos:  
 Temperatura  $-20^{\circ}\text{C}$ ~ $60^{\circ}\text{C}$  ; Drėgmė: 10%RH~90%RH  
 Atmosferos slėgis : 70KPa~106KPa



Paskirtosios įstaigos spaudas,  
nurodantis pavadinimą, adresą, telefoną  
(spaudas dedamas išrašius sertifikata)

## PATIKROS SERTIFIKATAS

Nr. 1848508

2023m. spalio 12d.

Data .....

1 (vienas)

Puslapių skaičius .....

Maksimalus elektrinis termometras KFT-01 Nr.MG.15440.1  
matavimo ribos (+32.0...+42.9)°C, II kl.

MP pavadinimas, tipas, Nr.,  
matavimo ribos, tikslumas .....

Skirgesa, UAB, Energetikų g. 8, LT-52461 Kaunas, Įm.kodas

MP savininko pavadinimas,  
įmonės kodas arba adresas .....

234449420,

BPM 8871101-103:2006

Patikros metodo žymuo .....

Skaitmeninis termometras XP100 Nr.028.1116.1212.006 dav. Pt100  
Nr.362111-1

Naudotos etaloninės MP  
pavadinimas, tipas, Nr. ....

ML90137 PR23-02588 2023m. spalio 12d. lipd. Nr.: 420698

Patikros protokolo  
registracijos Nr., data .....

UAB "MILGEDA" laboratorija

Matavimo priemonė patikrinta  
(laboratorijos pavadinimas) .....


Tinka naudoti

Išvada .....

2024m. spalio 11d.

Patikros sertifikatas galioja iki .....

(data)



Metrologas .....

(pareigų pavadinimas)




Pavel Purpurovič

(parašas)

(vardas ir pavardė)

Laboratorijos vadovas .....

(pareigų pavadinimas)




Sergej Voicechovskij

(parašas)

(vardas ir pavardė)

## Cure-Aid® - pirmos pagalbos pleistrai



### Aprašymas

Cure Aid – individualiai supakuotas pirmos pagalbos žaizdos pleistras su sugeriančia pagalvėle. Pleistro pagrindas gali būti pagamintas iš skirtingų medžiagų: polietileno plėvelės, poliuretano plėvelės, audinio, elastingo audinio, elastingos neaustinės medžiagos, putų pagrindo, spalvingomis dekoracijomis vaikams ir kt.

### Privalumai

- Itin sugerianti, nelimpanti prie žaizdos pagalvėlė;
- Itin hipoalergiški klijai;
- Itin gerai prilimpa;
- Individualiai supakuoti, gali būti sterilizuoti etileno oksidu pagal pageidavimą;
- Polietileno pagrindo pleistrai atsparūs skysčiams;
- Neaustinės medžiagos pleistrai turi mikroporas ir yra pralaidūs orui
- Perforuoti;
- Elastingi.

### Naudojimas


- Odos įbrėžimams ir nedidelių žaizdelių uždengimui;
- Tinka naudoti ant lanksčių kūno vietų.

### Dydžiai

- 9x38mm, 19x38mm, 18x40mm, 20 x 40mm, **19x72mm** ir kt.;
- Skirtingų išmatavimai atsižvelgiant į kliento pageidavimą.



## Cure-Aid® first aid dressing



### Description

Cure Aid - Individually packaged first aid wound patch with absorbent pad. The base of the patch can be made of different materials: polyethylene film, polyurethane film, fabric, elastic fabric, elastic non-woven material, foam base, colorful decorations for children, etc.

### Advantages


- Highly absorbent, non-sticky pad;
- Extremely hypoallergenic glue;
- Very good adhesion;
- Individually packed, sterilized with ethylene oxide upon request;
- Polyethylene-based patches are resistant to liquids;
- Non-woven based are has micropores and are air permeable
- Perforated;
- Elastic.

### Use

- For covering skin scratches and small wounds;
- Suitable for use on flexible areas of the body.

### Sizes

- 9x38mm, 19x38mm, 18x40mm, 20 x 40mm, **19x72mm** etc.
- Various sizes according to customer's needs.



### Technical data

-1. product name: Nasal Oxygen Cannula Disposable

vienkartinės

**2. Structure:** Nasal Oxygen Cannula is composed of nasal oxygen cannula and oxygen tube. The nasal oxygen cannula is consisted of the nose cannula, nasal oxygen tube, the adjusting ring and three way connector. The oxygen tube is composed of oxygen tube connectors and tube.

**3. Intended use and how to use:** It can be used in providing oxygen to emergency and seriously ill patients .

Before use, the user must read the usage manual.

- Attach the oxygen supply tubing to the gas source and set the oxygen to the prescribed flow.
- Check for gas flow through the device.
- Insert the nasal tips into the nostrils with the two plastic tubes over the ears and under the chin.
- Gently adjust the plastic slide until the cannula is secure.


NOTE: Nasal tips may be trimmed with scissors to fit smaller patients.

**4. Raw materials:** The materials is polyvinyl chloride (PVC).

**5. Product Specification:** L (large), M (middle), S (small) , XS (extra small).

### 6. drawings of the product:

atšakėles fiksuojanti atramėlė



	<b>DEGUONIES KANIULĖ</b> <b>sterili</b> <b>TD-24</b>	<b>Priedas Nr..3</b> <b>taškas 2 pp. 2.3</b>
		<b>Patikrinimas I</b> 2016.11.18
		Puslapis 1/1

**GAMINTOJAS:** ZARYS International Group sp. z o.o. sp.k.

**ADRESAS:** g. Pod Borem 18, 41-808 Zabrze, Lenkija, tel. +48 32 376 07 00-02, fax. +48 32 370 38 08

**INDIKACIJA:**

- Priemonė, skirta pernešti deguoniui per nosį

**SAVYBĖS:**

- Pagaminta iš PVC
- Labai minkšti antgaliai lygiais kraštais
- Tvirtinimas po smakru su minkštu, paslankiu reguliavimo žiedu
- Įvairios nosies dalies tarpų kaniulės, priklausomai nuo dydžio
- Universalus konektorius, suderinamas su bet kokių deguonies šaltiniu
- Žvaigždės skerspjūvio tipo vamzdelis, atsparus lenkimui
- Įmanomi 3 dydžiai: suaugusiems, vaikams, kūdikiams su 140 cm, 200 cm, 300 cm ir 500 cm vamzdeliais
- Be latekso
- EO sterilizuota
- Individualiai įpakuota: folija



**DYDŽIAI IR STANDARTAI**

	KATALOGO NUMERIS	VAMZDELIO ILGIS	DYDIS	TARPINĖ PAKUOTĖ	DIDELĖ PAKUOTĖ
80.2 p.d.	CTNN-200	200 cm	kūdikiams	25 vnt.	10 x 25 vnt.
80.2 p.d.	CTNZ-200	200 cm	vaikams	25 vnt.	10 x 25 vnt.
	CTND-140	140 cm	suaugusiems	30 vnt.	10 x 30 vnt.
	CTND-200	200 cm	suaugusiems	25 vnt.	10 x 25 vnt.
	CTND-300	300 cm	suaugusiems	15 vnt.	10 x 15 vnt.
	CTND-500	500 cm	suaugusiems	10 vnt.	10 x 10 vnt.