

Convincing Precision!



## **marSeal IQ** and **SealSafe® IQ**

*The system for bipolar vessel sealing  
for open and laparoscopic surgery*





Convincing Precision!


The *marSeal 5* instrument set at a glance

**marSeal 5**  
**marSeal 5 IQ**  
**marSeal 5 IQ short**

*marSeal 5 system*

Instrument, complete	
80-630-50-04	Version "IQ", shaft length 370 mm
80-630-90-04	Version "IQ" short, shaft length 200 mm
80-630-54-04	Version "m", shaft length 370 mm
80-630-55-04	Version "i"/"e", shaft length 370 mm
Components	
80-630-52-04	Handle "IQ"
80-630-56-04	Handle "m" version
80-630-57-04	Handle "i"/"e" version
80-630-51-04	Tubular shaft, 370 mm long, cleaning tool and flushing adapter
80-630-92-04	Tubular shaft for marSeal 5 IQ short, 200 mm long, cleaning tool and flushing adapter
80-630-53-04	Disposable blade insert, 370 mm long, SU = 6/pack
80-630-91-04	Disposable blade insert for marSeal 5 IQ short, 200 mm long, SU = 6/pack
80-630-83-04	Disposable blade insert, 370 mm long, sterile, SU = 10/pack
80-630-58-04	Cleaning tool
80-630-59-04	Flushing adapter

2. Elektrochirurginio instrumento darbinės dalies ilgis 200 mm



2. Elektrochirurginio instrumento darbinės dalies ilgis 200 mm



1. Siūloma elektrochirurginio instrumento darbinė dalis yra techniškai suderinama su LSMU ligoninėje Kauno klinikose naudojamu gamintojo „KLS Martin“ elektrochirurgijos prietaisu „MAXIUM-m“

Can be activated also via finger switch operation



*marSeal 5 instrument components*

	Handle, IQ-connector (HW 06) 80-630-52-04	Handle, KLS Martin connector 80-630-56-04	Handle, international connector 80-630-57-04	Tubular shaft* 370 mm 80-630-51-04	Tubular shaft* 200 mm 80-630-92-04	Blade shaft with blade 370 mm 80-630-53-04	Blade shaft with blade 200 mm 80-630-91-04
marSeal 5 cpl., 370 mm, IQ-connector (HW 06), 80-630-50-04	X			X		X	
marSeal 5 short cpl., 200 mm IQ-connector (HW 06), 80-630-90-04	X				X		X
marSeal 5 cpl., 370 mm, KLS Martin connector, 80-630-54-04		X		X		X	
marSeal 5 cpl., 370 mm, international connector, 80-630-55-04			X	X		X	

\*incl. blade shaft, cleaning tool and flushing adapter

