

Temporary Cements

TempBond

TempBond® is a self-curing zinc-oxide eugenol-based temporary cement indicated for temporary crowns, bridges or splints and for trial cementing permanent restorations. Its smooth flow means exceptional handling and results in effortless and complete restoration seating.

TempBond Introductory Pack (with modifier) 60058

Contents: 1 tube base (50 g), 1 tube accelerator (15 g), 1 tube modifier (13 g), 1 mixing pad

TempBond Standard Pack (without modifier) 61087

Contents: 1 tube base (50 g), 1 tube accelerator (15 g), 1 mixing pad

TempBond Unidose 31377

Contents: 50 unidose, 2.4 g each (total of 120 g)

TempBond Original Automix Syringe 33684

Contents: 2 Syringes (11.8 g each), 20 Mixing Tips

TempBond Refill Packs

Base material only (50 g) 60254

Accelerator only (15 g) 60255

Modifier only (13 g) 61081

TempBond® NE

This non-eugenol temporary cement will not inhibit the polymerisation of resin cements and acrylic temporaries and offers an option for patients allergic to eugenol. It provides the same flow and retentive properties of TempBond.

TempBond NE 60256

Non-Eugenol temporary zinc-oxide cement. Contents: 1 tube base (50 g), 1 tube accelerator (15 g), 1 mixing pad

TempBond NE Unidose 31973

Contents: 50 unidose, 2.4 g each (total of 120 g)

TempBond NE Automix Syringe 33685

Contents: 2 Syringes (11.7 g each), 20 Mixing Tips

TempBond Clear

A translucent, dual-cure, resin-based cement for temporary and provisional restorations. TempBond® Clear is eugenol free allowing for its use with resin luting cements and is easy to mix and clean up, leaving no residue. As it is resin based, a much stronger bond is achieved, but removal, when required, is fast and trouble-free.

TempBond Clear with Triclosan Automix Syringe 33351

Contents: 1 Syringe (6 g), 10 Mixing Tips

TempBond Accessories 33361

Automix Tip Refill (50 tips)



TempBond Automix Syringes



TempBond Unidose



TempBond Introductory Pack (with modifier)



TempBond NE