



**Poz. 19 Balionu išplečiama transfemoralinė
introduserio sistema**

SoloPath™

Balloon Expandable TransFemoral Access System

Micro to macro solutions
for vascular access



 **TERUMO**

SoloPath™

Provides the smallest insertion profile of any sheath in its class

- ▶ Designed to minimize insertion forces
- ▶ Reduced profile enables easier access even in difficult anatomy
- ▶ Designed to track vessels similar to a guidewire

Sheath is dilated using Controlled Deployment Technology (CDT™)

- ▶ Provides a protective conduit
- ▶ Uses natural elasticity of tissue & vessel
- ▶ Facilitates more predictable closure

Unique malleable design

- ▶ Designed to protect the vessel
- ▶ Provides enhanced maneuverability

Order details

SoloPath™ Balloon Expandable TransFemoral Introducer

Model	Description	I.D.	O.D.	Expandable Length	Tip Diameter	Working Length
STFI-1825	SoloPath TransFemoral 18 Fr x 25 cm	18 Fr / 6.0 mm	21 Fr / 7.0 mm	29 cm	5.3 Fr / 1.8 mm	25 cm
STFI-1835	SoloPath TransFemoral 18 Fr x 35 cm	18 Fr / 6.0 mm	21 Fr / 7.0 mm	39 cm		35 cm
STFI-1925	SoloPath TransFemoral 19 Fr x 25 cm	19 Fr / 6.3 mm	22 Fr / 7.3 mm	29 cm		25 cm
STFI-1935	SoloPath TransFemoral 19 Fr x 35 cm	19 Fr / 6.3 mm	22 Fr / 7.3 mm	39 cm		35 cm
STFI-2125	SoloPath TransFemoral 21 Fr x 25 cm	21 Fr / 7.0 mm	24 Fr / 8.0 mm	29 cm		25 cm
STFI-2135	SoloPath TransFemoral 21 Fr x 35 cm	21 Fr / 7.0 mm	24 Fr / 8.0 mm	39 cm		35 cm



Terumo Corporation
44-1, 2-chome, Hataegaya
Shibuya-ku, Tokyo 151-8572,
Japan
Tel.: +81 3 3374 8111
Fax: +81 3 3374 8196

European Headquarters

Terumo Europe N.V.
Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 38 12 11
Fax: +32 16 40 02 49

**Terumo Europe
Interventional Systems
Liaison Office**
Immenhölz Kröbbs
8, rue Jean-Pierre Timpaud
CS 16720
78057 Saint Quentin
en Yvelines Cedex, France
Tel.: +33 139 30 95 75
Fax: +33 139 30 95 70

European Offices

**Terumo Europe N.V.
Benelux Sales Division**
Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 38 25 80
Fax: +32 16 39 25 99
The Netherlands:
Tel.: 0800 0220396
Fax: 0800 0220414

Laboratoires Terumo France S.A.
Parc Ariane - Bât. Uranus
1, rue Hélène Boucher
78284 Guyancourt Cedex, France
Tel.: +33 139 95 13 00
Fax: +33 139 43 60 85

Terumo Italia S.r.l.
Via Simone Martini, 143/145
00142 Roma RM, Italy
Tel.: +39 0651 95 14 20
Fax: +39 0650 30 407

Terumo UK Ltd
Tamesis, The Causeway
Egham, Surrey
TW20 9AW, United Kingdom
Tel.: +44 1784 476 200
Fax: +44 1784 476 201

Terumo Deutschland GmbH
Ludwig-Erhard-Strasse 6
65760 Eschborn, Germany
Tel.: +49 6196 80 230
Fax: +49 6196 80 23 200

**Terumo Deutschland GmbH
Zweigniederlassung Österreich**
Liebmannsstrasse 110 301
2345 Brunn am Gebirge, Österreich
Tel.: +43 2236 378 020
Fax: +43 2236 378 020 100

**Terumo Deutschland GmbH
Zweigniederlassung Soreitenbach**
Bodenackerstrasse 3
8957 Spreitenbach, Switzerland
Tel.: +41 56 418 10 10
Fax: +41 56 419 10 11

Terumo Europe España S.L.
Avda. Juan Carlos I, Nº 19 - 7ª Planta
Edificio Torre La Garena
28806 Alcalá de Henares
Madrid, Spain
Tel.: +34 9021 01 298
Fax: +34 9021 01 358

Terumo Denmark
Ejllal of Terumo Sweden AB
Tel.: +45 7020 93 80
Fax: +45 7020 94 80

Terumo Sweden AB
Nya Varvet, Byggnad 80
Sven Kallfelts Gäta 18
426 71 Västra Frölunda, Sweden
Tel.: +46 3174 85 880
Fax: +46 3174 85 890

**Terumo Europe N.V.
Moscow Representative Office**
90 Northern Tower
13Fr., 19 Teatovskaya Street
Moscow 123317, Russia
Tel.: +7 495 888 4740
Fax: +7 495 888 4739

**Terumo Corporation Dubai
Branch
(Middle East and Africa)**
Al Masraf Tower, 22nd Floor
P.O. Box 20281 Dubai
U.A.E.
Tel.: + 971 4 221 2220
Fax: + 971 4 221 3330

Addresses subject to change.
Please consult our Website.

© Registered Trademark
Published by Terumo Europe N.V.