



# Hi-Torque Steelcore

.018 SUPPORTIVE DELIVERY GUIDE WIRE

.018 | GUIDE WIRE

10 p.d Katalogas

## THE .018 SUPPORT YOU NEED.

Support and flexibility for a successful crossing and device delivery.



### Extra support and flexibility

- Supportive stainless steel shaft for excellent tactile feedback
- 1:1 torque for navigation around vessel bends
- Provides proximal through distal support

### Multiple options

- Available in two transition lengths; Steelcore LT has a longer core taper than Steelcore for increased flexibility and steering
- Shapeable platinum/nickel tip coil
- Wire lengths of 190 cm and 300 cm



Information contained herein for DISTRIBUTION outside of the U.S. ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.



# Hi-Torque Steelcore

## .018 SUPPORTIVE DELIVERY GUIDE WIRE

.018 | GUIDE WIRE

### ORDERING INFORMATION

DESCRIPTION	STOCK NUMBER	WIRE SIZE (in)	WIRE LENGTH (cm)	TIP SHAPE	TIP TAPER
Hi-Torque 018 Steelcore 190 cm	1003281	0.018	190	Straight	Straight
Hi-Torque 018 Steelcore 300 cm	1003282	0.018	300	Straight	Straight
Hi-Torque 018 Steelcore LT 190 cm	1007709	0.018	190	Straight	Straight
Hi-Torque 018 Steelcore LT 300 cm	1007710	0.018	300	Straight	Straight
Hi-Torque 018 Steelcore LT J 190 cm	1007709-J	0.018	190	J	Straight
Hi-Torque 018 Steelcore LT J 300 cm	1007710-J	0.018	300	J	Straight

#### Abbott Vascular International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use inside the product carton (when available) or at [efu.abbottvascular.com](http://efu.abbottvascular.com) for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. Information contained herein for DISTRIBUTION outside of the U.S. ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Tests performed by and data on file at Abbott Vascular unless otherwise noted. Illustrations are artist's representations only and should not be considered as engineering drawings or photographs.

Hi-Torque Steelcore is a trademark of the Abbott Group of Companies.

[www.AbbottVascular.com](http://www.AbbottVascular.com)

©2016 Abbott. All rights reserved. AP2942010-OUS Rev. A

