




10 p.d Katalogas

.035 STEERABILITY AND SUPPORT.

Place and exchange devices with support and confidence.



Purpose and support

- Facilitates the placement and exchange of interventional devices during procedures
- Supportive .035 wire with soft atraumatic tip
- 1:1 torque response provides exceptional steering
- Hydrophobic coating to reduce friction
- Wire lengths of 145 cm, 190 cm, and 300 cm

Efficiency

- Reduces multiple .035 wire exchanges during procedure
- Facilitates catheter placement for diagnostic and contralateral approach interventions





Hi-Torque Supra Core

.035 WORKHORSE & SPECIALTY GUIDE WIRES

.035 | GUIDE WIRE

ORDERING INFORMATION

DESCRIPTION	STOCK NUMBER	WIRE SIZE (in)	WIRE LENGTH (cm)	TIP TAPER
Hi-Torque Supra Core .035 145 cm	1002703	0.035	145	Straight
Hi-Torque Supra Core .035 190 cm	1002703-01	0.035	190	Straight
Hi-Torque Supra Core .035 300 cm	1002703-02	0.035	300	Straight

Abbott Vascular International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use inside the product carton (when available) or at efu.abbottvascular.com for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. Information contained herein for DISTRIBUTION outside of the U.S. ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Tests performed by and data on file at Abbott Vascular unless otherwise noted. Illustrations are artist's representations only and should not be considered as engineering drawings or photographs.

Hi-Torque Supra Core is a trademark of the Abbott Group of Companies.

www.AbbottVascular.com

©2016 Abbott. All rights reserved. AP2942014-OUS Rev. A

