

 Search for products, documents, pages, or articles



[Home](#) / [PRODUCTS](#) / [PERFORMER™ INTRODUCER AND SETS](#)

[Share this page](#)

# Performer™ Introducer and Sets

[SPECIFICATIONS](#)

[IMAGES](#)

Used to introduce balloons, closed and non-tapered end catheters or other diagnostic and interventional devices.

**This product line is serviced by the following clinical division(s):**

[Aortic Intervention »](#)

[Peripheral Intervention »](#)

Instructions for Use (IFUs) linked below are the most current version associated with the product. To find an IFU associated with a specific Lot Number, please visit <https://ifu.cookmedical.com/>.

Feedback

[Cookie Settings](#)

**LEARN HOW TO ORDER »**

Set up a customer account and order this or other Cook products

**ASK A SUPPORT QUESTION »**

Ask Customer Support about our products or services

**FIND PRODUCT TRAINING »**

Connect with our Vista® Education & Training program

**VIEW YOUR ACCOUNT INFORMATION »**

Register your customer number on our Customer Portal

 **Specifications**

COLLAPSED EXPANDED

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Fr	Length cm	Accepts Wire Guide Diameter inch
<b>Introducer Only</b>						
G07427	RCFW-5.0-35	–	N/A	5	13	.035
G07600	RCFW-10.0-38		N/A	10	13	.038
G07601	RCFW-12.0-38		N/A	12	13	

[Cookie Settings](#) 

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Fr	Length cm	Accepts Wire Guide Diameter inch
<b>Introducer with a Radiopaque Band</b>						
G07919	RCFW-6.0-35-RB		N/A	6	13	.035
G08550	RCFW-6.0-38-30-RB		N/A	6	30	.038
G07707	RCFW-7.0-38-30-RB		N/A	7	30	.038
G08139	RCFW-8.0-38-30-RB		N/A	8	30	.038
G08686	RCFW-10.0-38-30-RB		N/A	10	30	.038
G08955	RCFW-10.0-38-45-RB		N/A	10	45	.038
G08956	RCFW-12.0-38-30-RB		N/A	12	30	

[Cookie Settings](#)

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Fr	Length cm	Accepts Wire Guide Diameter inch
G28433	RCFW-12.0-38-45-RB		N/A	12	45	.038
G10407	RCFW-12.0-38-80-RB		N/A	12	80	.038
G08957	RCFW-14.0-38-30-RB		N/A	14	30	.038
G28438	RCFW-14.0-38-45-RB		N/A	14	45	.038
G09843	RCFW-14.0-38-80-RB		N/A	14	80	.038
G09691	RCFW-16.0P-38-30-RB		N/A	16	30	.038
G28445	RCFW-16.0P-38-45-RB		N/A	16	45	.038
G09582	RCFW-16.0P-38-70-RB		N/A	16	70	.038

Cookie Settings

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Fr	Length cm	Accepts Wire Guide Diameter inch
G09695	RCFW-16.0P-38-80-RB		N/A	16	80	.038
G11672	RCFW-18.0P-38-30-RB		N/A	18	30	.038
G27025	RCFW-18.0P-38-85-RB		N/A	18	85	.038

#### Introducer Set that Includes a Wire Guide

G07605	RCF-5.0-35-J	—	N/A	5	13	.035
G07487	RCF-5.0-38-J	—	N/A	5	13	.038
G27242	RCF-5.0-35-45-J		N/A	5	45	.035
G07376	RCF-6.0-35-J	—	N/A	6	13	.035

[Cookie Settings](#)

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Fr	Length cm	Accepts Wire Guide Diameter inch
G07338	RCF-7.0-38-J	–	N/A	7	13	.038
G07490	RCF-8.0-38-J	–	N/A	8	13	.038
G07564	RCF-9.0-38-J		N/A	9	13	.038
G07417	RCF-10.0-38-J		N/A	10	13	.038
G07369	RCF-12.0-38-J		N/A	12	13	.038
G08024	RCF-14.0-38-J		N/A	14	13	.038

### Micropuncture® Introducer Access Set that Includes a Wire Guide and a Needle

G12664	RCFN-4.0-18-MPIS-NT		N/A	4	13	.018
--------	---------------------	---	-----	---	----	------

[Cookie Settings](#)



Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Fr	Length cm	Accepts Wire Guide Diameter inch
G12424	RCFN-5.0-18-MPIS-NT		N/A	5	13	.018
G12423	RCFN-6.0-18-MPIS-NT		N/A	6	13	.018
<b>Introducer Access Set that Includes a Wire Guide and a Needle</b>						
G10025	RCFN-5.0-35-J		N/A	5	13	.035
<b>Radial Introducer Access Set that Includes a Wire Guide and a Needle</b>						
G27022	RCFN-4.0-18-5.5-RA1.5		N/A	4	5.5	.018
G21107	RCFN-4.0-18-13-RA2.5		N/A	4	13	.018
G10692	RCFN-4.0-18-23-RA1.5		N/A	4	23	.018

[Cookie Settings](#)

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Fr	Length cm	Accepts Wire Guide Diameter inch
G12018	RCFN-5.0-18-5.5-RA1.5		N/A	5	5.5	.018
G10268	RCFN-5.0-18-7-RA2.5		N/A	5	7	.018
G10266	RCFN-5.0-18-13-RA2.5		N/A	5	13	.018
G10267	RCFN-5.0-18-23-RA2.5		N/A	5	23	.018
G10387	RCFN-6.0-18-7-RA2.5		N/A	6	7	.018
G10270	RCFN-6.0-18-13-RA2.5		N/A	6	13	.018
G10271	RCFN-6.0-18-23-RA2.5		N/A	6	23	.018

PI-D56839-EN

## Image Library



Not all products shown on this website may be approved or available in all jurisdictions. Consult with your local Cook representative, distribution company or customer support centre for details.

© Cook 2024

 **+353 61 334 440**

[Products](#)  
[Patient Resources](#)  
[Support](#)  
[About](#)  
[Divisions](#)  
[Contact](#)

[Privacy Statement](#)  
[Terms of Use](#)  
[Human Rights](#)  
[Responsibility](#)  
[Statement](#)  
[Cookie](#)  
[Compliance](#)

**Sign Up. Get Email about News, Products, and Events from Cook Medical.**

[SIGN UP](#)

[Cookie Settings](#)

[Education & Training](#)  
[Customer Portal](#)

[Sitemap](#)

[Cookie Settings](#)

