

economical applicator

[Product Information](#)[Instructions](#)[Reorder Codes](#)[Download](#)

A basic applicator for all single dose capsule systems

- Autoclavable
- Economical

Economical Complet Applicator
8100120

14 pirkimo dalis



Shop by Category

Compomers

15 pirkimo dalis

Shop by Procedure

Shop by Brand

▶ Preventive

▼ Restorative

▼ Direct Restoration

▼ Composites & Flowables

Universal Composites

Flowables

Compomers

Bulk Fill Composites

▼ Refine by

Brand

Shade (14)

Sort by: Product Title

▼ Dyract A2 Refill Package - 20pk

Buy Now



Dyract, The Evidence-Based Restorative

Read more →

SKU: 60604202

▼ Dyract A3 Refill Package - 20pk

Buy Now



Dyract, The Evidence-Based Restorative

Read more →

Top Dent Mixing sheets

TOP DENT

18 pirkimo dalis

Polyethylene-coated paper pad.

		Art nr
4 x 5 cm	st	415 084
7 x 9 cm	st	415 085
14 x 20 cm	st	415 086



Apsaugoti

3M™ Clinpro™ XT ilgalaikė danga



24 pirkimo dalis

Ilgalaikė apsauga

Apsauginis barjeras, iškart sumažinantis jautrumą, atpalaiduojantis fluoridą ir išliekantis šešis mėnesius.

Clinpro XT ilgalaikė danga yra pažeistoms sritims skirta šviesa kietinama patvari medžiaga iškart suteikianti apsauginį sluoksnį padidėjusiam dantino jautrumui palengvinti.


- ▶ Puikiai tinka atviriems, jautriems šaknų paviršiams apsaugoti, ypač tik ką nuvalytiems paviršiams po profesionalios burnos higienos.
- ▶ Lieka ant dantų šešis mėnesius arba ilgiau
- ▶ Atsparus drėgmei: nereikia džiovinti ar paruošti šaknies paviršiaus
- ▶ Atpalaiduoja fluoridą, kalcį ir fosfatą
- ▶ Iškart sumažina jautrumą

Clinpro XT lakas apsaugo dantų paviršių nuo rūgščių poveikio

Clinpro XT lako in vitro tyrimai rodo, kad danga:

- ▶ Sukuria barjerą, saugantį nuo demineralizacijos ir rūgščių sukeltos erozijos
- ▶ Ne tik padengia, bet ir apsaugo
- ▶ Padeda atkurti demineralizuotą emalį šalia dangos

Clinpro XT ilgalaikė danga užsandarina ir apsaugo dantino kanalėlius




Atvirų dantino kanalėlių (dešinėje) ir Clinpro XT ilgalaikės dangos (kairėje) SEM vaizdai, padidinti 3000x, parodo barjero efektą.



Clinpro XT ilgalaikės dangos skerspjūvio SEM vaizde, padidintame 1000x, matomas padengtas dantinas ir dervos įsiskverbimas į kanalėlius.


Kaupiamasis fluoridų atpalaidavimas* [Vidurkis ± SD] (trumpalaikis)



—◆— 3M™ Clinpro™ XT ilgalaikė danga
—■— Įprasti fluoro lakai

In vitro tyrimai rodo, kad Clinpro XT ilgalaikė danga per 24 val. atpalaiduoja daugiau fluoridų, nei įprastas lakas.

Kaupiamasis fluoridų atpalaidavimas* [Vidurkis ± SD] (ilgalaikis)



—◆— 3M™ Clinpro™ XT ilgalaikė danga
—■— Silantai stiklo jonomerų pagrindu
—▲— Derviniai fluoro silantai

In vitro tyrimai rodo, kad Clinpro XT ilgalaikė danga per šešis mėnesius atpalaiduoja daugiau fluoridų, nei kitos dangos.

*3M Oral Care vidiniai duomenys.

27 pirkimo dalis

ViscoStat™ Clear
25% ALUMINUM CHLORIDE



2018¹
REALITY
Four Star Award
★★★★



Metal Dento-Infusor™ Tip

- Does not discolor the gingiva
- Stops minor bleeding in the esthetic zone
- Rinses easily
- Viscous gel
- Does not interfere with bonding²

Procedure



1. Subgingival preparation and bleeding sulcus.



2. Rub hemostatic firmly against bleeding tissues with the Metal Dento-Infusor™ tip. The clear gel allows easy visibility and rinses away quickly.



3. Press soaked Ultrapak™ cord into the sulcus. Leave 4–5 minutes.



4. Remove cord. Firm air/water spray. Air dry. Rub hemostatic into the sulcus again. Leave 1 minute.



5. Finished restoration 2 weeks post-op. Facilitates great control in esthetic zone with no gingival stain.

33 pirkimo dalis

File-Eze™ EDTA Lubricant

FILE LUBRICANT



Certified
GF
Gluten-Free

NaviTip 30 ga




- Peroxide free; will not affect the set of resin sealers

File-Eze file lubricant is an effective 19% EDTA in a water-soluble, viscous solution for chelating, lubricating, and debriding root canal preparations.

Ultradent™ Citric Acid 20% Solution



35 pirkimo dalis



NaviTip FX 30 ga

- Recommended as a cleanser/conditioner of prepared root canals
- Removes mineral and smear layers
- Slightly viscous formula facilitates lubrication
- Removes calcium hydroxide paste

Ultradent Citric Acid is a mild acidic material that is effective at dissolving/cleaning calcium hydroxide from canals (e.g. UltraCal XS). It is also recommended as a cleanser/conditioner to remove smear layer from dentinal walls. Deliver with the NaviTip or the NaviTip FX.

Dyract flow



36 pirkimo dalis



Get Free Shipping on Online Orders of \$500+ *Exclusions Apply

Discover

Shop

Academy

Dentsply Sirona

USA

Contact

Sign In For Price €

Click Image to Open Expanded View

Dyract
[View all Dyract Produc](#)

Product Description

Dyract flow flowable compomer restorative offers sustained, long-term fluoride release and low net volumetric shrinkage for exceptional marginal integrity.

- It is suitable for use as a base/liner under Class I and II restorations and can be used to block out undercuts and restore Class V preps.
- It is formulated to provide easy handling, high surface hardness, exceptional flexural modulus and high elasticity as well as optimal opacity for natural esthetics.

[Product Information](#)
[Instructions](#)
[Reorder Codes](#)
[Download](#)

The tooth remineralizing restorative

Ideal for minimally invasive dentistry


Riva Self Cure is the ideal restorative material for use in minimally invasive dentistry (MID) as it is a bioactive material which prevents caries from occurring. Natural tooth structure is preserved, so large cavities and undercuts are not necessary.

Great marginal adaptation

Riva Self Cure is the answer to long term microleakage prevention. Use it as a restorative or a base under permanent restorations. "...self-cured glass ionomer bases were more effective in reducing microleakage in both the occlusal and proximal cavo surfaces than a flowable resin." (1)



51, 54 pirkimo dalys



Riva Self Cure Capsules Regular Set

- 50 Riva Self Cure Capsules

8600001 A1

8600002 A2

8600003 A3

8600004 A3.5

8600005 A4

8600007 B2

8600008 B3



- 6.9mL (8g) Riva Self Cure Liquid bottle
- 15g Riva Self Cure Powder jar
- accessories

8610501 A1

8610502 A2

8610503 A3

8610504 A3.5

8610505 A4

8610507 B2

8610508 B3

Top Dent Activa BioActive Restorative

TOP DENT by PULPDENT

Triple-curing **bioactive** restorative material with physical properties including shock-absorbing capacity as well as excellent wear resistance and strength. The material is free of bisphenol A (and derivatives), radiopaque, hydrophilic and displays ion (fluoride, phosphate, calcium) release/uptake properties mimicking tooth substance (physiological situation)

ACTIVA BioACTIVE RESTORATIVE facilitates clinical work through independence of separate bonding procedures.

The material stimulates remineralization and apatite formation which seals the marginal gap between restorative material and tooth substance and thereby prevents Micro-leakage and secondary caries.

Can be used for Bulk-fill technology.

Two-component material in automix syringe á 5 ml (8 g), extruded with dispenser.

Shelf-life: 2 years. Can be stored in room temp.



55, 56 pirkimo dalys

	Art nr
Start kit A1 , 1 syringe, dispenser, 20 automix tips with bendable 20G metal tip1 st	fp 681 272
Start kit A2 , 1 syringe, dispenser, 20 automix tips with bendable 20G metal tip1 st	fp 681 273
Start kit A3 , 1 syringe, dispenser, 20 automix tips with bendable 20G metal tip1 st	fp 681 274
Refill A1 , 1 syringe, 20 automix tips with bendable 20G metal tip1 st	fp 681 275
Refill A2 , 1 syringe, 20 automix tips with bendable 20G metal tip1 st	fp 681 276
Refill A3 , 1 syringe, 20 automix tips with bendable 20G metal tip1 st	fp 681 277
Refill A3,5 , 1 syringe, 20 automix tips with bendable 20G metal tip1 st	fp 681 278
Valuepack A1 , 2 syringe, 40 automix tips with bendable 20G metal tip2 st	fp 681 279
Valuepack A2 , 2 syringes, 40 automix tips with bendable 20G metal tip2 st	fp 681 280
Valuepack A3 , 2 syringes, 40 automix tips with bendable 20G metal tip2 st	fp 681 281
Valuepack A3,5 , 2 syringes, 40 automix tips with bendable 20G metal tip2 st	fp 681 282

Top Dent Activa BioActive Base/Liner

TOP DENT by PULPDENT

Triple-curing bioactive base/liner restorative material with physical properties including shock-absorbing capacity as well as excellent wear resistance and strength characteristics. The material is free of bisphenol A (and derivatives), radiopaque, hydrophilic while displaying ion (fluoride, phosphate, calcium) release/uptake properties mimicking tooth substance (physiological situation)

ACTIVA BioACTIVE Base/Liner facilitates clinical work through independence of separate bonding procedures and display mild intrinsic self-etching properties. The material stimulates remineralization and apatite formation which seals the marginal gap between restorative material and tooth substance and thereby prevents micro-leakage and secondary caries.

Dentin-colored two-component material in automix syringe á 5 ml (7 g).

Shelf-life: 2 years. Can be stored in room temperature



57 pirkimo dalis

		Art nr
Syringe 4g 2 st	fp	679 307

Ultra-Blend™ plus

LINER



58 pirkimo dalis



- Bioactive¹ liner and pulp-capping material
- Superior calcium release²
- Light curable
- Controlled, precise syringe delivery
- No mixing necessary
- Will not dissolve over time
- Radiopaque
- Highly filled

Ultra-Blend plus liner is a light-activated radiopaque material with calcium hydroxide in a urethane dimethacrylate (UDMA) base. It's perfect for pulp capping and will not dissolve over time. Ultra-Blend plus liner is highly filled for minimal shrinkage.

UltraSeal XT™ hydro™

HYDROPHILIC PIT AND FISSURE SEALANT



59 pirkimo dalis

UltraSeal XT™ hydro™
Inspiral™ Brush Tip

- It is hydrophilic before it is cured, hydrophobic once cured, and has a self-adhesive quality
- Advanced adhesive technology
- Fluoresces under black light to ensure sealant is still in place
- Highly filled resin (53%)
- Thixotropic/ideal viscosity
- Two shades: Opaque White and Natural

UltraSeal XT hydro hydrophilic pit and fissure sealant is a light-cured, radiopaque, fluoride-releasing composite sealant. It is stronger and more wear resistant because it is a 53%-filled resin and has less polymerization shrinkage than competitive products. The spiral brush action of the tip causes shear thinning of the thixotropic UltraSeal XT hydro sealant. The resin stops-flowing when shear thinning ceases and placement is complete, preventing it from running before it can be light cured. The advanced hydrophilic chemistry works when all visible moisture has been removed, UltraSeal XT hydro sealant is more forgiving of moisture deep inside pits and fissures.

Adper™ Single Bond 2 Adhesive

61 pirkimo dalis

General Information

Adper Single Bond 2 Adhesive, manufactured by 3M ESPE, is a simple, moist bonding adhesive containing 10%, 5nm colloidal filler.

Adper Single Bond 2 adhesive offers the dental practitioner a wide range of applications. These include bonding to all classes of direct composite restorations as well as procedures involving porcelain, composite, metal repair, set amalgam, root surface desensitization and bonding of porcelain veneers with RelyX™ Veneer Cement System and RelyX™ Ceramic Primer, manufactured by 3M ESPE.

After light curing Adper Single Bond 2 adhesive, it may also be used for amalgam and indirect bonding procedures when combined with RelyX™ ARC Adhesive Resin Cement, manufactured by 3M ESPE. Compatibility with indirect bonding procedures is due to the low film thickness (approximately 10 µm) of cured Adper Single Bond 2 adhesive.

Adper Single Bond 2 adhesive is available in two delivery systems, a unit dose delivery and a multi-use vial dispenser.

Use of etchant is critical for both enamel and dentin surfaces.




riva light cure

[Product Information](#)
[Instructions](#)
[Reorder Codes](#)
[Download](#)

The glass ionomer restorative with true esthetics

Great translucency

*all shades A2



(Riva Light Cure HV)



ionglass

Riva Light Cure Powder / Liquid Kits

- 7.2mL (8g) Riva Light Cure Liquid bottle
- 15g Riva Light Cure Powder jar
- accessories

8700501 A1

8700502 A2

8700503 A3

8700504 A3.5

8700505 A4

8700507 B2

8700508 B3

8700510 B4

8700512 C2

8700514 C4

8700516 Bleach

62 pirkimo dalis

