

www.acteongroup.com



To visualize our Youtube channel:
www.youtube.com/user/acteonsatelec

SATELEC® S.A.S I A Company of ACTEON® Group
17 av. Gustave Eiffel BP 302161 33708 MERIGNAC cedex FRANCE
Tel +33 (0) 556 340 607 Fax +33 (0) 556 349 232
E-mail: info@acteongroup.com www.acteongroup.com



EN

Non contractual document - Réf. D00601 - V8 - 06/2016 - Copyright © 2014 SATELEC - All rights reserved. No information or part of this document may be reproduced or transmitted in any form without the prior permission of ACTEON.

I AM GARING

DENTAL TIPS
Preserve teeth





AS3D
Ref. F00065

Diamond-coated universal tip 30µm, length 3mm, taper 9%

The AS3D tip is intended for apical surgery of anterior teeth. It should be used without pressure, at the lowest possible effective power.



AS9D
Ref. F00067

Diamond-coated tip 30µm, length 9mm, taper 8%

The AS9D tip is for complex cases and allows preparation of the root canal up to the coronal third. The diamond coating of the AS9D is only present on the extremity of the instrument not to over-prepare the canal.



The AS9D tip should first be introduced into the canal and oriented in the root axis before being activated to prevent the creation of a "false route".



ASLD
Ref. F00080

Left-oriented tip, diamond-coated 30µm, length 3mm, taper 10%

Recommended for apical surgery of premolars and molars. This instrument should be used with very light pressure.



AS6D
Ref. F00079

Diamond-coated tip 30µm, length 6mm, taper 9%

Second instrument in the sequence, it should be used to obtain a preparation length of 5mm at least.



ASRD
Ref. F00081

Right-oriented tip, diamond-coated 30µm, length 3mm, taper 10%

The ASRD tip is recommended for apical surgery of premolars and molars.



The Apical Surgery kit, with its unique 3-c-9mm concept, offers a controlled retrograde endodontic treatment with greater preservation of bone and dental tissue.



EndoSuccess Apical Surgery Kit (Ref. F000059): supplied with AS3D, AS6D, AS9D, ASLD and ASRD tips, a metal holder and an autoclavable universal wrench.



S12-70D
Ref. F00118

Retro surgery tip angled at 70°, diamond-coated 30µm, length 5mm, taper 9%

The S12-70D tip is recommended for the treatment of posterior areas, in canals that are difficult to access or roots with specific orientations.



P14D
Ref. F00106

Universal retro surgery tip, diamond-coated 30µm, length 5mm, taper 7%

The P14D tip is recommended for the preparation of canals in anterior teeth.



The micro-retro tips make minimum treatment possible providing fast healing.



P15LD
Ref. F00107

Left-oriented retro surgery tip, diamond-coated 30µm, length 5mm, taper 7%

The P15LD is recommended for the preparation of premolar and molar canals.



P15RD
Ref. F00108

Right-oriented retro surgery tip, diamond-coated 30µm, length 5mm, taper 7%

The P15RD is recommended for the preparation of premolar and molar canals.



Endodontija

Apikalinė chirurgija

„Endosuccess Apical Surgery“

1. Kodas F00065
AS3D

30 µm deimantu dengtas universalus instrumentas, ilgis 3 mm, kūgis 9%

Instrumentas AS3D skirtas priekinių dantų apikalinei chirurgijai. Jis turi būti naudojamas be spaudimo, su mažiausiu galimu efektyviu galingumu.