

# Rapid Rhino<sup>®</sup> for Epistaxis

Controls epistaxis in minutes



## Rapid Rhino is the epistaxis solution

Combines gentle tamponade with platelet aggregation for fast, enduring epistaxis relief that maximizes patient comfort

- Easy to insert and remove due to self-lubrication and ultra-low profile
  - Minimizes rebleeds because it does not stick to clots upon removal
  - Quickly conforms to nasal anatomy upon inflation
  - Allows healing in a moist environment with carboxymethylcellulose (CMC) fabric
- anatomiskai pūsitakankis karboksimetri lėliulozė*

Available in a range of sizes:

### Unilateral Models



#### Rapid Rhino 450

- 4.5 cm anterior inflatable
- Successfully implemented in pediatric epistaxis cases
- Also for use in patients with small nasal anatomy



#### Rapid Rhino 551

- 5.5 cm anterior nasal tamponade that contains an airway
- For patients with COPD




#### Rapid Rhino 550

- 5.5 cm anterior tamponade
- Highly effective in stopping 90% of nose bleeds



#### Rapid Rhino 751

- 7.5 cm anterior/posterior device
- For patients with COPD



#### Rapid Rhino 750

- 7.5 cm device is also a successful tool against anterior and many posterior bleeds

*su lateksine pūšle su "druogelio" tipo apsauga*

### Posterior Model

*50.*




#### Rapid Rhino 900 NEW!

- 9 cm inflatable epistaxis device
- Contains 2 balloons to control intranasal hemorrhage originating from the sphenopalatine artery

*2 balionai - galiniam ir pūšikiam karyjavimui stablę*

### Bilateral Models




#### Rapid Rhino 552

- 5.5 cm bilateral device for more severe epistaxis cases
- Stops bleeding vessels in the anterior region within both nostrils




#### Rapid Rhino 752

- 7.5 cm anterior/posterior device features an airway lumen made of medical grade PVC



Rapid Rhino anatomy

## Platelet aggregation and tamponade working together



Rapid Rhino with Carboxymethylcellulose (CMC) combines platelet aggregation with tamponade.

## Available in a range of sizes:

Catalog no.	Epistaxis solution	Size	Units per box
<b>Rapid Rhino Anterior Device</b> Ultra low profile provides for unsurpassed ease of insertion			
RR 450	for small anatomy/pediatrics	4.5 cm	10
RR 550	for anterior adult epistaxis	5.5 cm	10
RR 551	for COPC patients (Airway included)	5.5 cm	10
<b>Rapid Rhino Posterior/Anterior Device</b> Cuff will conform to nasal anatomy and provide gentle tamponade			
RR 750	for challenging posterior and anterior epistaxis	7.5 cm	10
RR 751	for posterior and anterior epistaxis including COPD patients (Airway included)	7.5 cm	10
RR 900	<b>Rapid Rhino Posterior Device</b> Two independent balloons for severe posterior epistaxis	9 cm	10
<b>Rapid Rhino Bilateral Device</b> Offers simultaneous inflation via single Pilot Cuff			
RR 552	for anterior/bilateral epistaxis ideal for counter packing	5.5 cm	10
RR 752	for anterior/posterior epistaxis	7.5 cm	10



### International

**ArthroCare ENT**  
Skeppsbron 2, 111 30 Stockholm,  
Sweden  
Tel: +46 8 546 172 00  
Fax: +46 8 546 172 39  
info@arthrocare.se

### Australia

**ArthroCare (Australasia Pty Ltd)**  
PO Box 2186, Normanhurst NSW 2076  
Australia  
Tel: +61 2 9479 9200 Fax: +61 2 9479 9299  
Email: service@arthrocare.net.au

### Austria

**Atlantech Med. Prod. Vertriebs-AG**  
Kaufmannstrasse 15, 6020 Innsbruck, Austria  
Tel: +43 512 361 538 Fax: +43 512 361 548

### France

**ArthroCare France**  
9, rue Noirrot, 70000 Vesoul, France  
Tel: +33 84 76 81 38 Fax: +33 3 84 76 74 63

### Germany

**ArthroCare (Deutschland) AG**  
Berghauser Str. 62, 42859 Remscheid, Germany  
Tel: +49 2191 93342 0 Fax: +49 2191 93342 49  
info@arthrocare.de

### Scandinavia

**ArthroCare Scandinavia ApS**  
Science Center Skejby  
Brendstrupgårdsvej 21  
8200 Århus N.  
Denmark  
Tel: +45 7023 2850  
Fax: +45 7023 2851  
Email: info@arthrocare.dk

### UK

**ArthroCare UK Ltd**  
St James Business Park  
9, Grimbald Crag Court  
Knaresborough  
North Yorkshire HG5 8QB  
United Kingdom  
Tel: +44 (0) 1423 888806  
Fax: +44 (0) 1423 888816  
Email: info@arthrocareuk.com

### USA

**ArthroCare ENT**  
680 Vaqueros Avenue  
Sunnyvale, CA 94085-3523  
Tel: +1 408-736-0224

Rapid Rhino products are latex free.

**Not for distribution in the United States**

© 2009 ArthroCare Corporation. ArthroCare, Rapid-Pac and Rapid Rhino are registered trademarks of ArthroCare Corporation.

P/N 16674.gb Rev C