

Septoplasty

	Catalog Number
 ENTACT® Septal Stapler 3-Pack	601-00100

Sinus Dressings and Implants

Removables

	Catalog Number
 RAPID RHINO 300 Dressing - Riemann	RR 300
 RAPID RHINO 400 Dressing - Riemann	RR 400
 RAPID RHINO 500 Dressing - Goodman	RR 500
 RAPID RHINO 800 Dressing - Mannheim	RR 800

Dissolvables

	Catalog Number
 RAPID RHINO SINU-KNIT® Dissolvable Dressing	RR 600
 SINU-FOAM® Dissolvable Dressing	RR 650
 NASASTENT Dissolvable Dressing	RR1000

RAPID RHINO® Epistaxis

Unilateral Models

	Catalog Number
 RAPID RHINO 450 Tamponade 4.5cm, inflatable	RR 450
 RAPID RHINO 550 Tamponade 5.5cm, inflatable	RR 550
 RAPID RHINO 551 Tamponade 5.5cm, inflatable with airway	RR 551
 RAPID RHINO 750 Tamponade 7.5cm, inflatable	RR 750
 RAPID RHINO 751 Tamponade 7.5cm, inflatable with airway	RR 751
 RAPID RHINO 900 Tamponade 9cm, inflatable	RR 900
 RAPID RHINO 530 Tamponade - Rapid-Pac®, 5.3cm, non-inflatable	RR 530

Bilateral Models

	Catalog Number
 RAPID RHINO 552 Tamponade 5.5cm, inflatable	RR 552
 RAPID RHINO 555 Tamponade 5.5cm, inflatable with airway	RR 555
 RAPID RHINO 752 Tamponade 7.5cm, inflatable	RR 752
 RAPID RHINO 755 Tamponade 7.5cm, inflatable with airway	RR 755

Learn more at smith-nephew.com

Manufactured by:
ArthroCare Corporation
7000 West William
Cannon Drive
Austin, TX 78735, USA

®Trademark of Smith+Nephew.
All trademarks acknowledged.
©2020 Smith+Nephew. All rights reserved. 18338-EMEA V3 1120

Carboxymethylcellulose (CMC) is a highly absorbent hydrocolloid-polymer polysaccharide material that forms a gel when hydrated.

Anterior *Prüelins Neureybas*

RR 450 4.5cm

RR 550 5.5cm


RR 750 7.5cm

Anterior with Airway

RR 551 5.5cm

RR 751 7.5cm

12 poe



RR 750 Model



The central airway of the **RR 551** & **RR 751** allows the patient some continued nasal airflow with the device properly positioned.

18. p.d. Melingauhs uocies tarupouai

ANTERIOR/ANTERIOR WITH AIRWAY

subanferis durējus balionus

Rapid Rhino 900 contains two balloons to control intranasal hemorrhage originating from the sphenopalatine artery for posterior epistaxis.


13202

Posterior Galiniam kraujavimui

RR 900 9.0cm

Naivovo padėjimas su kuboksinetifikacija - hidrodinaminė veikė

Tampovai kraujavimui iš nosies stabdys, sterilūs



labesini pūslėli → mechan - AIR inflation provides a delicate and lightweight pressure mechanism.

Green-lined swallow guard provides indication for posterior balloon.

POSTERIOR