


„Lunderquist" tipo didelio stangrumo viela



# Lunderquist®

EXTRA-STIFF WIRE GUIDE



## Guide with optimal support.

### Support

A **PTFE-coated, stainless steel mandril** is the key to Lunderquist providing maximum support and low surface friction during endovascular aortic repair and interventional cardiovascular procedures.

### Safety

The **smooth transition** from the stiff mandril to the floppy gold tip (in 260 and 300 cm lengths) and the **flexibility** of the tips aid Lunderquist's atraumatic vessel engagement. A radiopaque gold inner coil allows **precise tip visibility and location** during wire-guide manipulation.

### Compatibility

Lunderquist is compatible with a wide variety of anatomies and needs. **Various lengths and three tip configurations** are available, including a double-curved tip design that mimics the anatomy with a precurved distal tip.

## Cross-Reference Table

| Boston Scientific Meier Wire Guides      |                | Cook Medical Lunderquist Wire Guides |                      |               |           |                        |              |
|------------------------------------------|----------------|--------------------------------------|----------------------|---------------|-----------|------------------------|--------------|
| Product Description                      | Catalog Number | Global Product Number                | Order Number         | Diameter inch | Length cm | Flexible Tip Length cm | Curve        |
| Back-up Meier .035 in, 185 cm Flex J Tip | SCH-30600      | G46729                               | TSMG-35-180-4-LES    | .035          | 180       | 4                      | straight     |
|                                          |                | G49228                               | TSMG-35-180-7-LES    | .035          | 180       | 7                      | straight     |
|                                          |                | G49230                               | TSCMG-35-180-4-3-LES | .035          | 180       | 4                      | single curve |
| Back-up Meier .035 in, 300 cm Ex C Tip   | SCH-30601      | G45352                               | TSMG-35-300-LES      | .035          | 300       | 4                      | straight     |
|                                          |                | G45407                               | TSCMG-35-300-7-LES   | .035          | 300       | 4                      | single curve |
|                                          |                | G45408                               | TSCMG-35-300-LESDC   | .035          | 300       | 4                      | double curve |
| Back-up Meier .035 in, 260 cm w/C Tip    | 30-602         | G31453                               | TSMG-35-260-LES      | .035          | 260       | 4                      | straight     |
|                                          |                | G45353                               | TSCMG-35-260-7-LES   | .035          | 260       | 4                      | single curve |
|                                          |                | G45208                               | TSCMG-35-260-LESDC   | .035          | 260       | 4                      | double curve |



[www.cookmedical.com](http://www.cookmedical.com)

**COOK MEDICAL INCORPORATED**

P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.  
Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

**COOK (CANADA) INC.**

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA  
Phone: 905.640.7110, Toll Free: 800.668.0300

**WILLIAM A. COOK AUSTRALIA PTY. LTD.**

95 Brandl Street, Brisbane Technology Park, Eight Mile Plains  
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

**WILLIAM COOK EUROPE ApS**

Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2010 AI-BOUSM-LUNDMEIR-EN-201005

AORTIC  
INTERVENTION

CRITICAL  
CARE

ENDOSCOPY

INTERVENTIONAL  
RADIOLOGY

LEAD  
MANAGEMENT

PERIPHERAL  
INTERVENTION

SURGERY

UROLOGY

WOMEN'S  
HEALTH