

BoneScalpel Tip Selection Guide



Tip	Part #	Approach	Purpose	Procedures
 <p>10mm Blade, Blunt .5 x 5 mm</p>	MXB-10 Working Length* = 34 mm 15.4	Open	Osteotomy En-Bloc Resection Bone Cutting	Laminoplasty Laminectomy Facetectomy PSO TLIF/PLIF Corpectomy Laminotomy Foraminotomy ACDF
	MXB-10LS Working Length* = 145 mm 15.6	MIS		
	MXB-10LC Working Length* = 145 mm 15.7	MIS		
	MXB-MIS-10 Working Length* = 100 mm 15.12	Mini-Open/MIS		
 <p>20mm Blade, Blunt 1.0 x 6.4 mm</p>	MXB-20 Working Length* = 40 mm 15.3	Open		
	MXB-20LC Working Length* = 145 mm 15.8	MIS		
	MXB-MIS-20LCRS Working Length* = 100 mm 15.11	Mini-Open/MIS		
 <p>25mm Blade, Blunt 1.0 x 6.4 mm</p>	MXB-25 Working Length* = 45 mm 15.5	Open		
 <p>20mm Blade, Serrated 1.0 x 6.4 mm</p>	MXB-B1 Working Length* = 40 mm 15.9	Open		
 <p>Micro Hook Shaver 1.3 x 1.8 mm</p>	MXB-S1 Working Length* = 59 mm 15.1	Open		
	MXB-MIS-S1 Working Length* = 140 mm 15.13	Mini-Open/MIS		
 <p>Macro Hook Shaver 2.2 x 3.0 mm</p>	MXB-S2 Working Length* = 50 mm 15.2	Open		
 <p>Diamond Shaver 4.4 mm</p>	MXB-S3 Working Length* = 52 mm 15.10	Open	360° Bone Shaving & Bone Shaping Decompression Gross Resection	Corpectomy VCR Decortication Laminotomy

*Approximate distance from typical index finger placement to end of probe using BoneScalpel Handpiece (BCM-HP)

If you would like further information or would like to evaluate BoneScalpel, please contact us at +1.631.694.9555



1938 NEW HWY, FARMINGDALE, N.Y. 11735 | +1.631.694.9555 +1.631.694.3285 FAX
MISONIX, INC. | NASDAQ SYMBOL: MSON | MISONIX.COM

Ordering Information

BoneScalpel® System

BCM-SY	<p>Misonix BoneScalpel System</p> <ul style="list-style-type: none"> • Ultrasonic Console • Two Handpieces with Wrenches • Footswitch and System Accessories <p>Console can be configured for 110-130V, 60Hz and 200-240V, 50Hz.</p>
--------	---

BoneScalpel Reusable Items - Autoclavable

BCM-HP	BoneScalpel Universal Handpiece
BCM-SS	Probe Cover, Hard Tissue
BCM-CW	Handpiece Wrench
BCM-2W	Blade T-Wrench
BCM-CBS	Small Cleaning Brush Set
BCM-CBL	Large Cleaning Brush Set
BCM-PR	Replacement Pump Head
E-SYSCART-WS	Misonix System Cart with Shelf

BoneScalpel Tubing – Single Use, Sterile

MXB-T	Irrigation Tube Set 15.14
-------	----------------------------------

BoneScalpel Tips For Hard Tissue Removal – Single Use, Sterile

BLADES	
MXB-10	10 mm Blade, Blunt – includes 10 mm blade, short extension and silicone sleeve
MXB-20	20 mm Blade, Blunt – includes 20 mm blade, short extension and silicone sleeve
MXB-B1	20mm Blade, Serrated – includes 20 mm blade, short extension and silicone sleeve
MXB-25	25 mm Blade, Blunt - includes 25 mm blade, short extension and silicone sleeve
SHIVERS	
MXB-S1	Micro Hook Shaver – includes shaver tip, short extension and silicone sleeve
MXB-S3	Ø4.4 Diamond Shaver – includes shaver tip, short extension and silicone sleeve
MIS	
MXB-10LS	10 mm Blade, Long Straight – includes 10 mm blade, long straight extension and silicone sleeve
MXB-10LC	10 mm Blade, Long Curved – includes 10 mm blade, long curved extension, and silicone sleeve
MXB-20LC	20 mm Blade, Long Curved – includes 20 mm blade, long curved extension, and silicone sleeve
MXB-MIS-20 LCRS	20 mm Blade, Long Curved, Rigid Sheath – includes 20 mm blade, long curved extension, and rigid sheath