




Komet Medical. Powered Instruments.



Rotierende Instrumente

Rotary instruments



p.d.1

- 1.Silomi gržtai techniškai suderinami su LSMU ligoninje Kauno klinikose naudojamu gamintoju „Bien Air“ antgaliu „ChiroPro 980“.
- 2.Gržtas su J formos jungtimi; 3.Pagamintas iš kietmetalio;4.Gržto ilgis 44,5 mm;
- 5.Fisrins galvuts diametras 2,1 mm.



Hartmetall · Tungsten carbide

Hartmetall · Rund



X

Tungsten carbide · Round

volframo karbidas- kietmetalīs

J formos jungtis

X

NC01...	Ø (mm)	# S/FL	Round Shank					J-Notch
			(L = mm) 44,5	(L = mm) 51	(L = mm) 70	(L = mm) 95	(L = mm) 125	
	0,5	6			NC0192.009 STE	NC0192.012 STE		
	0,5	6			NC0192.009	NC0192.012		
	0,8	6			NC0149.009 STE	NC0149.012 STE		
	0,8	6			NC0149.009	NC0149.012		
	1,0	6	NC0101.002 STE	NC0101.004 STE	NC0101.009 STE	NC0101.012 STE		
	1,0	6	NC0101.002	NC0101.004	NC0101.009	NC0101.012		
	1,2	6			NC0193.009 STE	NC0193.012 STE		
	1,2	6			NC0193.009	NC0193.012		
	1,4	6	NC0102.002 STE		NC0102.009 STE			
	1,4	6	NC0102.002		NC0102.009			
	1,6	6		NC0103.004 STE	NC0103.009 STE	NC0103.012 STE		NC0103.014 STE
	1,6	6		NC0103.004	NC0103.009	NC0103.012		NC0103.014
	1,8	6	NC0104.002 STE	NC0104.004 STE	NC0104.009 STE	NC0104.012 STE		NC0104.014 STE
	1,8	6	NC0104.002	NC0104.004	NC0104.009	NC0104.012		NC0104.014
	2,1	6		NC0105.004 STE	NC0105.009 STE	NC0105.012 STE		NC0105.014 STE
	2,1	6		NC0105.004	NC0105.009	NC0105.012		NC0105.014
	2,3	6		NC0106.004 STE	NC0106.009 STE	NC0106.012 STE	NC0106.026 STE	NC0106.014 STE
	2,3	6		NC0106.004	NC0106.009	NC0106.012	NC0106.026	NC0106.014
	2,7	8	NC0128.002 STE		NC0128.009 STE			
	2,7	8	NC0128.002		NC0128.009			
	3,1	6		NC0107.004 STE		NC0107.012 STE		NC0107.014 STE
	3,1	6		NC0107.004		NC0107.012		NC0107.014
	3,1	8	NC0117.002 STE		NC0117.009 STE		NC0117.026 STE	
	3,1	8	NC0117.002		NC0117.009		NC0117.026	
	3,5	8	NC0194.002 STE		NC0194.009 STE	NC0194.012 STE		
	3,5	8	NC0194.002		NC0194.009	NC0194.012		
	4,0	8		NC0108.004 STE	NC0108.009 STE	NC0108.012 STE		NC0108.014 STE
	4,0	8		NC0108.004	NC0108.009	NC0108.012		NC0108.014
	4,0	10	NC0195.002 STE		NC0195.009 STE	NC0195.012 STE	NC0195.026 STE	
	4,0	10	NC0195.002		NC0195.009	NC0195.012	NC0195.026	
	4,5	10	NC0196.002 STE		NC0196.009 STE	NC0196.012 STE		
	4,5	10	NC0196.002		NC0196.009	NC0196.012		
	5,0	8		NC0109.004 STE	NC0109.009 STE	NC0109.012 STE		NC0109.014 STE
	5,0	8		NC0109.004	NC0109.009	NC0109.012		NC0109.014
	5,0	12	NC0155.002 STE		NC0155.009 STE	NC0155.012 STE	NC0155.026 STE	
	5,0	12	NC0155.002		NC0155.009	NC0155.012	NC0155.026	
	6,0	8		NC0110.004 STE	NC0110.009 STE			NC0110.014 STE
	6,0	8		NC0110.004	NC0110.009			NC0110.014
	6,0	14			NC0156.009 STE	NC0156.012 STE	NC0156.026 STE	
	6,0	14			NC0156.009	NC0156.012	NC0156.026	
	6,5	8		NC0170.004 STE	NC0170.009 STE	NC0170.012 STE		
	6,5	8		NC0170.004	NC0170.009	NC0170.012		
	7,0	16	NC0125.002 STE		NC0125.009 STE	NC0125.012 STE	NC0125.026 STE	
	7,0	16	NC0125.002		NC0125.009	NC0125.012	NC0125.026	

p.d.1

p.d.2 sterilus

p.d. 2

1. Silomi gržtai techniškai suderinami su LSMU ligoninje Kauno klinikose naudojamu gamintoju „Bien Air“ antgaliu „ChiroPro 980“.
2. Gržtas, sterilus, su J formos jungtimi; 3. Pagamintas iš kietmetalio ;
4. Gržto ilgis 44,5 mm; 5. Fisrins galvuts diametras 6,0 mm.

inches

1

2

3

4

5

6

151