



YANGZHOU WUYANG MEDICAL EQUIPMENT CO., LTD.

Add: Tongda Road, Touqiao Town, Yangzhou City China
Tel: 0514-87481431 Fax: 0514-87480500
P.C.: 225109 Http://www.wuyangmedical.com

CE Technical Documentation

Cotton Tipped Applicators

Product Name: Cotton Tipped Applicators

Scope: It is mainly used in protecting, disinfecting skin mucous membranes and cleaning wounds.

Product classification: medical equipment/disposable medical supplies

The product is classified as Class I* according to Rule 5 of MDD93/42/EEC Annex IX.

RULE 5	VERDICT
All invasive devices with respect to body orifices, other than surgically invasive devices and which are not intended for connection to an active medical device or which are intended for connection to an active medical device in Class I:	Class I*

Specifications: wooden sticks/PP sticks/bamboo sticks and so on

Product performance and structure:

The product is mainly composed of absorbent cotton ball and wooden cotton bar, its packing is PE or paper-plastic packaging, sterilization treatment, and after various performance testing qualified, supply for use in medical units. The product is used as a sterile product after sterilization and should be destroyed after use.

The product is mainly consists of cotton wool ball and wooden cotton club.

Instructions:

1. A medical staff decides to use a q-tip before using it;
2. Check the package without damage before use and confirm whether it is within the validity period;
3. Open the small package, when using, twist the cotton swab and press the cotton swab on the wound to the wound.
4. The product can only be used at once, and can not be used repeatedly

Disinfection method: ethylene oxide gas disinfection.

YANGZHOU WUYANG MEDICAL EQUIPMENT CO., LTD.
Tongda Road, Touqiao Town, Yangzhou 225109, China



YANGZHOU WUYANG MEDICAL EQUIPMENT CO., LTD.

Add: Tongda Road, Touqiao Town, Yangzhou City China

Tel: 0514-87481431

Fax: 0514-87480500

P.C.: 225109

Http://www.wuyangmedical.com

Use method:

1. The product needs to be operated by a trained physician;
2. Follow the hospital's medical specifications and product instructions when used;
3. This product is only used for one-time use, should be used immediately after opening the package, and should be destroyed after use.

Packing: Mainly paper-plastic bag

Warning:

1. Use this product at one time and do not reuse it;
2. Package damaged, do not use;
3. Pay attention to the sterilization failure period on the packaging bag, which is not allowed to be used.

Storage requirements:

This product should be stored in the relative humidity less than 80%, no corrosive gas, good ventilation, dry, cool clean indoor, avoid extrusion.

Manufacturer:

YANGZHOU WUYANG MEDICAL EQUIPMENT CO., LTD.

Tongda Road, Touqiao Town, Yangzhou 225109, China

Tel.: 86-514-87481431 15152749638 Fax: 86-514-87480500

mail: alex@wuyangmedical.com

Contact Person: Mona Chen

EU Authorizing Representative:

SUNGO Europe B.V.

Fascinatio Boulevard 522, Unit 1.7, 2909VA Capelle aan den IJssel, The Netherlands

Notified Body:

TÜV Rheinland LGA Products GmbH

Tillystraße 2, 90431 Nürnberg, Deutschland

Identification no.: 0197

YANGZHOU WUYANG MEDICAL EQUIPMENT CO., LTD.
Tongda Road, Touqiao Town, Yangzhou 225109, China