


Grow your discovery, fall in love with big savings Save now >



Home > Shop All Products > Gel Electrophoresis Equipment And Supplies > Gel Electrophoresis Reagents And Buffers > Electrophoresis Gel Stains > SYBR™ Green I Nucleic Acid Gel ...



[Certificates](#) [SDS](#)

Invitrogen™

SYBR™ Green I Nucleic Acid Gel Stain, 10,000X concentrate in DMSO

SYBR™ Green I Nucleic Acid Gel Stain is one of the most sensitive stains available for detecting double-stranded DNA

[Read more](#)

Have Questions? [Contact Us](#)

Change view



Catalog Number	Quantity
S7563	500 µL
S7585	20 x 50 µL
S7567	1 mL

3 Options

Catalog number S7563



.e (EUR) / Each

Grow your discovery, fall in love with big savings Save now >

ThermoFisher
SCIENTIFIC



Add to cart

Quantity: 500 µL

[Request bulk or custom format](#)

Product Overview

[Figures](#)

[Recommendations](#)

[Documents](#)

[FAQ](#)

[Citations & References](#)

SYBR™ Green I Nucleic Acid Gel Stain is one of the most sensitive stains available for detecting double-stranded DNA (dsDNA) in agarose and polyacrylamide gels. Because SYBR™ Green I has greater sensitivity for dsDNA, it is especially useful for assays where the presence of contaminating RNA or ssDNA might obscure results. With exceptionally low background fluorescence and spectral characteristics that closely match light sources and filter sets in existing instruments, SYBR™ Green I stain is ideal for use with laser scanners. SYBR™ Green I stain provides:

- **Sensitivity**—at least four times greater than ethidium bromide for DNA in agarose gels
- **Ease of use**—gels soaked in diluted stain can be visualized without desalting
- **Compatibility**—with UV transilluminators, gel documentation systems, and laser scanners
- **Flexibility**—for use with a broad range of applications, including DNA typing, PCR-based assays, DNA damage assays, analysis of complex samples, and real-time PCR detection

Product use

One mL stains 100 minigels. Various pack sizes are available.

For Research Use Only. Not for use in diagnostic procedures.

Specifications


Detection Location	In-Gel Detection
Detection Method	Fluorescence
Label or Dye	SYBR Green I
Product Type	Nucleic Acid Gel Stain
Quantity	500 µL
Shipping Condition	Room Temperature
Target Molecule	DNA
Unit Size	Each

Grow your discovery, fall in love with big savings Save now >




Store at -20°C, protected from light in a desiccator.


Figures



Electrophoretic bandshift assay using SYBR® Green I Nucleic Acid Gel Stain.




DNA molecular weight ladders stained with SYBR® Green I Nucleic Acid Gel



Comparison of single-stranded oligonucleotide detection using SYBR®

Fluorescence spectra

[Open in SpectraViewer](#)



Customers who viewed this item also viewed

Grow your discovery, fall in love with big savings Save now >

ThermoFisher
SCIENTIFIC



SYBR™ Safe DNA Gel Stain

Catalog number: S33102

128,00 / Each

[Add to cart](#)

SYBR™ Gold Nucleic Acid Gel Stain (10,000X...

Catalog number: S11494

282,00 / Each

[Add to cart](#)

SYBR™ Green PCR Master Mix, 1 x 5 mL

Catalog number: 4309

582,00 / Each

[Add to cart](#)

Documents & Downloads

Certificates

[Search](#)

Lot #	Certificate Type	Date	Catalog Number(s)
2994111	Certificate of Analysis	Sep 09, 2024	S7585
2980704	Certificate of Analysis	Aug 12, 2024	S7567
2806848	Certificate of Analysis	Dec 31, 2023	S7585
2720336	Certificate of Analysis	Nov 06, 2023	S7563
2734942	Certificate of Analysis	Sep 25, 2023	S7567

5 results displayed, search above for a specific certificate

[Request a Certificate](#)

Safety Data Sheets



[SDS](#)

Product Information

nuals



Grow your discovery, fall in love with big savings Save now >

ThermoFisher
SCIENTIFIC



[Usage Example: MitoTracker Orange CMTMRos & SYBR Green I \(Japanese\)](#)

Molecular Probes® Handbook

[Nucleic Acid Detection on Gels, Blots and Arrays—Section 8.4](#)

[Nucleic Acid Quantitation in Solution—Section 8.3](#)

[Nucleic Acid Stains—Section 8.1](#)

Frequently asked questions (FAQs)

What are the differences between all the SYBR dyes?

What is the recommended filter for my gel documentation system?

Is SYBR Safe DNA Gel Stain the same as SYBR Green I dye?

Can I stain my TBE gel? How?

Can I use the ethidium bromide filters on my camera to image SYBR dyes?

[View more](#)

Citations & References (438)

Search citations by name, author, journal title or abstract text

[Search](#)

Citations & References

Abstract

[Rapid analysis of the DNA-binding specificities of transcription factors with DNA microarrays.](#)

Authors: Mukherjee S, Berger MF, Jona G, Wang XS, Muzzey D, Snyder M, Young RA, Bulyk ML

Journal: Nat Genet

PubMed ID: 15543148

'We developed a new DNA microarray-based technology, called protein binding microarrays (PBMs), that allows rapid, high-throughput characterization of the in vitro DNA binding-site sequence specificities of transcription factors in a single day. Using PBMs, we identified the DNA binding-site sequence specificities of the yeast transcription factors Abf1, Rap1 and Mig1. ... [More](#)

[Involvement of endonuclease G in nucleosomal DNA quantitation under sustained endogenous oxidative stress.](#)

'We have previously shown that inhibition of catalase and glutathione peroxidase activities by 3-amino-1,2,4-triazole (ATZ) and mercaptosuccinic acid (MS), respectively, in rat primary hepatocytes

Grow your discovery, fall in love with big savings Save now >

ThermoFisher
SCIENTIFIC



[Role of diacylglycerol kinase alpha in the attenuation of receptor signaling.](#) [↗](#)

Authors: Sanjuán MA, Jones DR, Izquierdo M, Mérida I

Journal: J Cell Biol

PubMed ID: 11285286

'Diacylglycerol kinase (DGK) is suggested to attenuate diacylglycerol-induced cell responses through the phosphorylation of this second messenger to phosphatidic acid. Here, we show that DGKalpha, an isoform highly expressed in T lymphocytes, translocates from cytosol to the plasma membrane in response to two different receptors known to elicit T cell ... [More](#)

[Conserved role of a complement-like protein in phagocytosis revealed by dsRNA knockout in cultured cells of the mosquito, Anopheles gambiae.](#) [↗](#)

Authors: Levashina EA, Moita LF, Blandin S, Vriend G, Lagueux M, Kafatos FC

Journal: Cell

PubMed ID: 11257225

'We characterize a novel hemocyte-specific acute phase glycoprotein from the malaria vector, Anopheles gambiae. It shows substantial structural and functional similarities, including the highly conserved thioester motif, to both a central component of mammalian complement system, factor C3, and to a pan-protease inhibitor, alpha2-macroglobulin. Most importantly, this protein serves as ... [More](#)

[Molecular cloning of mouse type 2 and type 3 inositol 1,4,5-trisphosphate receptors and identification of a novel type 2 receptor splice variant.](#) [↗](#)

Authors: Iwai M, Tateishi Y, Hattori M, Mizutani A, Nakamura T, Futatsugi A, Inoue T, Furuichi T, Michikawa T, Mikoshiba K

Journal: J Biol Chem

PubMed ID: 15632133

'We isolated cDNAs encoding type 2 and type 3 inositol 1,4,5-trisphosphate (IP(3)) receptors (IP(3)R2 and IP(3)R3, respectively) from mouse lung and found a novel alternative splicing segment, SI(m2), at 176-208 of IP(3)R2. The long form (IP(3)R2 SI(m2)(+)) was dominant, but the short form (IP(3)R2 SI(m2)(-)) was detected in all tissues ... [More](#)

438 total citations

1 2 3 4 5 ... 88 Next >

Other products to consider

UltraPure™ Low Melting Point Agarose, 100 g

Catalog number: 16520100

692,00 / 100 g

UltraPure™ TAE Buffer, 10X, 1 L

Catalog number: 15558042

31,20 / Each

Grow your discovery, fall in love with big savings Save now >

