



- Features**
- ✓ ✓ • Stainless steel semi sharp cutting edge for soft or weak tissue
 - ✓ ✓ • Size embossed on each curette for easy identification
 - ✓ ✓ • Wide range of sizes
 - ✓ ✓ • Sterile and ready to use
 - ✓ ✓ • Disposable

Used in the following minor skin surgical procedures


- Curettage of a basal cell carcinoma
- Curettage of a verruca
- Curettage of other benign neoplasms of the skin

Clinical areas

- Dermatology
- Surgery
- Podiatry
- Veterinary

Package

- 20 curettes in a dispenser box
- 20 boxes per shipping carton

REF	SIZE	Actual size
MK402	2 mm	
MK403	3 mm	
MK404	4 mm ✓	
MK405	5 mm	
MK407	7 mm ✓	

p.d.10

p.d.10



A dermal curette biopsy technique for various dermatological diseases



A punch (Trepan) biopsy technique for various dermatological diseases

References : Manabu MAEDA, MD & Ph D *Department of Dermatology, Prefectural Gifu Hospital, 4-6-1 Noishiki, Gifu City, 500-8717, Japan*

Product specifications may be changed without notice



kai industries co., ltd.
 1110, Oyana, Seki City, Gifu Pref. 501-3992 Japan
 Phone: +81 (0) 575 28 6600 Fax: +81 (0) 575 28 6611
 Email: info-ind@kai-group.com http://www.kai-ind.co.jp

kai Europe GmbH
 Kottendorfer Str. 5, 42697 Solingen, Germany
 Phone: +49 (0) 212 232380 Fax: +49 (0) 212 2323899

- p.d. 10. Kitretės (diametras 4 mm; 7 mm):
- Sterilios (pažymėta simboliu);
 - Vienkartinės (simbolis ant pakuotės);
 - Laikiklis – plastikinis;
 - Antgalis – iš nerūdijančio medicininio plieno;
 - Ant pakuotės pažymėtas produkto galiojimo laikas;
 - Su numatyta pakuotės atidarymo vieta;
 - Įpakuota po 1 vnt.