



AU

# Instructions For Use

© 2018 Beckman Coulter, Inc. All rights reserved.

## AU WASH SOL. Wash Solution

**REF**

OSR0001 6 x 2 L  
ODR2000 4 x 5 L

### ANNUAL REVIEW

Reviewed by	Date	Reviewed by	Date

### PRINCIPLE

#### INTENDED USE

Cleaning Concentrate.

### REAGENTS

#### WARNING AND PRECAUTIONS

Exercise the normal precautions required for handling all laboratory reagents.

Dispose of all waste material in accordance with local guidelines.

#### GHS HAZARD CLASSIFICATION

Wash Solution

DANGER



H314

Causes severe skin burns and eye damage.

P280

Wear protective gloves, protective clothing and eye/face protection.

P301+P330+P331

IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P303+P361+P353

IF ON SKIN (or hair): Rinse skin with water.

P305+P351+P338

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310

Immediately call a POISON CENTER or doctor/physician.

Sodium Hydroxide 1 - 2%

Genapol X080 1 - 2%

Sulfonic acids, C14-17-sec-alkane, sodium salts 1 - 5%

SDS

Safety Data Sheet is available at [techdocs.beckmancoulter.com](http://techdocs.beckmancoulter.com)

## STORAGE AND STABILITY

The product is stable up to the stated expiry date when stored at 2...25°C.

## ADDITIONAL INFORMATION

### REVISION HISTORY

IFU updated to add Vietnamese language.

Updated Warning and Precautions section

### Preceding version revision history

Added Storage and Stability section.

EC REP

Beckman Coulter Ireland Inc., Lismeehan, O'Callaghan's Mills Co. Clare, Ireland (001) 703-527-3887



Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.