

OPHTHALMIC DIAGNOSTIC STRIPS

In Medical Grade Tyvek Pouches

Schirmer Strips

CE
0434



46

Sterile Ophthalmic strips for Dry Eye test

Packaging:

Each box contains 100 sterile strips.



Indications:

- Used for measurement of tear production
- Gauged for easy reading

Fluorescein Sodium Ophthalmic Strips

CE
0434




47

Composition (per strip):

Each strip is impregnated with 1.0 mg of Fluorescein Sodium

Packaging:

Each box contains 100 sterile strips.



Indications:

- For Applanation Tonometry
- For evaluating Contact Lens fitting

High Molecular Fluorescein Sodium Ophthalmic Strips




Composition (per strip):

Each strip is impregnated with 0.5 mg of high molecular Fluorescein Sodium

Packaging:

Each box contains 100 sterile strips.



Indications:

- For Applanation Tonometry
- For evaluating Contact Lens fitting