

OPHTHALMIC DIAGNOSTIC STRIPS

In Medical Grade Tyvek Pouches

Schirmer Strips



Sterile Ophthalmic strips for Dry Eye test

Packaging:
Each box contains 100 sterile strips.



Indications:


- Used for measurement of tear production
- Gauged for easy reading

Fluorescein Sodium Ophthalmic Strips



Composition (per strip):
Each strip is impregnated with 1.0 mg of Fluorescein Sodium

Packaging:
Each box contains 100 sterile strips.



Indications:


- For Applanation Tonometry
- For evaluating Contact Lens fitting

High Molecular Fluorescein Sodium Ophthalmic Strips



Composition (per strip):
Each strip is impregnated with 0.5 mg of high molecular Fluorescein Sodium

Packaging:
Each box contains 100 sterile strips.



Indications:

- For Applanation Tonometry
- For evaluating Contact Lens fitting