

[Home](#) — MassCheck® Antimycotic Drugs/EXTENDED Plasma Controls

MassCheck® Antimycotic Drugs/EXTENDED Plasma Controls

Order No.: 0252/XT; 0253/XT; 0254/XT

Parameters:

5-Flucytosine, Anidulafungin, Caspofungin, Fluconazole, Hydroxyitraconazole, Isavuconazole, Itraconazole, Ketoconazole, Miconazole, Posaconazole, Voriconazole

MassTox® TDM Series A Antimycotic Drugs/EXTENDED in Serum/Plasma – LC-MS/MS

Product Name	Qty
--------------	-----

Product Name	Qty
MassCheck® Antimycotic Drugs/EXTENDED Plasma Control, Level I Order No.: 0253/XT	0
MassCheck® Antimycotic Drugs/EXTENDED Plasma Control, Level II Order No.: 0254/XT	0
MassCheck® Antimycotic Drugs/EXTENDED Plasma Control, Bi-Level (I+II) Order No.: 0252/XT	0

Quote / Order

Parameters ▼

Description ▲

Technical Data ▼

Downloads ▼

MassCheck® Antimycotic Drugs/EXTENDED Plasma Controls

For MassTox® TDM Series A Parameter Set Antimycotic Drugs/EXTENDED in Serum/Plasma by LC-MS/MS (92922/XT)

0252/XT MassCheck® Antimycotic Drugs/EXTENDED Plasma Control, Bi-Level I + II (lyoph.), 2 x 5 x 1 ml

0253/XT MassCheck® Antimycotic Drugs/EXTENDED Plasma Control, Level I (lyoph.), 5 x 1 ml

0254/XT MassCheck® Antimycotic Drugs/EXTENDED Plasma Control, Level II (lyoph.), 5 x 1 ml

These lyophilised quality controls from Chromsystems are based on human plasma and should be used according to the instruction manual. They are designed to monitor the accuracy and precision of analytical procedures in clinical diagnostics for the quantitative determination of antimycotic drugs in serum/plasma. While following the instructions the controls should be handled and measured in the same manner as a patient specimen.