

[Home](#) — 3PLUS1® Multilevel Plasma Calibrator Set Antimycotic Drugs/EXTENDED



3PLUS1® Multilevel Plasma Calibrator Set Antimycotic Drugs/EXTENDED

Order No.: 92051/XT

Parameters:

5-Flucytosine, Anidulafungin, Caspofungin, Fluconazole, Hydroxyitraconazole, Isavuconazole, Itraconazole, Ketoconazole, Micafungin, Posaconazole, Voriconazole

MassTox® TDM Series A Antimycotic Drugs/EXTENDED in Serum/Plasma – LC-MS/MS

Quantity

Quote / Order

Parameters



Description



Technical Data



Downloads



3PLUS1® Multilevel Plasma Calibrator Set Antimycotic Drugs/EXTENDED

For MassTox® TDM Series A Parameter Set Antimycotic Drugs/EXTENDED in Serum/Plasma by LC-MS/MS (92922/XT)

These lyophilised calibrators from Chromsystems are based on human plasma and should be used according to the instruction manual. They are designed to calibrate the measurement method in clinical diagnostics for the quantitative determination of antimycotic drugs in serum/plasma. While following the instructions the calibrators should be handled and measured in the same manner as a patient specimen. These calibrators are only available as multilevel calibrator set with 3 different concentration levels plus 1 blank value.