

Subject: R: Question about HHV6 Real-TM Quant kit
Date: 2023 January 4, Wednesday 14:00:20 Eastern European Standard Time
From: Sacace Info
To: Aurelija Bilkauskienė
CC: 'Specialists Sacace'

Dear Aurelija,

Simone and Mariarosaria are in holidays but I want to inform you that you can use also serum or plasma for HHV6 detection but the sensibility will be lower.

Best regards,
Dr. Alexander Safoian
Sacace Biotechnologies s.r.l.
via Scalabrini, 44
22100 Como, Italy
Tel. [+390314892927](tel:+390314892927)
Fax [+390314892926](tel:+390314892926)
www.sacace.com

Da: Aurelija Bilkauskienė [mailto:aurelija@interautomatika.com]
Inviato: giovedì 29 dicembre 2022 15:23
A: sacace.ceo@gmail.com
Cc: Sacace Info <info@sacace.com>
Oggetto: Question about HHV6 Real-TM Quant kit
Priorità: Alta

Dear all,

Since Simone is out of office until January 9th and we have a public tender deadline on January 6th, maybe you could help answering this question?

In the IFU of HHV6 Real-TM Quant kit (ref. V10-100FRT) it is indicated that this kit can analyze DNA extracted from whole blood but serum and plasma samples are not mentioned. Could you please confirm that this kit can analyze DNA extracted from serum and plasma samples?

Pagarbiai / Best Regards,

Aurelija Bilkauskienė
Produktų specialistė | Product specialist
Mob. +370 699 33418
El. p. aurelija@interautomatika.com

INTERAUTOMATIKA, UAB
Tel. +370 5 2607810
Faks. +370 5 2411464
Spaudos g. 6, 05131, Vilnius, Lietuva
interautomatika.eu
interlab.lt