



No.DOC-DQ010- 0245/DOC- DC-0140969
Rev.24

DECLARATION OF CONFORMITY

We, **TERUMO CORPORATION**

44-1, 2-chome, Hatagaya, Shibuya-ku, Tokyo 151-0072, Japan

being the manufacturer of:

RADIFOCUS Introducer II

Product : Catheter Introducer

declare that the above products of **Class IIa** are in conformity with the provisions of the EC Council Directive 93/42/EEC of 14 June 1993, as amended, concerning medical devices, and have been subject to the conformity assessment procedure laid down in Article 11, 2 and 11, 3(a) of the Directive, relating to the "Full quality assurance" set out in Annex II, and by certification of Annex II, excluding Section 4 under the supervision of TÜV Rheinland LGA Products GmbH (Registration No.: HD 60145252 0001), Tillystraße 2, 90431 Nürnberg Germany, as Notified Body authorized by the German Competent Authority and carrying the Notified Body No. 0197.

Authorized European Representative :

TERUMO EUROPE N.V.

Interleuvenlaan 40, 3001 Leuven, Belgium

Object of the declaration: see Appendix A

Tokyo, 2024-07-25

(place and date of issue)

署名者:
hiroki sasagawa
署名者名: hiroki sasagawa
署名理由: この文書を承認する
署名時刻: 2024-07-25 | 4:42:35 午後 JST
3F7B1A275BA74182B04CDC2E69B7202D

Hiroki Sasagawa
General Manager
Quality Assurance Department
TERUMO CORPORATION



Appendix A - List of Code Number Structure

□ □ * □ □ □ □ □ □ □ □ □ □
1 2 3 4 5 6 7 8 9 10 11 12

Character number	Characters & Meaning
1, 2	Product name RS: Introducer kit RM: Introducer kit containing hydrophilic polymer-coated sheath.
3	Destination *: for export
4	Kit contents A : Sheath, Dilator, Mini guide wire, Syringe, Entry needle (SurfloFlash), Scalpel* ¹ , (Guide inserter* ²) B : Sheath, Dilator, Mini guide wire, (Guide inserter* ²) C : Sheath, Dilator E : Sheath, Dilator, Mini guide wire, Entry needle (SurfloFlash), (Guide inserter* ²) G : Sheath, Dilator, Mini guide wire, Scalpel* ¹ (Guide inserter* ²) H : Dilator J : Sheath, Dilator, Scalpel* ¹ K : Sheath, Dilator, Mini guide wire, Syringe, Entry needle (metallic needle), Scalpel* ¹ , (Guide inserter* ²) L : Sheath, Dilator, Syringe, Entry needle (SurfloFlash), Scalpel* ¹ M : Sheath, Dilator, Entry needle (SurfloFlash), Scalpel* ¹ N : Sheath, Dilator, Mini guide wire, Syringe, Scalpel* ¹ , (Guide inserter* ²) P : Sheath, Dilator, Mini guide wire, Syringe, Entry needle (SurfloFlash), (Guide inserter* ²) Q : Dilator, Mini guide wire, (Guide inserter* ²) R : Sheath, Dilator, Mini guide wire, Entry needle (metallic needle / Metallic Entry Needle improved product), (Guide inserter* ²) S : Sheath, Dilator, Mini guide wire, Entry needle (SurfloFlash), Scalpel* ¹ , (Guide inserter* ²) W : Mini guide wire * ¹ : not contained in the export specifications * ² : contained when the mini guide wire has an angled tip or a J tip.



Character number	Characters & Meaning																																																																																																														
5-6	Sheath Size (w/o hydrophilic polymer coating) Characters: 40 50 55 60 65 70 75 80 85 90 10 11 00 Size: 4.0 5.0 5.5 6.0 6.5 7.0 7.5 8.0 8.5 9.0 10.0 11.0 no sheath Sheath Size (with hydrophilic polymer coating) Characters: F4 F5 F6 F7 Size: 4.0 5.0 6.0 7.0																																																																																																														
7	Mini guide wire OD, Dilator ID, Size of Entry needle (length of projecting portion of dilator is 25mm) Standard Type (the items with their product code starting with RS*) <table border="1" data-bbox="481 676 1449 1066"> <thead> <tr> <th colspan="2" rowspan="2"></th> <th colspan="5">Mini guide wire diameter/ Dilator inner diameter</th> <th rowspan="2">Type of Surflo</th> </tr> <tr> <th>0.018"</th> <th>0.021"</th> <th>0.025"</th> <th>0.035"</th> <th>0.038"</th> </tr> </thead> <tbody> <tr> <td rowspan="4">Type of scalpel and tray</td> <td>b-type</td> <td>A</td> <td>D</td> <td>G</td> <td>K</td> <td>N</td> <td rowspan="2">Standard</td> </tr> <tr> <td>a-type</td> <td>B</td> <td>E</td> <td>H</td> <td>L</td> <td>P</td> </tr> <tr> <td>b-type</td> <td>C</td> <td>F</td> <td>J</td> <td>M</td> <td>Q</td> <td rowspan="2">With adapter</td> </tr> <tr> <td>a-type</td> <td>V</td> <td>W</td> <td>X</td> <td>Y</td> <td>Z</td> </tr> <tr> <td colspan="2">Entry needle size</td> <td>22G×1"</td> <td>22G×1"</td> <td>20G×2"</td> <td>18G×2 1/2"</td> <td>16G×2 1/2"</td> <td rowspan="2"></td> </tr> <tr> <td colspan="2">metallic needle size</td> <td>---</td> <td>---</td> <td>---</td> <td>18G×2 3/4"</td> <td>18G×2 3/4"</td> </tr> </tbody> </table> <p>*Kit containing a hydrophilic polymer-coated sheath (the items with their product code starting with RM*):</p> <table border="1" data-bbox="379 1133 1449 1635"> <thead> <tr> <th colspan="2" rowspan="2"></th> <th colspan="5">Mini guide wire diameter/ Dilator inner diameter</th> <th rowspan="2">Type of Surflo</th> </tr> <tr> <th>0.018"</th> <th>0.021"</th> <th>0.025"</th> <th>0.035"</th> <th>0.038"</th> </tr> </thead> <tbody> <tr> <td rowspan="4">Type of scalpel and tray</td> <td>b-type</td> <td>A</td> <td>D</td> <td>G</td> <td>K</td> <td>N</td> <td rowspan="2">With adapter</td> </tr> <tr> <td>a-type</td> <td>B</td> <td>E</td> <td>H</td> <td>L</td> <td>P</td> </tr> <tr> <td>b-type</td> <td>C</td> <td>F</td> <td>J</td> <td>M</td> <td>Q</td> <td rowspan="2">Standard</td> </tr> <tr> <td>a-type</td> <td>V</td> <td>W</td> <td>X</td> <td>Y</td> <td>Z</td> </tr> <tr> <td colspan="2">Entry needle size</td> <td>22G×1"</td> <td>22G×1"</td> <td>20G×2"</td> <td>18G×2 1/2"</td> <td>16G×2 1/2"</td> <td rowspan="2"></td> </tr> <tr> <td colspan="2">metallic needle size</td> <td>---</td> <td>Metallic Entry Needle improved product 21G×1 2/5"</td> <td>Metallic Entry Needle improved product 20G×1 2/5"</td> <td>18G×2 3/4"</td> <td>18G×2 3/4"</td> </tr> </tbody> </table>			Mini guide wire diameter/ Dilator inner diameter					Type of Surflo	0.018"	0.021"	0.025"	0.035"	0.038"	Type of scalpel and tray	b-type	A	D	G	K	N	Standard	a-type	B	E	H	L	P	b-type	C	F	J	M	Q	With adapter	a-type	V	W	X	Y	Z	Entry needle size		22G×1"	22G×1"	20G×2"	18G×2 1/2"	16G×2 1/2"		metallic needle size		---	---	---	18G×2 3/4"	18G×2 3/4"			Mini guide wire diameter/ Dilator inner diameter					Type of Surflo	0.018"	0.021"	0.025"	0.035"	0.038"	Type of scalpel and tray	b-type	A	D	G	K	N	With adapter	a-type	B	E	H	L	P	b-type	C	F	J	M	Q	Standard	a-type	V	W	X	Y	Z	Entry needle size		22G×1"	22G×1"	20G×2"	18G×2 1/2"	16G×2 1/2"		metallic needle size		---	Metallic Entry Needle improved product 21G×1 2/5"	Metallic Entry Needle improved product 20G×1 2/5"	18G×2 3/4"	18G×2 3/4"
				Mini guide wire diameter/ Dilator inner diameter						Type of Surflo																																																																																																					
		0.018"	0.021"	0.025"	0.035"	0.038"																																																																																																									
Type of scalpel and tray	b-type	A	D	G	K	N	Standard																																																																																																								
	a-type	B	E	H	L	P																																																																																																									
	b-type	C	F	J	M	Q	With adapter																																																																																																								
	a-type	V	W	X	Y	Z																																																																																																									
Entry needle size		22G×1"	22G×1"	20G×2"	18G×2 1/2"	16G×2 1/2"																																																																																																									
metallic needle size		---	---	---	18G×2 3/4"	18G×2 3/4"																																																																																																									
		Mini guide wire diameter/ Dilator inner diameter					Type of Surflo																																																																																																								
		0.018"	0.021"	0.025"	0.035"	0.038"																																																																																																									
Type of scalpel and tray	b-type	A	D	G	K	N	With adapter																																																																																																								
	a-type	B	E	H	L	P																																																																																																									
	b-type	C	F	J	M	Q	Standard																																																																																																								
	a-type	V	W	X	Y	Z																																																																																																									
Entry needle size		22G×1"	22G×1"	20G×2"	18G×2 1/2"	16G×2 1/2"																																																																																																									
metallic needle size		---	Metallic Entry Needle improved product 21G×1 2/5"	Metallic Entry Needle improved product 20G×1 2/5"	18G×2 3/4"	18G×2 3/4"																																																																																																									



No.DOC-DQ010- 0245/DOC- DC-0140969

Rev.24

Character number	Characters & Meaning
8-9	Length of sheath 05 ~: 50 mm ~
10	Type of mini guide wire A: Plastic, Angled B: Plastic, Angled, 4mm of tip cut E: Straight E F: Flex, Angled H: Short tapered, short angle J: Plastic, 3mm J K: Stiff, Angled M: Spring, J N: No mini guide wire contained. P: Spring, Straight S: Plastic, Straight V: Plastic, 1.5mm J-angle Y: Flex, Straight
11	Packaging Q: Tray package (Multi-language) R: Pouch package (Multi-language)
12	Reserved 1: With scalpel. 5: Inner diameter of dilator at distal end: 0.038", Outer diameter of mini guide wire: 0.035" Z: Entry needle: 20Gx2"→20Gx1 1/4" Length of mini guide wire: 80cm→45cm, scalpel contained. W: Entry needle: 20Gx2"→20Gx1 1/4", scalpel contained.



Appendix A - List of Code Number Structure

□ □ □ □ □ □ □ □ □ □ □ □
1 2 3 4 5 6 7 8 9 10 11 12

Character number	Characters & Meaning																																	
1, 2	Product name RS: Introducer kit RM: Introducer kit containing hydrophilic polymer-coated sheath																																	
3	Destination + / *: Manufactured by TVC for worldwide excluding Japan																																	
4	Kit contents A : Sheath, Dilator, Mini guide wire, Syringe, Entry needle (SurfloFlash), Scalpel* ¹ , (Guide inserter* ²) B : Sheath, Dilator, Mini guide wire, (Guide inserter* ²) C : Sheath, Dilator R : Sheath, Dilator, Mini guide wire, Entry needle (metallic needle / Metallic Entry Needle improved product), (Guide inserter* ²) *1: not contained in the export specification *2: contained when the mini guide wire has an angled tip or a J tip.																																	
5-6	Sheath Size (Fr) (w/o hydrophilic polymer coating) <table border="1"> <tr> <td>Characters</td> <td>40</td> <td>50</td> <td>60</td> <td>70</td> <td>80</td> </tr> <tr> <td>Size</td> <td>4.0</td> <td>5.0</td> <td>6.0</td> <td>7.0</td> <td>8.0</td> </tr> </table> Sheath Size (Fr) (with hydrophilic polymer coating) <table border="1"> <tr> <td>Characters</td> <td>F5</td> <td>F6</td> </tr> <tr> <td>Size</td> <td>5.0</td> <td>6.0</td> </tr> </table>	Characters	40	50	60	70	80	Size	4.0	5.0	6.0	7.0	8.0	Characters	F5	F6	Size	5.0	6.0															
Characters	40	50	60	70	80																													
Size	4.0	5.0	6.0	7.0	8.0																													
Characters	F5	F6																																
Size	5.0	6.0																																
7	Mini guide wire OD, Dilator ID, Size of Enrty needle (length of projecting portion of dilator is 25mm) Standard Type (the items with their product code starting with RS*) <table border="1"> <thead> <tr> <th rowspan="2"></th> <th colspan="5">Mini guide wire diameter/ Dilator inner diameter</th> <th rowspan="2">Type of Surflo</th> </tr> <tr> <th>0.018"</th> <th>0.021"</th> <th>0.025"</th> <th>0.035"</th> <th>0.038"</th> </tr> </thead> <tbody> <tr> <td>Type of scalpel and tray</td> <td>A</td> <td>D</td> <td>G</td> <td>K</td> <td>N</td> <td>Standard</td> </tr> <tr> <td>Entry needle size</td> <td>---</td> <td>---</td> <td>20G×2"</td> <td>18G×2 1/2"</td> <td>16G×2 1/2"</td> <td></td> </tr> <tr> <td>Metallic needle size</td> <td>Metallic Entry Needle improved product 22G×1 2/5"</td> <td>Metallic Entry Needle improved product 21G×1 2/5"</td> <td>Metallic Entry Needle improved product 20G×1 2/5"</td> <td>18G×2 3/4"</td> <td>18G×2 3/4"</td> <td></td> </tr> </tbody> </table>		Mini guide wire diameter/ Dilator inner diameter					Type of Surflo	0.018"	0.021"	0.025"	0.035"	0.038"	Type of scalpel and tray	A	D	G	K	N	Standard	Entry needle size	---	---	20G×2"	18G×2 1/2"	16G×2 1/2"		Metallic needle size	Metallic Entry Needle improved product 22G×1 2/5"	Metallic Entry Needle improved product 21G×1 2/5"	Metallic Entry Needle improved product 20G×1 2/5"	18G×2 3/4"	18G×2 3/4"	
	Mini guide wire diameter/ Dilator inner diameter					Type of Surflo																												
	0.018"	0.021"	0.025"	0.035"	0.038"																													
Type of scalpel and tray	A	D	G	K	N	Standard																												
Entry needle size	---	---	20G×2"	18G×2 1/2"	16G×2 1/2"																													
Metallic needle size	Metallic Entry Needle improved product 22G×1 2/5"	Metallic Entry Needle improved product 21G×1 2/5"	Metallic Entry Needle improved product 20G×1 2/5"	18G×2 3/4"	18G×2 3/4"																													



No.DOC-DQ010- 0245/DOC- DC-0140969

Rev.24

Kit containing a hydrophilic polymer-coated sheath (the items with their product code starting with RM*)				
	Mini guide wire diameter/ Dilator inner diameter			Type of Surflo
	0.021"	0.025"	0.035"	
Type of scalpel and tray	---	G	---	With adapter
	F	J	M	Standard
Entry needle size	---	20G×2"	---	
Metallic needle size	Metallic Entry Needle improved product 21G×1 2/5"	Metallic Entry Needle improved product 20G×1 2/5"	---	
8-9	Length of sheath 00 : no sheath 05~ : 50mm~ ~			
10	Type of mini guide wire A: Plastic, Angled M: Spring, J angled N: No mini guide wire contained. P: Spring, Straight S: Plastic, Straight H: Short tapered, short angle			
11	Packaging Q: Tray package (Multi-language #, Chinese) R: Pouch package (Multi-language #, Chinese)			
12	Reserved 1: With scalpel. 5: Inner diameter of dilator at distal end: 0.038", Outer diameter of mini guide wire: 0.035" Z: Entry needle: 20Gx2"→20Gx1 1/4" Length of mini guide wire: 80cm→45cm, scalpel contained. W: Entry needle: 20Gx2" →20Gx1 1/4" , scalpel contained.			