

To whom it may concern

Statement overheat protection

01 June 2023

Dear Sir / Madam,

With this statement we confirm:

Temperature sensors are placed in the lid seal, in the centrifuge chamber, on the condenser and at the suction pipe; the latter in order to prevent overheating of the motor and electronic parts. This is confirmed by the CB report and described in more detail in the service manual.

Due to changing IVD regulations the technical documentation is in process.

Yours faithfully
Andreas Hettich GmbH & Co. KG

**Andreas Hettich
GmbH & Co. KG**
Föhrenstr. 12, 78532 Tuttlingen
Postfach 2 60, 78503 Tuttlingen
Telefon (07461) 705-0
Telefax (07461) 705-1125



Dennis Maurer
Senior Area Sales Manager Europe