




Other packing units upon request!

**Original Herenz Steristoppers®**

Item	No.	for glasses with inner-Ø	PU = MOQ with appr.
1012205	22 k	21.5–23.0 mm	500 pcs/carton
1012208	22 l	22.0–24.0 mm	2,000 pcs/carton
1012350	23,5 P	25.0–27.0 mm	2,000 pcs/carton
1012600	26	26.0–27.0 mm	2,000 pcs/carton
1012700	27	27.0–32.5 mm	1,000 pcs/carton
1012705	27 k	26,0–26,5 mm	2,500 pcs/carton
1012900	29 ****	29.5–31.0 mm	850 pcs/carton
1012919	29 A (< No.29)	29.0–30.0 mm	1,500 pcs/carton
1012905	29 k	26.5–27.5 mm	2,500 pcs/carton
1013200	32	32.5–35.0 mm	800 pcs/carton
1013201	32 P	32.0–33.0 mm	750 pcs/carton
1013225	32 PD	30.0–32.5 mm	750 pcs/carton
1013400	34	33.5–36.5 mm	750 pcs/carton
1013600	36	35.5–39.5 mm	600 pcs/carton
1013700	37	37.5–41.5 mm	600 pcs/carton
1013800	38	40.0–43.5 mm	500 pcs/carton
1014001	40 P	42.0–46.5 mm	500 pcs/carton
1014501	45 P	47.0–48.0 mm	300 pcs/carton
1014504	45 PDL	47.0–48.0 mm	300 pcs/carton
1014800	48 *****	47.0–49.0 mm	400 pcs/carton
1015003	50 Z	47.0–49.0 mm	200 pcs/carton
1016001	60 P	57.0–60.0 mm	200 pcs/carton
1017005	70 k	69.0–71.0 mm	100 pcs/carton



- \* = for thick-walled test tubes
- \*\* = for thin-walled test tubes
- \*\*\* = for small Erlenmeyerflasks
- \*\*\*\* = for big Erlenmeyerflasks
- \*\*\*\*\* = for Erlenmeyerflasks wide neck

**Original Herenz Steristoppers® Magnum**

For glasses with high tolerances, 2 executions.

Item	Description	for glasses w/ inner-Ø of neck	PU = MOQ with appr.
1012721	Magnum	26.0–30.0 mm	800 pcs/carton
1012723	MA 2/3 short	25.5–29.0 mm	250 pcs/carton

**Original Herenz Steristopper®-caps**

Small and big Steristopper®-Caps for glasses without rim.

Item	outer-Ø mm	PU = MOQ with appr.
1011411	appr. 14.0	5,000 pcs/carton

**Original Steristoppers® with cotton-swab**

ABS ←

For sealing of test tubes (available with cork-stoppers upon request).

Item	Description	PU = MOQ
1020006	Steristoppers® with cotton-swab	100 pcs/bag