

05 June 2020



Fisher & Paykel Healthcare Limited  
Paykel Building  
15 Maurice Paykel Place, East Tamaki  
P O Box 14 348, Panmure  
Auckland, New Zealand  
Telephone: +64 9 574 0100  
Facsimile: +64 9 574 0158  
Website: [www.fphcare.com](http://www.fphcare.com)

### Material Composition of F&P AirSpiral™+ Circuit Kits

To Whom It May Concern:

Please find below material component and composition information for F&P AirSpiral Circuit Kits.

<b>Item code</b>	900PT561 (including MR290 chamber)
<b>Materials</b>	Ethylene-octene copolymer (POE), insulated copper, polypropylene (PP), thermoplastic elastomer (TPE), acetal copolymer (POM), liquid silicone rubber (LSR), bronze, thermoplastic vulcanizate (TPV), acrylonitrile-butadiene-styrene copolymer (ABS), polyethylene (PE), thermoplastic elastomer (TPE), aluminum, ink, acrylate adhesive, polycarbonate (PC), polyvinyl chloride (PVC)
<b>Materials not present</b>	Not made with natural rubber latex, or phthalates (DEHP, DBP, BBP)

Please contact your local Fisher & Paykel Healthcare sales representative if you have any further questions regarding these products.

Yours sincerely,

**Emma Beacham**  
Product Manager  
AIRVO / Optiflow Marketing