


Larident Clin Grip



1 pirkimo dalis

Description

For quick temporary storage of root canal instruments. Can be filled with disinfectant solution, autoclavable, made of aluminum.



M+W

M+W Select foam inserts

variant *

Height 0.5 cm, Ø 5 cm

Package contents

Foam inserts, height 0.5 cm, diameter 5 cm, suitable for Clin Grip and Interim Stand, 50 pieces.



Ionoseal

2 pirkimo dalis



One component light-cured resin modified glass ionomer liner

Indications

Lining under all types of filling material

Extended fissure sealing

Restoration of small lesions

Advantages

- Saves time: Ready-to-use one component material
- Fast 20 sec. light-cure (2 mm)
- High compressive strength (226 MPa)
- Acid resistant allows for total-etch bonding technique
- Patented NDT^{®1} Syringe saves material
- High fluoride release protects against secondary caries
- High biocompatibility
- Radiopaque

Presentation

REF 1326 pack 3 × 2.5 g syringes, intraoral tips type 41

REF 2146 intraoral tips type 41 - 100 pack

7 Pirkimo dalis

SDI

riva coat

[Product Information](#)

[Instructions](#)

[Reorder Codes](#)

[Download](#)

Riva Coat is designed for the prevention of saliva and water contamination to exposed surfaces of glass ionomer cements during the first 24 hours immediately after placement.



8 Pirkimo dalis



SIMPLEE®

All you need. No more. But no less.


[Simplee Products for Praxis](#) ▾

[Distributors](#)

[Download](#)

[Gallery](#)

[Warehouse Info](#)



SIMPLEE FLOSS

Waxed, menthol flavoured in a handy dispenser box.

Simple Floss, 40 m

156 820

9 Pirkimo dalis

Essentia - the aesthetic restorative with duo-layering and a reduced shade inventory



Syringe Refill

- 2 mL syringe refill with touch cap: LD; MD; DD; LE; DE; U; OM
- 2 mL syringe with 20 Dispensing Tips: ML
- 2.4 mL syringe: RBM; BM; WM

Unitip Refill

- Refills of 15 unitips: LD; MD; DD; LE; DE; U

Starter Kit Syringe

- 7 syringes, 1 of each shade: LD; MD; DD; LE; DE; U; ML
- 20 Dispensing Tips
- Accessories

Starter Kit Unitip

- 30 unitips (5 unitips each in 6 shades: LD; MD; DD; LE; DE; U)
- 1 syringe ML; 20 Dispensing Tips
- Accessories

Modifier Kit


- 4 Syringes: RBM, BM, WM, OM



Essentia Shades:



Essentia Modifiers:



OptiDisc™ General Assorted Kit

To buy this product, please create an online account/login.

SKU#: 4200

♥ ADD TO WISH LIST

Made of durable polyester impregnated with aluminium oxide particles, the four grit system achieves a long lasting result. A colour-coded system enable immediate identification of the grit type. The flush-mounted mandrel ensures correct positioning of the disc and protects against scarring of the restoration. OptiDisc also features the choice of a short size mandrel to facilitate access to posterior and other restricted areas.

Content:

120 x ø 9.6 mm Assorted OptiDisc [30 x Extra-coarse, 30 x Coarse, 30 x Fine, 30 x Extra-fine],

120 x ø 12.6 mm Assorted OptiDisc [30 x Extra-coarse, 30 x Coarse, 30 x Fine, 30 x Extra-fine],

3 x Standard Mandrels, 2 x Short Mandrels, 1 x OptiShine



11.1 Pirkimo dalis



SAFE

OptiDisc's active side is easy to recognise and the patented mandrel locks below the surface of the disc, protecting your restoration against scarring and scratching. It also guarantees optimal torque transmission, so the disc does not slide and is not rpm sensitive, giving you better and safer handling.



4 COLOURS 4 STAGES OF ABRASIVENESS

The OptiDisc system offers the choice of four grit level discs for contouring, finishing and high gloss polishing. It is suitable for composite, glass ionomers, amalgams, semiprecious and precious metals. Discs are also available in three different diameter sizes depending on your needs: 9.6mm, 12.6mm and 15.9mm.



The colour-coded discs help you identify the grit size you are looking for easily and clearly.



11.2-5 Pirkimo dalys


SAFE

OptiDisc's active side is easy to recognise and the patented mandrel locks below the surface of the disc, protecting your restoration against scarring and scratching. It also guarantees optimal torque transmission, so the disc does not slide and is not rpm sensitive, giving you better and safer handling.



4 COLOURS 4 STAGES OF ABRASIVENESS

The OptiDisc system offers the choice of four grit level discs for contouring, finishing and high gloss polishing. It is suitable for composite, glass ionomers, amalgams, semiprecious and precious metals. Discs are also available in three different diameter sizes depending on your needs: 9.6mm, 12.6mm and 15.9mm.



Refill packs

Ø 9.6mm

4197 OptiDisc Extra-Coarse 9.6mm (80 pcs)

4181 OptiDisc Coarse/Medium 9.6mm (100 pcs)

4182 OptiDisc Fine 9.6mm (100 pcs)

4183 OptiDisc Extra-Fine 9.6mm (100 pcs)

Ø 12.6mm

4198 OptiDisc Extra-Coarse 12.6mm (80 pcs)

4184 OptiDisc Coarse/Medium 12.6mm (100 pcs)

4185 OptiDisc Fine 12.6mm (100 pcs)

4186 OptiDisc Extra-Fine 12.6mm (100 pcs)

Ø 15.9mm

4199 OptiDisc Extra-Coarse 15.9mm (80 pcs)

4191 OptiDisc Coarse/Medium 15.9mm (100 pcs)

4192 OptiDisc Fine 15.9mm (100 pcs)

4193 OptiDisc Extra-Fine 15.9mm (100 pcs)

SIMPLEE UNIVERSAL

surišiklis

5 ml,

- Visoms ėsdinimo technikoms
- Hidrofiliškas
- Be cheminio aktyvatoriaus suderinamas su visomis cheminio ir dvigubo kietėjimo bei šviesoje kietėjančiomis medžiagomis
- Su MDP monomeru
- Bepalvis
- Be bisfenolio A
- Tinka metalo ir keramikos paviršių paruošimui



Pirkimo dalis 12

Product description

Simplee Universal Bond is a single bottle self-etch light-curing adhesive system (Seventh Generation) responsible for the union between tooth structure (enamel and dentin) and restorative materials. The primer and bonding agent are combined in a single bottle. Prior use of acid etchant is optional on both enamel and dentin, it can be used, therefore, in three protocols: self-etch, selective-etch and total-etch.

It can be used for direct restorations made of composite resins and light-curing glass ionomer or for indirect procedures when associated with resin cement. The product contains MDP, which allows it to be used as a metallic and ceramic primer.

Simplee Universal Bond is made with a combination of different photoinitiators. These react with each other and enhance the curing effect of the light from the polymerisation lamps. This in turn increases the bond strength and improves the mechanical properties of the adhesive film. Simplee Universal Bond is colourless and thus contributes to high aesthetics.

Product Presentation

Package containing a bottle with 5ml of the product.

Basic Composition

Active Ingredients: MDP (10-methacryloyloxydecyl dihydrogen phosphate), methacrylate monomers, photoinitiators complex, co-initiators, stabilizers.

Inactive Ingredients: inert load (silica particles) and vehicle (ethanol).

LM

feel the
difference

Pirkimo dalis 15



Columbia 4L-4R SD

- With sharpen free LM Sharp Diamond™ coated tips

LM ErgoSense®

219-220SDES

OptraGate®

Moisture control made easy

OptraGate® is a latex-free lip and cheek retractor. The aid's unique flexibility and elasticity make it comfortable to wear and assist patients to keep their mouth open.

OptraGate
for popular relative isolation

16-17 pirkimo dalys

OptraDam Plus
for absolute isolation



Photos: courtesy of Ivoclar Vivadent AG, Liechtenstein

ivoclar | 100 years

- OptraGate White -
 - Assortment (40 x regular, 40 x small)
 - Refills x 80 (regular, small or junior)
 - Refills x 40 (regular, small or junior)

ProTaper Universal Hand



Files

Pirkimo dalys 18,19

Product Highlights

- Nickel titanium construction
- ISO-colored handles
- Sold as pack of 6



Taper

Assorted .02v-.09v

.07v

.08v

.09v

.06v

.04v

.02v

HAND INSTRUMENT STERILIZATION:


ProFile®, ProFile® Series 29®, ProFile® GT®, ProTaper®, Lexicon® Hand Instruments

Endo Instrument Sterilization:

Files must be sterilized before use. ANSI/ADA Specification 28 recommends:

- Scrub the instruments with soap and warm water.
- Rinse thoroughly with distilled or deionized water.
- Allow to air dry.
- Place the instruments, unwrapped, in the autoclave tray.
- Use fresh distilled or deionized water.
- Steam Autoclave at 136° C (plus or minus 2° C) for 20 minutes.

ROUND EXCAVATORS



Cariology:

- For the removal of softened, carious dentin
- For the removal of temporary fillings

Surgery:

- To sever the periodontal ligament prior to extraction
- For the removal of granulation tissue after extraction



Pirkimo dalis 20

LMErgoSense®

ES = LM-ErgoSense® handle (Ø 13.7 mm)



Handle options.



LMErgoSense®



LMErgoMax™



LMErgoNorm™

Round Excavator, 1.0 mm

LMErgoSense®


61-62ES

Design:

- Round working end
- Curved cutting edge

Pirkimo dalis 21

The innovative lamella shape and the individual adaptation of the material composition to the respective materials are the great strengths of this system. The two-layer lamellae interlock during polishing. This makes TWIST polishers very thin, enabling precise polishing even in hard-to-reach areas such as the occlusion and interdental areas. They adapt to any tooth surface, preserve the original surface structure, and are characterized by low wear. With WST shank. Speed 10-12,000, max. 20,000 rpm. Sterilizable and reusable.



Diacomp Plus Twist composite polisher set, Ø 11 mm




Diacomp Plus Twist composite polisher set, Ø 14 mm

SDI | RIVA STAR

TOOTH DESENSITIZING AGENT

NEW GENERATION DESENSITIZER WHICH IMMEDIATELY REMOVES TOOTH SENSITIVITY

A STAR PRODUCT CONTAINING THREE POWERFUL COMPONENTS:



21 pirkimo dalis

TWO YEARS AND BEYOND DESENSITIZING EFFECT

Clinical studies support the immediate effect of Riva Star and its long term effect of over two years.

INCREASING BOND STRENGTHS OF GICS

Pre-treating tooth surfaces with Riva Star is a great way to enhance bonding of glass ionomers to dentin.

PATENTED PROCEDURE

Unlike other silver fluoride systems, Riva Star's two step patented procedure minimizes the risk of staining. By applying the potassium iodine solution over the silver fluoride solution, a silver iodide precipitate is formed.

ASSISTS IN THE ELIMINATION OF BIOFILM¹⁻³

INHIBITS BIOFILM FORMATION¹⁻³

MAY ASSIST IN THE ALLEVIATION OF CARIES⁴

NEXT GENERATION SDF SYSTEM

SDI continues to innovate SDF solutions by releasing Riva Star Aqua, a next generation water-based silver fluoride solution that is also effective as an SDF solution, with additional patient benefits.

Similar to Riva Star (SDF), Riva Star Aqua (AgF) is also a non-invasive 2 step patented system used to desensitise tooth pain in a few drops. Riva Star Aqua is an aqueous silver fluoride solution without the Ammonia base.

Riva Star Aqua Bottle kit

- 1 x 1.5mL Riva Star Aqua Bottle Step 1 (Blue label)
- 1 x 3mL Riva Star Aqua Bottle Step 2 (Green label)

8800531

