

OptraGate®

Moisture control made easy

OptraGate® is a latex-free lip and cheek retractor. The aid's unique flexibility and elasticity make it comfortable to wear and assist patients to keep their mouth open.

OptraGate
for popular relative isolation

16-17 pirkimo dalys

OptraDam Plus
for absolute isolation



Photos: courtesy of Ivoclar Vivadent AG, Liechtenstein

ivoclar | 100 years

- OptraGate White -
 - Assortment (40 x regular, 40 x small)
 - Refills x 80 (regular, small or junior)
 - Refills x 40 (regular, small or junior)

ProTaper Universal Hand



Files

Pirkimo dalys 18,19

Product Highlights

- Nickel titanium construction
- ISO-colored handles
- Sold as pack of 6



Taper

Assorted .02v-.09v

.07v

.08v

.09v

.06v

.04v

.02v

HAND INSTRUMENT STERILIZATION:

ProFile®, ProFile® Series 29®, ProFile® GT®, ProTaper®, Lexicon® Hand Instruments

Endo Instrument Sterilization:

Files must be sterilized before use. ANSI/ADA Specification 28 recommends:

- Scrub the instruments with soap and warm water.
- Rinse thoroughly with distilled or deionized water.
- Allow to air dry.
- Place the instruments, unwrapped, in the autoclave tray.
- Use fresh distilled or deionized water.
- Steam Autoclave at 136° C (plus or minus 2° C) for 20 minutes.