


Nuflow[®]

Nasal cannula set



Available sizes:

| | 1 Nasal cannula set small | 2 Nasal cannula set medium | 3 Nasal cannula set large | 4 Nasal cannula set xlarge |
|----------------------------|----------------------------------|-----------------------------------|----------------------------------|-----------------------------------|
| REF | 1320 | 1330 | 1340 | 1350 |
| Flow limits ¹ | 10,0 l/min | 14 l/min | 23 l/min | 27 l/min |
| Outer diameter prong | 2,3 mm | 2,7 mm | 3,2 mm | 3,7 mm |
| Septum distance | 3,7 mm | 4,3 mm | 4,8 mm | 5,3 mm |
| Patient group ² | Premature | Premature, Newborns | Infants | Children |
| Packing Unit | 10 pieces | 10 pieces | 10 pieces | 10 pieces |
| Fixation patches | 1380 small | 1380 small | 1381 large | 1381 large |

1 - See user manual for details

2 - Depends on size of patient nose

Connectors

Tube connection: OD 10 mm or OD 15 mm available

Material

- Hydrocolloid, Silicone, Zylar, TPE, PP,
- free of PVC, DEHP, natural rubber latex



Consumables

| | |
|----------|---|
| REF 1380 | Fixation patch, small (for use with Nuflow size REF 1320 or REF 1330) |
| REF 1381 | Fixation patch, large (for use with Nuflow size REF 1340 or REF 1350) |

Combinable devices

medinSINDI®, medinCNO®mini, medin-NC3®, medinBLENDER
Ventilators with high-flow mode and overpressure protection

Regulations

| | |
|---------------------------|---|
| Classification | EU IIa |
| Country of origin | Germany |
| Intended use | single patient use, non invasive, max 7 days |
| Manufacturing surrounding | clean room |
| Disposal | The Nuflow Nasal Cannula and its packaging can be discarded with regular household waste. |

Manufacturer:
medin Medical Innovations GmbH
Adam-Geisler-Straße 1
☎ +49 (0)8142 448 46 -0
info@medin-medical.com
www.medin-medical.com