



Recumbent wheelchair for transporting patients JORDAN A2h



Technical parameters:

1. Purpose:

wheelchair with increased maneuverability (adapted for narrow rooms and narrow passageways), suitable for transporting children and adults in a lying position.

2. Stroller frame:

1. Metal, covered with a coating that protects against corrosion, is resistant to mechanical effects and the chemical effects of disinfectants;

2. With protection ("bumpers") in all four corners.

3. Chassis:

1. Mobile, with four maneuverable wheels with a diameter of 125 mm;

2. Central wheel brake.

4. Adjusting the height of the lying part:

1. Hydraulic;

2. Height adjustment limits from 55 cm to 95 cm;

3. Height adjustment pedals are installed on both sides of the stroller.

5. Protective frames on the sides of the stroller:

1. Made of steel, covered with a coating that protects against corrosion, is resistant to mechanical effects and the chemical effects of disinfectants (surfaces are polished, welding seams are mechanically cleaned);

2. Can be lowered/raised, fixed in the upper position.

6. Stroller handles:

1. Equipped in both the head and foot end of the stroller, suitable for pushing/pulling the stroller; (one handle is visible in the catalog, but at the request of the buyer, the handles will be made on both sides of the cart).

2. Made of steel, covered with a coating that protects against corrosion and is resistant to the mechanical and chemical effects of disinfectants (surfaces are polished, welding seams are mechanically cleaned).

7. Mattress:

1. The length and width of the mattress correspond to the base of the mattress;

2. Mattress height 5 cm thick;

3. The mattress is covered with antistatic, wet cleaning and disinfectant resistant, non-absorbing artificial leather.

8. Dimensions of the lying part of the stroller (mattress base): Length 200 cm, width 70 cm;

9. External (overall) dimensions of the stroller: Length 217 cm, width 83 cm;

10. Carrying capacity of the trolley (the maximum operating load of the trolley provided by the manufacturer) 250 kg;

11. Warranty period 24 months.

The logo for RESI, featuring the word "RESI" in a bold, sans-serif font with a stylized 'R'.

RESI Třeboň spol. s r.o. , Novohradská 1153, 379 01 Třeboň, Czech Republic
Tel.: +420 384 724 292, www.resi.cz

