

Peties sąnario pozicionavimo sistema su paklos laikikliu

Patient Positioning

T-MAX System (beach chair positioner)



The T-MAX Shoulder Positioner is a patient positioning system that provides arthroscopy surgeons with the full range of features they need.

Exposure – shoulder pad slides laterally, providing maximum exposure to the operative shoulder

Access – patient can be positioned closer to the side of the OR table to allow maximum accessibility

Adjustability – hydraulic assist mechanism allows any angle from sitting to supine (0-90°)

Versatility – the T-MAX attaches to most OR tables and is designed for intraoperative use of radiographs and fluoroscopy

4. Universalus sistemos pajungimo mechanizmas, leidžiantis ję prijungti prie bet kurio operacinio stalo.

Reference #	Description
7210996	T-MAX System

System includes

✓ 7210551	T-MAX Shoulder Positioner	Pozicionavimo blokas.
✓ 7210552	Adjustable Arm Board	Reguliuojamas rankos laikiklis.
✓ 7210553	T-MAX Cart	Pozicionavimo bloko vežimas.
✓ 7210554	Square Rail Clamps (2 per box)	Bėgų ventžiklis
✓ 7210555	Intubation Pad	Laikiklis palvairi
7210560	Knee Bolster	
72201184	Body Strap	

2. Sistema turi hidraulinį mechanizmą, kuris gali kisti paciento padėtį nuo gulimos iki sėdimos.



Tech Support, Endoscopy <Endoscopy.TechSupport@smith-nephew.com>
 Kam Goralska, Katarzyna Tech Support, Endoscopy
 Kopija Greta Vilkaite

11615 PROTIS
 ↓
 AUKŠTIS

S&N #	TENET #	PRODUCT DESCRIPTION	Weight (lbs.)	L (in.)	W (in.)	H (in.)	Weight (kg)	L (cm)	W (cm)	H (cm)	Description
77210551	7401	T-MAX II Shoulder Exposure Positioner	49.2	14	18.75	14.5	22.4	77.5	47.5	37	Chair like device / attach to OR Table
77210552	7411	Adjustable Arm Board	5.6	21.25	6	3	2.5	54	15.5	7.5	Board to Position non-operative arm
77210553	7412	T-MAX Cart	8.4	42.75	7.75	6.5	3.8	109	20	17	Cart for transport and storage of T-MAX

SmithNephew

Michael Oxner | Technical Support
 Michael.Oxner@smith-nephew.com
 Office: +1.800.343.5717
 www.smith-nephew.com

From: Goralska, Katarzyna <Katarzyna.Goralska@smith-nephew.com>
 Sent: Monday, December 22, 2025 7:52 AM
 To: Tech Support, Endoscopy <Endoscopy.TechSupport@smith-nephew.com>
 Cc: Greta Vilkaite <RV@ostecal.lt>
 Subject: T-MAX - dimenzions

Dear Colleagues,

I would like to kindly request your support in relation to our distributor's in Lithuania ongoing tender process. The currently available documentation does not include these measurements which are necessary for our technical review. We would greatly appreciate your assistance and thank you in advance for your prompt support.

Patient Positioning

SPIDER2 Limb Positioner and Accessories

SPIDER2 System (T-MAX)

Reference #	Description
72203565	T-MAX System

System includes

✓ 7210551	T-MAX Shoulder Positioner
✓ 7210552	Adjustable Arm Board
✓ 7210553	T-MAX Cart
✓ 7210554	Square Rail Clamps (2 per box)
✓ 7210555	Intubation Pad
72203299	SPIDER2 Limb Positioner
7210570	Piggy Back Connector (2 per box)
7210571	Shoulder Connection Bar (2 per box)
7210576	Lateral Shoulder Connection Bar (2 per box)
7210577	Traction Accessory
72200918	SPIDER Cart
72203300	Switch Drape (20 per box)
72203301	Battery Pack

SPIDER2 System (Beach Chair)

Reference #	Description
72203516	Beach Chair System

System includes

72203299	SPIDER2 Limb Positioner
7210570	Piggy Back Connector (2 per box)
72200918	SPIDER Cart
72203300	Switch Drape (20 per box)
7210571	Shoulder Connection Bar (2 per box)

SPIDER2 System (Lateral Decubitus)

Reference #	Description
72203517	Lateral Decubitus System

System includes

72203299	SPIDER2 Limb Positioner
7210570	Piggy Back Connector (2 per box)
7210576	Lateral Shoulder Connection Bar (2 per box)
7210577	Traction Accessory
72200918	SPIDER Cart
72201428	Lateral JACK
72202466	Lateral JACK Sterile Cover (25 per box)
72203300	Switch Drape (20 per box)

Miscellaneous

Reference #	Description
72200918	SPIDER Cart
7210553	T-MAX Cart

Disposables

Reference #	Description
7210573	Shoulder Stabilization Kit, disposable, sterile (12 per case)
7210559	Disposable Face Mask, disposable (8 per case)

Patient Positioning

T-MAX Shoulder Positioners and Accessories

Combination System

Reference #	Description
7210997	T-MAX and SPIDER Combination System

System includes

7210551	T-MAX Shoulder Positioner
7210552	Adjustable Arm Board
7210553	T-MAX Cart
7210554	Square Rail Clamps (2 per box)
7210555	Intubation Pad
7210569	SPIDER Limb Positioner
7210570	Piggy Back Connector (2 per box, required for 7210569)
7210571	Shoulder Connection Bar (2 per box)
7210576	Lateral Shoulder Connection Bar (2 per box)
7210577	Traction Accessory, two included in system
72200918	SPIDER Cart

T-MAX Shoulder Positioner and Accessories

Reference #	Description
7210551	T-MAX Shoulder Positioner
7210552	Adjustable Arm Board
7210553	T-MAX Cart
7210554	Square Rail Clamps (2 per box)
7210555	Intubation Pad
7210556	Gel Face Mask <i>gelinei bouke</i>
7210557	Extended Lateral Support
7210558	Mask Straps (pair)
7210559	Disposable Face Mask, disposable (8 per box)
7210560	Knee Bolster
7210563	Intubation Drape Bar
7210564	Lateral Support Knobs, t-shaped (2 per box)
72200481	T-MAX, Japanese/Eschman
72200482	T-MAX, OPT
72200668	T-MAX, DE
72201185	T-MAX, MA
72203318*	Adjustable Lateral Support
72203415*	Non-operative Arm Board




* Item shown

T-MAX Beach Chair Positioner

Operating with confidence, the T-MAX Beach Chair Positioner affords the entire operating room staff the tools needed to position patients for beach chair shoulder arthroscopy and open procedures.⁹ Offering a dual ball and socket head rest, adjustable lateral supports, sliding shoulder pads and more, the T-MAX is a unique beach chair positioner that benefits surgical staff, anesthesiologists and surgeons.^{9, 24, 25} Equipped with adjustable table legs and a variety of rail clamps, the T-Max fits a wide array of OR table offerings.

Hidnoulinis
mechanizmas

Hydraulic assist
mechanism



Access

Optimize access to the repair site. T-MAX slides laterally on the frame to either side of the OR table bringing the patient closer to the surgeon. The sliding shoulder pad laterally adjusts with ease to maximize exposure or fully support the posterior aspect of the operative shoulder.⁹

Versatility

Suitable for arthroscopic and total shoulder procedures, our T-MAX Beach Chair Positioner attaches to a variety of OR tables and offers head rest movement and an intubation tree, well arm positioning and body support.

Stability

Securely holds the patient while positioning and re-positioning between 0 and 90°. Built-in lateral supports are designed to provide sturdy counter traction when distracting the operative arm.

Intubation tree

Used to hold intubation tubing and can simply be rotated out of the way if not required

Head rest

Offers a wide range of positions of tilt and rotation, maintaining the patient's head and airway in a static position throughout the procedure

Lateral support pads

Designed to ensure patient stability while providing counter traction when distracting the operative arm

Lateral adjustment

Allows the T-MAX to slide laterally on the frame, positioning patients close to the surgeon

Spider rail

The only beach chair positioner designed to work with the SPIDER Limb Positioner

Veido Rauts

Face mask

The disposable face mask allows attachment to head rest

su tvintinimo septukais

Shoulder pad

Sliding shoulder pad and posterior placement can provide stability and allow for unrestricted access to the surgical site⁹

Non-operative arm board (static)

Provides support and adjustability for the non-operative arm

Hidraulinis mechanizmas

Hydraulic assist mechanism

Assists in raising and lowering patient

Knee bolster

Assists in preventing the patient from sliding down the table

Mobilios jungtys, ipoline sistema lengvai keisti (papal pority) pozicionerimo z' laimis' an desinimi.