

ISE Cleaning Solution / Elecsys SysClean

cobas[®]

11298500 316

5 x 100 mL

English**Intended use**

For the cleaning of ISE units on Roche/Hitachi analyzers.
For the cleaning of Elecsys analyzers.

Summary

ISE Cleaning Solution / Elecsys SysClean is an alkaline cleaning solution with antimicrobial properties.

Reagents – contents and concentrations

5 bottles, each containing 100 mL

Components:

Sodium hydroxide 3 mol/L

Sodium hypochlorite solution (< 2% active chlorine)

Additive


Precautions and warnings

For in vitro diagnostic use

Exercise the normal precautions required for handling all laboratory reagents.

Safety data sheet available for professional user on request.

Warning. This kit contains components classified as follows according to the European Directive 88/379 EEC:

 C – Corrosive (sodium hydroxide) R35; S26 – 37/39 – 45:

Causes severe burns. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable gloves and eye/face protection. In case of accident, or if you feel unwell, seek medical advice immediately (show the label where possible).

Contact phone: all countries: +49-621-7590; USA: +1-800-428-2336.

Disposal of all waste material should be in accordance with local guidelines.

Handling

The solution is ready for use.

Storage and stability

Store protected from light.

The solution is stable up to the stated expiration date when stored at 2–8°C.

Materials provided

ISE Cleaning Solution / Elecsys SysClean

Materials required (but not provided)

- General laboratory equipment
- Roche/Hitachi or Elecsys analyzer

Instructions for use

Please refer to the operator's manual and maintenance instruction sheets.

ELECSYS and COBAS are trademarks of Roche.
Other brand or product names are trademarks of their respective holders.
Significant additions or changes are indicated by a change bar in the margin.
©2008 Roche Diagnostics



Roche Diagnostics GmbH, D-68298 Mannheim

