

Fix-clip

Fix-clip Cat. No. FS 0030

Intracutaneous suture with Fix-clip for an optimal cosmetic result.

Pierce the skin at two points only


No suture knots

Fixation with two clips

Easy to remove after wound healing: simply cut off one of the clips

Easy suture removal by pulling out the second clip

For threads EP 1-2 (USP 5/0-3/0)



UK03
201

Special Products

Minimal invasive surgery

| Material | EP | USP | Length | Cat. No |
|--|-----|-----|------------|---------|
| CHIRASORB violet, Loop + Applicator | 3 | 2/0 | 1 x 50 cm | MS 0143 |
| CHIRASORB violet, Loop + Applicator | 3.5 | 0 | 1 x 50 cm | MS 0144 |
| CHIRASORB violet, Loop + Applicator, KR 27 | 4 | 1 | 1 x 120 cm | MS 1125 |
| CHIRLAC violet, Loop + Applicator | 4 | 1 | 1 x 70 cm | MS 0125 |
| CHIRLAC violet, Loop + Applicator | 5 | 2 | 1 x 70 cm | MS 0126 |
| 5.1 POLYDOX violet, Loop + Applicator | 3.5 | 0 | 1 x 70 cm | MS 0134 |
| POLYDOX violet, Loop + Applicator | 4 | 1 | 1 x 70 cm | MS 0135 |
| SILON monofilament blue, Loop + Applicator | 3.5 | 0 | 1 x 70 cm | MS 1124 |