


Ruby™ Case Examples

38 mm Hepatic Artery Aneurysm



Pre-Treatment Post-Treatment

Dr. Jim Benenati
Baptist Cardiac and Vascular Institute, Miami, FL


Pulmonary AVM



Pre-Treatment Post-Treatment

Dr. J Moskovitz
Florida Hospital, Orlando, FL


Pre Y90 GDA Embolization



Pre-Treatment Post-Treatment

Dr. Jie Mao
Boulder Community Hospital, Boulder, CO


Hypogastric Embolization



Pre-Treatment Post-Treatment

Dr. Marcus D'Ayala
New York Methodist Hospital, Brooklyn, NY

Type 2 AAA Endoleak



Pre-Treatment Intra-Operative Post-Treatment

Dr. Frank Arko
Sanger Heart and Vascular Institute, Charlotte, NC

Bilateral Ovarian Vein Embolization



Post-Treatment


Dr. Anthony Venbrux
George Washington Hospital, Washington, DC

PERIPHERAL

Ruby™ Coil

Widest range of sizes to treat all anatomies.

Ordering Information




Frame			Fill		
COMPLEX STANDARD			COMPLEX SOFT		
Catalog Number	Secondary Diameter (mm)	Length (cm)	Catalog Number	Secondary Diameter (mm)	Length (cm)
RB2C0305	3	5	RB4C0201	2	1
RB2C0312	3	12	RB4C0202	2	2
RB2C0320	3	20	RB4C0204	2	4
RB2C0410	4	10	RB4C0305	3	5
RB2C0415	4	15	RB4C0315	3	15
RB2C0420	4	20	RB4C0406	4	6
RB2C0430	4	30	RB4C0415	4	15
RB2C0435	4	35	RB4C0620	6	20
RB2C0512	5	12	RB4C0630	6	30
RB2C0520	5	20	RB4C0835	8	35
RB2C0530	5	30	RB4C0860	8	60
RB2C0612	6	12	RB4C1035	10	35
RB2C0620	6	20	RB4C1240	12	40
RB2C0630	6	30	RB4C1650	16	50
RB2C0725	7	25	RB4C2060	20	60
RB2C0740	7	40			
RB2C0825	8	25			
RB2C0840	8	40			
RB2C0930	9	30			
RB2C1035	10	35			
RB2C1240	12	40			
RB2C1260	12	60			
RB2C1440	14	40			
RB2C1460	14	60			
RB2C1660	16	60			
RB2C1857	18	57			
RB2C2060	20	60			
RB2C2457	24	57			
RB2C2860	28	60			
RB2C3260	32	60			


FEATURED IN
Endovascular
TODAY

Technical Tips
and Clinical Benefits
of the Penumbra
Ruby™ Embolization Coil

April 2013



Catalog Number	Tip Shape	ID	Distal	Proximal
PXSLIMSTR	Straight Tip	.025"	2.6F	2.95F
PXSLIM045	45° Tip	.025"	2.6F	2.95F
PXSLIM090	90° Tip	.025"	2.6F	2.95F
PXSLIM130	130° Tip	.025"	2.6F	2.95F



Catalog Number	Description
RH1	Detachment Handle
RH5	Detachment Handle (5 pack)

1. Schönholz C, Benenati J, Arko F, Teigen C, Moskovitz JD. Ruby Detachable Coils—Case reports illustrating the utility of this novel technology. Endovascular Today. 2013;12(4):78-82.



Penumbra, Inc. USA
1351 Harbor Bay Pkwy
Alameda, CA 94502
USA
1.888.BRAIN.06
1.888.272.4606
T 1.510.748.3200
F 1.510.748.3232
order@penumbrainc.com
info@penumbrainc.com


Penumbra Europe GmbH
Am Borsigturm 44
13507 Berlin
Germany
T +49 30 2005 676-0
F +49 30 2005 676-10
order@penumbrainc.de
info@penumbrainc.de

Penumbra Neuro Australia PTY LTD
Suite 3, Level 5, 1 Oxford Street
Darlinghurst NSW 2010
Australia
T +61-1300 817 025
F +61-1300 817 026
order.anz@penumbrainc.com

PERIPHERAL

Ruby™ Coil

Large Volume Detachable Coils



High Flow Microcatheter
Compatible



Ruby™ Coil

Large Volume Detachable Coils

VOLUME ADVANTAGE

Reduces the number of coils and procedure time by maximizing length and coil thickness.¹

30 cm Ruby™ Coil


30 cm Conventional Detachable Coil



Only Coil to Provide 60 cm Length and .020" Thickness

Ruby™ Coil

High Flow Microcatheter
Penumbra PX SLIM™




60 cm
Max. Length

Coil Thickness

.020"

Conventional 35 Coil
035 Diagnostic Catheter



40 cm
Max. Length

.021"

Conventional 18 Coil
18 Microcatheter



30 cm
Max. Length

.012"

COMPLETE CONTROL

Fully retractable coil with active mechanical detachment delivers complete control.

Instant 1-Click Detachment



Standard
32 mm x 60 cm




Soft
2 mm x 1 cm

EXTREME DELIVERABILITY

Only coil designed to be delivered through a high flow microcatheter.

PX SLIM™ Microcatheter



SOFTNESS

Choice of Standard and Soft coils provide superior packing density while avoiding catheter kick-out.

Unique Tri-Layer Coil Design

