

Angio-Seal™

Vascular Closure Device

185 pirkimo dalis



The Angio-Seal™ Vascular Closure Device uses **three bioabsorbable components** to actively seal the arteriotomy

- **Anchor**

Intravascular: placed against the inside of the vessel wall.

- **Collagen**

Extravascular: placed on top of the arteriotomy in the tissue tract

- **Suture**

Cinches the anchor and collagen together to form a secure seal

COLLAGEN

*9-Hole Weave
V-Twist Folding
Pattern*

ANCHOR

*50:50 Blend of lactide
and glycolide polymers*

SUTURE

*Poly-glycolic Acid
(PGA) coated*

