

ETHICON

PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

CONTOUR™ Curved Cutter Stapler

Designed for deeper reach
and better visibility



Introducing ARC

Answering surgeons' needs

By aligning our devices to the surgical treatment pathway, we can ensure you have the right solutions for the right procedure.

Entry
Dissection
Ligation

Access



Transection
Resection
Anastomosis

Repair



Layer closure
Port closure

Closure



**COLORECTAL
THORACIC
SKELETAL MUSCLE**

The Contour™ Curved Cutter Stapler is a multifire, single patient use device with a curved head, that cuts and staples simultaneously. The device allows placement of a 40-mm cut line in the width of a 30-mm space

Key surgical concerns

- Placement below Tumor line in colorectal surgery
- Access in hard-to-reach tissues low within the rectum
- Placement in tight confines of the pelvic region

Main device benefits


- Unique curved head design for lower pelvic access & enhanced visibility of the anatomic structures¹
- Facilitates maneuverability and mobility to place staples on the colon just where you want them
- One-handed firing & closure, free the other hand for use in pelvis
- Intermediate locking position allows tissue manipulation and inspection before closing device. Manual pin placement option captures intended tissue before closing device





1 Unique curved head design

Designed for lower pelvic access and enhanced visibility¹; allows placement of a 40-mm cut line in the width of a 30-mm space²



2 Manual pin placement option

Captures intended tissue before closing device




3 Low pelvic access¹

Delivers low pelvic access with enhanced visibility



4 Intermediate locking position

Allows tissue manipulation and inspection before closing device



5 Reloadable cartridges

A single device may be fired up to 6 times in a single procedure

6 Simultaneous stapling and cutting

Delivers 4 curvilinear rows of staples with a single cut between—eliminates² procedure steps the need for a bowel clamp and a scalpel⁴

7 Ergonomic design

One-handed firing, pin placement and closure free the other hand for use in pelvis³

8 Audible and tactile feedback⁵

Breakthrough washer provides audible and tactile feedback with every firing and a new stainless steel anvil with every cartridge - no picture



Ordering CONTOUR™ Curved Cutter Stapler:

	Code	Description	Qty/ Sales Unit
Devices	CS40B	CONTOUR™ Curved Cutter Stapler with blue cartridge	3
	CS40G	CONTOUR™ Curved Cutter Stapler with green cartridge	3
Cartridges	CR40B	CONTOUR™ Reload Cartridge - blue	6
	CR40G	CONTOUR™ Reload Cartridge - green	6

¹ C0000001342. Data on file, Ethicon Endo-Surgery (Europe) GmbH ² C0000001404. Data on file, Ethicon Endo-Surgery (Europe) GmbH ³ C0000001347. Data on file, Ethicon Endo-Surgery (Europe) GmbH
⁴ C0000001341. Data on file, Ethicon Endo-Surgery (Europe) GmbH ⁵ C0000001348. Data on file, Ethicon Endo-Surgery (Europe) GmbH

Ethicon Endo-Surgery (Europe) GmbH
Hummelsbütteler Steindamm 71
22851 Norderstedt, Germany

www.ethicon.com

BR 596.1 ©2013 Ethicon Endo-Surgery (Europe) GmbH.
Please always refer to the Instructions for Use / Package Insert that come with the device for the most current and complete instructions.

ETHICON
PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES