

**Nuo:** SecondLevel, TM <secondlevelih@orthoclinicaldiagnostics.com>  
**Išsiųsta:** trečiadienis 2022 m. lapkritis 23 10:08  
**Kam:** VITROLAB; SecondLevel, TM  
**Tema:** RE: Reagents refilling

TS 8 p.

Dear

If by refill you mean loaded reagent vials/cassettes, yes it is possible at anytime to load reagent

Thanks & Regards

**Florence Ragué**  
Staff Engineer  
Ortho Care  
M: +33 604 65 49 67  
[florence.rague@OrthoClinicalDiagnostics.com](mailto:florence.rague@OrthoClinicalDiagnostics.com)

[OrthoClinicalDiagnostics.com](http://OrthoClinicalDiagnostics.com)

*Quidel and Ortho are becoming QuidelOrtho [www.quidelortho.com](http://www.quidelortho.com)*  
QDEL | Nasdaq

**Ortho Clinical Diagnostics**  
Because Every Test Is A Life

The information contained in this transmission is intended only for the person or entity to which it is addressed and may contain confidential, trade secret, proprietary, privileged, private and/or copyrighted material. If you are not the intended recipient(s) of this message, or the employee or agent responsible for delivery of this message to the intended recipient(s), do not review, retransmit, disclose, disseminate, use, or take any action in reliance upon, this information. If you received this transmission in error, please contact the sender and destroy all printed copies and delete the material from all computers. In connection with our business, QuidelOrtho may collect and process your personal data. For further information regarding how we use this data, please see our online [privacy notice](#).

**From:** : | VITROLAB <konkursai@vitrolab.lt>  
**Sent:** Wednesday, November 23, 2022 7:44 AM  
**To:** SecondLevel, TM <secondlevelih@orthoclinicaldiagnostics.com>  
**Subject:** Reagents refilling  
**Importance:** High

**EXTERNAL SENDER: Verify links, attachments and sender before taking action**

Dear All,

Could you confirm that the Ortho Vision Analyser' reagents can be refilled at **any time** during the system operation?

Public tender administrator



Baltai pr. 36-11, LT-48196 Kaunas, Lithuania | Phone/fax.: 8-37 33 33 17/29  
Mob. phone: +370 640 10809 | E-mail: [biuras@vitrolab.lt](mailto:biuras@vitrolab.lt) | [www.vitrolab.lt](http://www.vitrolab.lt)