



## CRYOFIX GEL

Embedding matrix for cryostat sectioning in different colors

Cryofix Gel matrix contains mixture of polymers soluble in water. The polymers' freezing point is ideal for cryostat sectioning. Stained matrices enable easier visualization of embedded samples and easier organization of work.

- CryoFix Gel
- CryoFix Gel pink
- CryoFix Gel blue
- CryoFix Gel green

packaging	product code
100 mL	CF-100
100 mL	CFR-100
100 mL	CFP-100
100 mL	CFZ-100

181 pirkimo dalis

## CRYO SPRAY

Medium in spray bottle for freezing histological samples

Cryo Spray enables quick freezing of histological samples prior to sectioning on cryostat and microtome.

- Cryo Spray

packaging	product code
150 mL	08-SPRAY

## DEMOLD SPRAY

Solution in spray bottle for easy removal of paraffin tissue blocks

DeMold Spray is applied to the mold prior to placing the samples and pouring the wax, allowing easy and quick removal of tissue blocks from the mold.

- DeMold Spray

packaging	product code
100 mL	DS-100

## DEMOLD SOLUTION

Solution for refill of DeMold Spray spray bottle

DeMold Spray is applied to the mold prior to placing the samples and pouring the wax. It allows easy and quick removal of tissue blocks from the mold.

- DeMold solution
- DeMold solution

packaging	product code
500 mL	DSO-500
1000 mL	DSO-1L

