

BioGram 4 kit

Kit for identification of bacteria according to Gram

Differential staining of bacteria according to Gram is widely known and used method in microbiology. BioGram 4 kit provides optimal staining and identification of Gram-positive and Gram-negative bacteria.

BioGram 4 kit

Four-reagent kit for identification of bacteria according to Gram. Kit contains Gram Crystal Violet 1% solution, stabilized Gram Lugol solution, double amount of Gram Decolorizer solution 2 and Gram Safranin solution as a counterstain.

5x100 mL	BGR4-K-100
5x250 mL	BGR4-K-250
5x500 mL	BGR4-K-500

Gram Crystal Violet 1% solution

Primary stain reagent for Gram-positive bacteria.

100 mL	GC1-OT-100
250 mL	GC1-OT-250
500 mL	GC1-OT-500
1000 mL	GC1-OT-1L
2500 mL	GC1-OT-2.5L

Gram Lugol solution, stabilized

Solution contains polyvinylpyrrolidone iodine (PVP-I).

100 mL	GLS-OT-100
250 mL	GLS-OT-250
500 mL	GLS-OT-500
1000 mL	GLS-OT-1L
2500 mL	GLS-OT-2.5L

Gram Decolorizer solution 2

Medium speed decolorizer used after stabilized iodine solution.

100 mL	GD2-OT-100
250 mL	GD2-OT-250
500 mL	GD2-OT-500
1000 mL	GD2-OT-1L
2500 mL	GD2-OT-2.5L

Gram Safranin solution

Counterstain for Gram-negative bacteria, contains Safranine O dye.

100 mL	GSF-OT-100
250 mL	GSF-OT-250
500 mL	GSF-OT-500
1000 mL	GSF-OT-1L
2500 mL	GSF-OT-2.5L

55.7 dalis

