

		Catheter								Introducer				Stylet
Code		Int. Ø mm	Ext. Ø mm	Gauge G	Length cm	Flow rate by gravity ml/min	Prim. vol. ml	Ø Fr	Tube Hardness	Length mm	G G	Type	Ext. Ø mm	-
	<u>1252.030</u>	0.3	0.6	24	30	5.0	0.12	2	Standard	18	20	Microflash	1.1	Yes
	<u>1252.15</u>	0.3	0.6	24	15	9.0	0.06	2	Standard	25	20	Splitting needle	1.0	No
	<u>1252.150</u>	0.3	0.6	24	15	9.0	0.06	2	Standard	18	20	Microflash	1.1	No
11.2.	<u>1252.30</u>	0.3	0.6	24	30	5.0	0.12	2	Standard	25	20	Splitting needle	1.0	No
	<u>1252.31</u>	0.3	0.6	24	30	5.0	0.12	2	Standard	-	-	-	-	No
	<u>1252.35</u>	0.3	0.6	24	30	5.0	0.12	2	Standard	18	20	Microflash	1.1	No
	<u>1253.15</u>	0.5	1.0	20	15	5.0	0.09	3	Standard	30	17	Splitting needle	1.5	No
	<u>1253.30</u>	0.5	1.0	20	30	2.8	0.16	3	Standard	30	17	Splitting needle	1.5	No
	<u>1254.30</u>	0.8	1.3	18	30	18	0.42	4	Standard	30	17	Splitting needle	1.5	No
	<u>1262.30</u>	0.3	0.6	24	30	5.0	0.12	2	Soft	25	20	Splitting needle	1.0	No
	<u>1262.35</u>	0.3	0.6	24	30	5.0	0.12	2	Soft	18	20	Microflash	1.1	No
	<u>1352.152</u>	0.3	0.6	24	15	9.0	0.06	2	Standard	18	20	Microflash	1.1	No

	Catheter								Introducer				Stylet
Code	Int. Ø mm	Ext. Ø mm	Gauge G	Length cm	Flow rate by gravity ml/min	Prim. vol. ml	Ø Fr	Tube Hardness	Length mm	G G	Type	Ext. Ø mm	-
<u>1352.302</u>	0.3	0.6	24	30	5.0	0.12	2	Standard	18	20	Microflash	1.1	No
<u>1353.302</u>	0.5	1.0	20	30	2.8	0.16	3	Standard	45	17	Peelable cannula	1.5	No

- **Contains Latex:** No
- **Contains DEHP:** No
- **Contains products of animal or biological origin:** No
- **Pyrogen-free:** Yes